

Medical + Optical

Caring about bringing quality to you



Surgistar

The Cutting Edge of Tomorrow™

Surgical Cutting Instruments

Surgistar

Quality. Performance. Value.

Expect this and more from Surgistar.

Surgistar has been manufacturing single-use Ophthalmic, ENT, Surgical and Dental cutting instruments for over twenty years. Their products are manufactured using the highest levels of care and precision at their manufacturing facility in Vista, California, and are shipped through their packaging center now also located in Vista, California.

Selection. Service. Innovation.

Their portfolio numbers over 600 items, featuring the widest range of microsurgical blades, knives, cannulas, trephines and accessories offered by any manufacturer.

They offer single-use products for all types of procedures, including but not limited to:

Cataract, Refractive, Vitreoretinal, Oculoplastic, Myringotomy, Cardiology, Dermatology, Arthroscopy, Hair Transplant, Gynecological and Periodontal.



Vista, CA Manufacturing Facility

Contents

Angled Incision Knives & Collets	Page 4
Mini Blades	Page 5
Other Knives & Blades	Page 5
Nasal Sinus Series Round Stock	Page 5
Handles	Page 6
Tympanplasty Flat Stock	Page 6
Tapered Series Round Stock	Page 6
Myringotomy Flat Blades	Page 7
Myringotomy Round Blades	Page 7
Shapeshifter Malleable Dental Range	Page 8

We can help you with your requirements for single-use ophthalmic, ENT, surgical and dental cutting instruments.


Call your Territory Manager to discuss the options and configurations available, alternatively call our main office.

Toll Free: 1800 225 065


Angled Incision Knives & Collets

Product Code	Product Description		Qty/Bx
38-961501	Stab & Paracentesis Knife Complete 15°		Box 6
38-961502	Collet 15°		Box 6
38-962201	Stab & Paracentesis Knife Complete 22.5°		Box 6
38-963001	Stab & Paracentesis Knife Complete 30.4°		Box 6
38-964501	Stab & Paracentesis Knife Complete 45°		Box 6

Angled Incision Knives & Collets Controlled Depth

Product Code	Product Description		Qty/Bx
38-7511	Collet 1.5mm Depth 15°		Box 6
38-7513	Collet 3.0mm Depth 15°		Box 6
38-7514	Knife 3.0mm Depth 15°		Box 6
38-7515	Collet 5.0mm Depth 15°		Box 6
38-7516	Knife 5.0mm Depth 15°		Box 6
38-7522	Collet 4.0mm Depth 22.5°		Box 6
38-7523	Knife 4.0mm Depth 22.5°		Box 6
38-7530	Collet 3.5mm Depth 30°		Box 6
38-7531	Knife 3.5mm Depth 30°		Box 6
38-7536	Knife 5.0mm Depth 30°		Box 6

Mini Blades

Product Code	Product Description		Qty/Bx
38-6100	Miniature 1.5mm Blade Chisel Tip		Box 12
38-6200	Miniature 2.7m Blade Chisel Tip		Box 12
38-6300	Miniature Blade Spear Tip		Box 6
38-6400	Miniature Blade Round Tip, Sharp 1 Side		Box 12
38-6500	Miniature Blade Lance Tip		Box 12
38-6600	Lamellar Blade Angled Bevel Up		Box 6
38-6700	Miniature Blade Curved Tip, Sharp 1 Side		Box 12
38-6900	Miniature Blade Round Tip, Sharp All Around, Straight		Box 6
38-6910	Miniature Blade Round Tip, Sharp All Around, Angled		Box 6
38-6961	Micro Blade Round Tip, Sharp All Around, Straight		Box 6


Other Knives and Blades

Product Code	Product Description		Qty/Bx
38-5910	Needle Blade 3.0mm		Box 6
38-5920	Needle Blade 5.0mm		Box 6
38-5700	Scleral Hockey Stock Blade, Medium		Box 6
38-7300	Sickle Blade		Box 6

Nasal Sinus Series - Round Stock

Product Code	Product Description		Qty/Bx
38-8920	Nasal Sinus Knife Sickle Tip Straight		Box 4


Handles

Product Code	Product Description		Qty/Bx
38-1305	Collet Cross-Cut Only		Box 1
38-1311	Handle 5cm Round		Box 1
38-1312	Handle 7.5cm Round		Box 1
38-1310	Handle 10cm Round		Box 1
38-1350	Handle 13cm Round		Box 1
38-1360	Handle 13cm Round, Heavy Duty		Box 1
38-1320	Handle 15.5cm Round		Box 1
38-1330	Handle 10cm Hexagonal		Box 1


Tympanoplasty - Flat Stock

Product Code	Product Description		Qty/Bx
38-7200	Tympanoplasty Blade Angled 60° 2.5mm Wide		Box 6
38-7210	Tympanoplasty Blade Straight 2.5mm Wide		Box 6



Tapered Series - Round Stock

Product Code	Product Description		Qty/Bx
38-8230	Tapered Straight Edge Blade Angled Upward Cut		Box 6
38-8231	Tapered Straight Edge Blade Angled Downward Cut		Box 6
38-8232	Tapered Straight Edge Blade Straight		Box 6
38-8233	Tapered Spear Tip Blade Angled		Box 6
38-3234	Tapered Spear Tip Blade Straight		Box 6
38-1382	Tapered Series Round Stock ENT Handle 10cm		Box 1

Myringotomy - Flat Blades

Product Code	Product Description		Qty/Bx
38-7100	Lance Blade - Adult		Box 6
38-7100H	Lance Knife - Adult		Box 6
38-7110	Lance Reverse Blade - Adult		Box 6
38-7121	Lance Blade - Juvenile		Box 6
38-7124	Lance Knife - Juvenile		Box 6
38-7120	Spear Blade - Juvenile		Box 6
38-7122	Spear Knife - Juvenile		Box 6
38-7300	Sickle Blade		Box 6
38-7300H	Sickle Knife		Box 6

Myringotomy - Round Blades

Product Code	Product Description		Qty/Bx
38-700715	Sickle Knife Straight - Flat Handle		Box 6
38-700716	Lance Knife Straight - Flat handle		Box 6
38-8124	Lance Knife Straight - Juvenile		Box 6
38-8125	Lance Knife Angled - Juvenile		Box 6
38-8130	Spear Micro-Point Blade Straight		Box 6


Product Code	Product Description		Qty/Bx
38-S6300	Miniature Blade Spear Tip, Straight		6/Box
38-S6310	Miniature Blade Spear Tip, Angled 10°		6/Box
38-S6400	Miniature Blade Round Tip, Sharp 1 Side, Straight		12/Box
38-S6500	Miniature Blade Lance Tip 15°, Straight		12/Box
38-S6700	Miniature Blade Curved Tip, Sharp 1 Side, Straight		12/Box
38-S6900	Miniature Blade Round Tip, Sharp All Around, Straight		6/Box
38-S6915	Miniature Blade Round Tip, Sharp All Around, Angled 10°		6/Box
38-S6945	Miniature Blade Round Tip, Full Radius, Angled 45°		6/Box




Specially designed for use in Endodontic Microsurgery.



Narrow blades are 60% smaller than Standard, improving visibility under a microscope and allowing for more precise incisions.




Malleable shank can be contoured to your preference. Negotiate turns and difficult angles with ease.



Size comparison of narrow, standard and traditional scalpel.

Surgistar

The Cutting Edge of Tomorrow
www.surgistar.com



Medical + Optical Head Office
81 Whiting St
Artarmon NSW 2064
info@medopt.com.au
TOLL FREE: 1800 225 065
www.medopt.com.au
ABN: 80 003 178 596
2024