




OrthoMedics

Caring about bringing quality to you

Sterilization Baskets & Accessories




OrthoMedics

Caring about bringing quality to you



Medical + Optical is a privately owned company. Our aim is to provide intelligent solutions for today's Healthcare needs. We work closely with our major international suppliers to ensure our product offering is of the highest quality and meets the requirements of the end user. Our sterilization baskets are manufactured in Italy by the global company CBM and Medical + Optical are delighted to be their Australian distributors.




CBM was founded by Cav. Busatti in 1947, based originally in Milan and moving to its new, much larger site in Torre de' Picenardi (Cremona) in 1962. CBM conforms to all the European and Global Standards for manufacturing and is a significant member of the European Technical Committee CEN and the International Standards Committee ISO. This involvement ensures that CBM remains one of the leading manufacturers of sterilization baskets in the world. The CBM range

of sterilization baskets use the highest standard of raw materials, such as, AISI-304 Medical Grade Stainless Steel and the most advanced machines available. All CBM products comply with the Australian Standard™ AS 5369:2023, and its Normative References for compliance in infection control, instrument transport and protection.



OrthoMedics

Sterilization Basket Solutions for
both wrapped and containerised instrument sets



Contents

Standard Sterilization Baskets	Page 4
Plastic Tray Replacement Baskets	Page 6
Mini Perforated Sterilization Baskets	Page 9
Micro Perforated Sterilization Baskets	Page 10
Rigid Scope Perforated Sterilization Baskets	Page 11
Perforated Sterilization Baskets - DIN Sizes	Page 12
Robotic Surgery & Laparoscopic Frame	Page 13
Sterilization Basket Accessories	Page 15
Useful Configurations of Sterilization Baskets and Accessories	Page 24

We can help you with your requirements for your Sterilization Baskets and accessories and can assist in developing your customised trays.

Call your Territory Manager to discuss the options and configurations available or alternatively call our main office.

Toll Free: 0800 333 909

Standard Sterilization Baskets

These stainless steel baskets are the most popular sizes used by hospitals. They can be wrapped or placed in a sterilization container, the choice is yours.

They have smooth surfaces and rounded corners to minimise the risk of tears and each perforated sheet basket has the option to take 2 different styles of lids, depending on your preference.

All of these baskets can accommodate a range of inserts and accessories that can help both organise and protect your surgical instruments.


Product Code	Product Description	Dimensions L x W x H (mm)
ST25035	Perforated Sheet Basket without Lid	250 x 255 x 35
ST25050	Perforated Sheet Basket without Lid	250 x 255 x 50
ST25070	Perforated Sheet Basket without Lid	250 x 255 x 70
STL250255	Lid for Perforated Sheet Basket with Central Indent Handle	250 x 255
STL250255-R1	Lid for Perforated Sheet Basket with Side Cut-Outs for Handles	250 x 255
ST42035	Perforated Sheet Basket without Lid	420 x 255 x 35
ST42050	Perforated Sheet Basket without Lid	420 x 255 x 50
ST42070	Perforated Sheet Basket without Lid	420 x 255 x 70
STL420255	Lid for Perforated Sheet Basket with Central Indent Handle	420 x 255
STL420255-R1	Lid for Perforated Sheet Basket with Side Cut-Outs for Handles	420 x 255
ST54035	Perforated Sheet Basket without Lid	540 x 255 x 35
ST54050	Perforated Sheet Basket without Lid	540 x 255 x 50
ST54070	Perforated Sheet Basket without Lid	540 x 255 x 70
STL540255	Lid for Perforated Sheet Basket with Central Indent Handle	540 x 255
STL540255-R1	Lid for Perforated Sheet Basket with Side Cut-Outs for Handles	540 x 255
ST23025230	Inlay Tray	230 x 252 x 30
ST36025230	Inlay Tray	360 x 252 x 30



ST54070 / ST23025230



ST25050
ST42050
ST54050



STL250255 Displayed
STL420255
STL540255



STL250255-R1 Displayed
STL420255-R1
STL540255-R1



ST25050 / STL250255



ST25050 / STL250255-R1

Plastic Tray Replacement Baskets

These 3 baskets sizes are the exact dimensions of the most commonly used plastic trays. These sizes facilitate a smooth transition from a plastic basket system to a stainless steel basket system. Unlike the plastic trays, these metal baskets can be used throughout the washing and sterilization processes without the need to double handle. However, these sizes are not suitable for use with a Rigid Sterilization Container as they are not International DIN Sizes

Product Code	Product Description	Dimensions L x W x H (mm)
ST30027060	Perforated Sheet Basket with Handles without Lid	300 x 270 x 60
ST40027060	Perforated Sheet Basket with Handles without Lid	400 x 270 x 60
ST54027060	Perforated Sheet Basket with Handles without Lid	540 x 270 x 60




ST30027060

Sterilization Baskets

These particular sizes are ideal for keeping smaller and/or delicate instruments safe. These baskets can be fitted with inserts and a range of accessories, such as a silicone mat in the bottom of the basket and silicone lips on the underside of the lid, to firmly secure the instruments inside.

Product Code	Product Description	Dimensions L x W x H (mm)	
STC1705025WOL	Perforated Sheet Basket without Lid	170 x 50 x 25	
STC1705025	Perforated Sheet Basket with Lid	170 x 50 x 25	
STC1705025COM	Perforated Sheet Basket with 2x Silicone Lips on Underside of Lid & Silicone Mat	170 x 50 x 25	Most Popular
ST1705025WOLIP	Perforated Sheet Basket with Lid. 90° Corners without Handle & Upper Fold	170 x 50 x 25	
STC17013535B	Perforated Sheet Basket without Lid	170 x 135 x 35	Most Popular
STC17013535COM	Perforated Sheet Basket with 2x Silicone Lips on Underside of Lid & Silicone Mat	170 x 135 x 35	



STC1705025COM



STC17013535B



ST1705025WOLIP

Sterilization Baskets


These baskets are slightly longer and are extremely useful in securing many different types of instrument sets.

These include baskets that feature pre-designed layouts (those ending in COM), and baskets that can be customised with inserts from our accessories range.

Product Code	Product Description	Dimensions L x W x H (mm)
ST25014035	Perforated Sheet Basket without Handle	250 x 140 x 35
ST27010040	Perforated Sheet Basket without Lid	270 x 100 x 40
STL270140	Lid for Perforated Sheet Basket	270 x 140
ST27014045RC	Perforated Sheet Basket without Lid & Handle	270 x 140 x 45 <small>Most Popular</small>
ST27017035RC	Perforated Sheet Basket with Lid	270 x 170 x 35 <small>Most Popular</small>
STC27017035	Perforated Sheet Basket without Lid & Handles	270 x 170 x 35 <small>Most Popular</small>
STC27017035COM	Perforated Sheet Basket with 2x Silicone Lips on underside of Lid & Silicone Mat (Horizontal)	270 x 170 x 35 <small>Most Popular</small>
STC27017035COMLFL	Perforated Sheet Basket with 2x Silicone Lips on underside of Lid & Silicone Mat (Vertical)	270 x 170 x 35 <small>Most Popular</small>
STC27017065	Perforated Sheet Basket without Lid & Handles	270 x 170 x 65 <small>Most Popular</small>
STLC270170	Lid for Perforated Sheet Basket with middle insert handle	270 x 170
STLC270170FLAT	Lid for Perforated Sheet Basket flat with handle on top	270 x 170
STLC270170LL	Lid for Perforated Sheet Basket with 2x Silicone Lips on underside of Lid (Vertical)	270 x 170
STLC270170WOH	Lid for Perforated Sheet Basket with No Handle - Lid is Smooth and Flat	270 x 170 <small>Most Popular</small>
STLC270170FLAT	Flat lid only for all baskets	270 x 170
ST28011035RC	Perforated Sheet Basket without Lid & without Handles with Rounded Corners	280 x 110 x 35 <small>Most Popular</small>
STLC280110RC	Rounded Lid for Perforated Sheet Basket	280 x 110

STC27017035COM

STC27017035COMLFL



ST27014045RC

ST28011040RC



STC27017065


Sterilization Baskets

These baskets are specifically designed to accommodate extra long instruments. These baskets have rounded edges to facilitate wrapping, as they are too long to fit inside Rigid Sterilization Containers.


Product Code	Product Description	Dimensions L x W x H (mm)	
STC30024035	Perforated Sheet Basket with Dog Bone	300 x 240 x 35	
STLC300240-P	Flat Lid with Handle for Perforated Sheet Basket	300 x 240	
ST32017035	Perforated Sheet Basket with Lid & without Handles	320 x 170 x 35	
STC36016060	Perforated Sheet Basket without Lid	360 x 160 x 60	
ST40012045RC	Perforated Sheet Basket Round Corners without Lid & Handles	400 x 120 x 45	Most Popular
ST55031070	Perforated Sheet Basket without Lid	550 x 310 x 70	
ST56016040	Perforated Sheet Basket without Lid	560 x 160 x 40	Most Popular
STL560160	Lid for Perforated Sheet Basket	560 x 160	Most Popular
ST60046090	Perforated Sheet Basket - Wire ring around top, flat bottom without Lid	600 x 460 x 90	
ST650240180	Perforated Sheet Basket with 2x 2mm Perforated Sides without Lid	650 x 240 x 180	
STLC650240	Lid for Perforated Sheet Basket	650 x 240	
ST70025080	Perforated Sheet Basket without Lid & Handles with 2x 2mm Perforated Sides	700 x 250 x 80	
STC70030060	Perforated Sheet Basket without Lid	700 x 300 x 60	
STLC700300	Lid for Perforated Sheet Basket	700 x 300	
ST75024080	Perforated Sheet Basket with Lid	750 x 240 x 80	



ST32017035



ST70025080




ST56016040

MINI Perforated Sterilization Baskets

MINI Baskets are suitable for the protection of delicate/sharp instruments, as they are constructed from 2 x 2mm Stainless Steel Mesh on all 4 sides with the base and removable lid being constructed from 4 x 4mm Stainless Steel Mesh. MINI Baskets are compatible with many inserts from our range of accessories and can be used inside larger baskets to protect delicate instruments. They are also designed to help protect the end users and other staff during all stages of reprocessing and storage.

Product Code	Product Description	Dimensions L x W x H (mm)	
STC606020MINI	Mini Perforated Sheet Basket with Lid	60 x 60 x 20	Most Popular
STC905030MINI	Mini Perforated Sheet Basket with Lid	90 x 50 x 30	
STC1008030MINI	Mini Perforated Sheet Basket with Lid	100 x 80 x 30	
ST1008060MINI-WOLIP	Mini Perforated Sheet Basket with Lid without Lip	100 x 80 x 60	
ST15012035	Battery Tray	150 x 120 x 35	
STC1705030MINI	Mini Perforated Sheet Basket with Lid	170 x 50 x 30	
STC17010025MINI	Mini Perforated Sheet Basket with Lid	170 x 100 x 25	
ST1807035	Mini Perforated Sheet Basket with Lid	180 x 70 x 35	
ST2008035	Mini Perforated Basket without Lid & Handles Round Corners	200 x 80 x 35	
STLC20080	Lid for Perforated Sheet Basket	200 x 80	
ST2008035MINIWOL	Mini Perforated Sheet Basket without Lid	200 x 80 x 35	Most Popular
STC2008035MINI	Mini Perforated Sheet Basket with Lid	200 x 80 x 35	Most Popular
STC2008035MINISL	Mini Perforated Sheet Basket with Lid & 1 x Silicone Lip (Horizontal) on underside of Lid	200 x 80 x 35	
STC2205035MINI	Mini Perforated Sheet Basket with Lid	220 x 50 x 35	
ST24012035	Mini Perforated Sheet Basket without Lid & Handles	240 x 120 x 35	Most Popular
STL24012035	Lid for Mini Perforated Sheet Basket	240 x 120	
STC26712549MINI	Mini Perforated Sheet Basket with Lid	267 x 125 x 49	
STC28011040-BOT	Mini Perforated Sheet Basket without Lid	280 x 110 x 40	
STC28011040MINI	Mini Perforated Sheet Basket with Lid	280 x 110 x 40	

Mini-Perforation (Top & Bottom)



STC606020MINI



ST2008035MINIWOL



STC2008035MINISL



MICRO Perforated Sterilization Baskets

These MICRO Baskets are ideal for delicate and sharp instruments, as they have only 2 x 2 mm perforations on all surface areas. They can be utilised with the silicone matting to not only protect the delicate instruments inside but also to protect the users during reprocessing, transportation and storage. Many of these Micro Baskets are used inside another larger basket to keep the delicate or sharp instruments away from the other instruments.

Please note inserts cannot be fitted to micro perforated sheet baskets due to perforation size.

Product Code	Product Description	Dimensions L x W x H (mm)	
STC508030MICRO	Micro Perforated Sheet Basket with Lid	50 x 80 x 30	
STC606020MICRO	Micro Perforated Sheet Basket with Lid	60 x 60 x 20	Most Popular
STC1008030MICRO	Micro Perforated Sheet Basket with Lid	100 x 80 x 30	
STC1209050MICRO	Micro Perforated Sheet Basket with Lid	120 x 90 x 50	
STC1705025MICRO	Micro Perforated Sheet Basket with Lid	170 x 50 x 25	
STC17013535MICRO	Micro Perforated Sheet Basket with Lid	170 x 135 x 35	
ST2008035MICRO	Micro Perforated Sheet Basket with Lid	200 x 80 x 35	
ST2008035MICRO-NF	Micro Perforated Sheet Basket without Lid & Handles	200 x 80 x 35	Most Popular
ST2008035MICROWOL	Micro Perforated Sheet Basket without Lid	200 x 80 x 35	Most Popular
STC21014550MICRO	Micro Perforated Sheet Basket with Lid	210 x 145 x 50	
ST24012035MICRO-NF	Micro Perforated Sheet Basket without Lid & Handles	240 x 120 x 35	Most Popular
STC24024050MICRO	Micro Perforated Sheet Basket with Lid	240 x 240 x 50	
STC27017035MICRO	Micro Perforated Sheet Basket with Lid	270 x 170 x 35	
STC32017058MBASE	Micro Perforated Sheet Basket without Lid	320 x 170 x 58	
STC32017058MICRO	Micro Perforated Sheet Basket with Lid	320 x 170 x 58	

Micro-Perforation (all sides)




Rigid Scope Perforated Sterilization Baskets

Scopes are an expensive asset and can often be damaged beyond repair if not adequately protected during safe removal at point of use, reprocessing, storage & transportation. By using these Rigid Scope Baskets and an appropriate selection of accessories, you can protect these valuable assets and help to minimise repair and replacement costs.

These Rigid Scope Baskets are suitable for use with both Steam Sterilization and Hydrogen Peroxide (high grade disinfection) processes.

Product Code	Product Description	Dimensions L x W x H (mm)
STC2607045RSBWOL	Rigid Scope Basket without Lid	260 x 70 x 45
STC2607045RSB	Rigid Scope Basket with Lid	260 x 70 x 45 <small>Most Popular</small>
STC2907060RSB	Rigid Scope Basket with Lid	290 x 70 x 60
ST2907060RSB-COM	Rigid Scope Basket with Lid & with Inserts for Smaller Scopes 3-5mm Diameter	290 x 70 x 60
STC35015045RSB	Rigid Scope Basket Double with Lid	350 x 150 x 45
STC4007045RSB	Rigid Scope Basket with Lid	400 x 70 x 45 <small>Most Popular</small>
STC40012045RSBWOL	Rigid Scope Basket Double without Lid	400 x 120 x 45
STC40012045RSB	Rigid Scope Basket Double with Lid	400 x 120 x 45 <small>Most Popular</small>
STC4507060RSB	Rigid Scope Basket with Lid	450 x 70 x 60 <small>Most Popular</small>
STC48015050RSB	Rigid Scope Basket Double with Flat Lid	480 x 150 x 50
STC48015080RSB	Rigid Scope Basket Double with Flat Lid	480 x 150 x 80
STC58015080RSB	Rigid Scope Basket Double with Lid	580 x 150 x 80
STC6007060RSB	Rigid Scope Basket with Lid	600 x 70 x 60 <small>Most Popular</small>
STC60015070RSB	Rigid Scope Basket Double with Lid	600 x 150 x 70
STC65015070RSB	Rigid Scope Basket Double with Lid	650 x 150 x 70
STC75015080RSB	Rigid Scope Basket Ideal for Robotic Surgical Instruments with Lid	750 x 150 x 80



STC2907060RSB


STC4007045RSB

Perforated Sterilization Baskets - DIN Sizes

Baskets recommended under DIN EN 868-8 and DIN 58952-3 German Standards Organisation.

Product Code	Product Description	Dimensions L x W x H (mm)
1 DIN		
STC480250150	Perforated Sheet Basket	480 X 250 X 150
STC480250100	Perforated Sheet Baskets	480 X 250 X 100
STC48025080	Perforated Sheet Basket	480 X 250 X 80
STC48025060	Perforated Sheet Basket	480 X 250 X 60
ST48025050	Perforated Sheet Basket	480 X 250 X 50
STLC480250	Lid for Perforated Sheet Basket	480 X 250
STLC480250-R1	Lid with cut outs for Perforated Sheet Basket	480 X 250
1/2 DIN		
STC25024060	Perforated Sheet Basket	240 X 250 X 60
STC240250100	Perforated Sheet Basket	240 X 250 X 100
STLC250240	Lid for Perforated Sheet Basket	250 X 240
1/4 DIN		
STC250120140	Perforated Sheet Basket	120 X 250 X 140
STC12025080	Perforated Sheet Basket	120 X 250 X 80
STC25012060	Perforated Sheet Basket	120 X 250 X 60
STC25012040	Perforated Sheet Basket	120 X 250 X 40
STLC120250	Lid for Perforated Sheet Basket	120 X 250

STC480250100
STLC480250





STC280250100
STLC480250-R1

Robotic Surgery

These baskets have been specifically designed to safely store and handle complex robotic surgical instrumentation. They have additional features to enhance their rigidity and strength, on both the base and the lid.

Product Code	Product Description	Dimensions L x W x H (mm)
STC683354101	Perforated basket with Double Closure Lid	683 X 354 X 101
STC683354101-XL	Perforated basket with Double Closure Lid & Inserts	683 X 354 X 101




Laparoscopic Frame

Laparoscopic frames are very useful for the creation of a user friendly Laparoscopic Instrument Set.

The removable Laparoscopic Frame is designed to be aseptically lifted out of either a base metal perforated basket and/or a rigid sterilization container and placed safely on the instrument trolley.

Product Code	Product Description	Dimensions L x W x H (mm)
STAL508221142-2A	Laparoscopic Frame 10 places	508 x 221 x 142mm
STAL508221142-14	Laparoscopic Frame 14 places	508 x 221 x 142mm



Customise your trays with Inserts & Accessories

These silicone inserts and accessories will help to both organise your instrument sets and help protect small and delicate instruments from damage. There are many different options available.

Your local Territory Manager can help you in the initial design and set-up of your instrument sets.



ST25050
STL250255
STAS23203/08SS
S23203/08XS
STAS232FS



STC27017035
STLC270170FLAT
STAS23203/08SS

Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAFB020	Fixing Bolts for Ring Handled Instruments	20mm high
STAFB032	Fixing Bolts for Ring Handled Instruments	32mm high
STAFB050	Fixing Bolts for Ring Handled Instruments	50mm high
STASB666	Stainless Steel Screws TOR x20	Pk/10
STASB667-1	M4 Fixing Screw - also available in packs of 10	1 each
STASB666P	Stainless Steel Flat Head Style Screw	Pk/10
STASSL1016S	Replacement Screw Head for Silicone Lid Lip Inserts	1 each
STASSLL90	Moveable 'Soft Lip' Silicone Bar for underside of Lid	90 x 16
STASSLL130	Moveable 'Soft Lip' Silicone Bar for underside of Lid	130 x 16
STASSLL220	Moveable 'Soft Lip' Silicone Bar for underside of Lid	220 x 16
STASSLL350	Moveable 'Soft Lip' Silicone Bar for underside of Lid	350 x 16
STF04	Set of 4 Plastic Feet	
STPTI-A4	Rigid Sterilization Container Document Holder for Set	150 x 80



Fixing Bolts
STAFB020
STAFB032
STAFB050



SS TORX Head Screw
STASB666




M4 Fixing Screw
STASB667-1



SS Flat Head Screw
STASB666P



Replacement Screw
Head for Lid Lip
STASSL1016S



Movable 'Soft Lip' Silicone Bar
for underside of Lid
STASSLL130








Set of 4 Plastic Feet
STF04



Document Holder
STPTI-A4


Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAS72FS	Comb Shape Instrument Holder 5 places	72 x 20
STAS100FS	Comb Shape Instrument Holder 7 places	100 x 20 <small>Most Popular</small>
STAS165FS	Comb Shape Instrument Holder 13 places	165 x 20
STAS178FS	Comb Shape Instrument Holder 14 places	178 x 20 <small>Most Popular</small>
STAS232FS	Comb Shape Instrument Holder 18 places	232 x 20 <small>Most Popular</small>




Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAS03/08	Star Shape Instrument Holder Open 1 places: 3 - 8mm diameter	38 x 20
STASD03/08	Star Shape Instrument Holder Open 2 places: 3 - 8mm diameter	38 x 20 Most Popular
STAS7203/08H	Star Shape Instrument Holder Open 3 places: 3 - 8mm diameter	72 x 20
STAS7208/15H	Star Shape Instrument Holder Open 3 places: 8 - 15mm diameter	72 x 30
STAS7203/08B	Star Shape Instrument Holder Open 4 places: 3 - 8mm diameter	72 x 20
STAS8903/08SS	Star Shape Instrument Holder Open 5 places: 3 - 8mm diameter	89 x 20 Most Popular
STAS7203/08SS	Star Shape Instrument Holder Open 6 places: 3 - 8mm diameter	72 x 20 Most Popular
STAS10003/08SS	Star Shape Instrument Holder Open 7 places: 3 - 8mm diameter	100 x 20 Most Popular
STAS17808/15H	Star Shape Instrument Holder Open 8 places: 8 - 15mm diameter	178 x 30
STAS17803/08SS	Star Shape Instrument Holder Open 10 places: 3 - 8mm diameter	178 x 20 Most Popular
STASH23203/08SS	Star Shape Instrument Holder Open 13 places: 3 - 8mm diameter	232 x 50
STAS23203/08SS	Star Shape Instrument Holder Open 18 places: 3 - 8mm diameter	232 x 20 Most Popular




Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAS7203S	Key Shape Instrument Holder 5 places 3mm diameter	72 x 20
STAS7205S	Key Shape Instrument Holder 5 places 5mm diameter	72 x 20 <small>Most Popular</small>
STAS7206S	Key Shape Instrument Holder XL 5 places 6mm diameter	72 x 20
STAS7203SB	Key Shape Instrument Holder 6 places 3mm diameter	72 x 20
STAS17803S	Key Shape Small Instrument Holder 12 places 3mm diameter	178 x 20
STAS17805S	Key Shape Large Instrument Holder 12 places 5mm diameter	178 x 20




Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAS11603/09XS	Christmas Tree Shape Instrument Holder 5 place 3 - 9mm diameter	116 x 30
STAS10003/08XS	Christmas Tree Shape Instrument Holder 7 place 3 - 5mm diameter	100 x 20
STAS23203/08XS	Christmas Tree Shape Instrument Holder 18 place 3 - 5mm diameter	232 x 20 <small>Most Popular</small>
STAS23203/09XS	Christmas Tree Shape Instrument Holder 11 place 3 - 9mm diameter	232 x 30









Product Code	Product Description	Dimensions L x H (mm)
STAS38CC	Silicone Separator Solid Insert	38 x 30
STAS72CC	Silicone Separator Solid Insert	72 x 30
STAS7243HRSB	Slotted Silicone Holder Small Hole	72 x 43
STAS7231HRSB	Slotted Silicone Holder Large Hole	72 x 31







Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAS3812K	Endoscope Holder Open End 12mm diameter	38 x 30
STAS3815K	Endoscope Holder Open End 15mm diameter	38 x 30
STAS3812RSB	Endoscope Holder Open End 15mm diameter	38 x 45
STAS3819K	Endoscope Holder Open End 19mm diameter	38 x 30 <small>Most Popular</small>
STAS5020RSB	Endoscope Holder Open End 20mm diameter	50 x 45
STAS5022RSB	Endoscope Holder Open End 22mm diameter	50 x 41 <small>Most Popular</small>
STAS3825RSB	Endoscope Holder Open End 25mm diameter	38 x 39
STAS5027RSB	Endoscope Holder Open End 27mm diameter	50 x 41 <small>Most Popular</small>
STAS7231LRB	Endoscope Holder Open End 31mm diameter	72 x 45 <small>Most Popular</small>
STAS7243LRB	Endoscope Holder Open End 43mm diameter	72 x 45



Product Code	Product Description	Dimensions L x H (mm)
STAS3804C	Endoscope Holder Closed End 4mm diameter	38 x 30
STAS3806C	Endoscope Holder Closed End 6mm diameter	38 x 30
STAS3808C	Endoscope Holder Closed End 8mm diameter	38 x 30 <small>Most Popular</small>
STAS3812C	Endoscope Holder Closed End 12mm diameter	38 x 30 <small>Most Popular</small>





Inserts & Accessories





Product Code	Product Description	Dimensions L x H (mm)	
STAS3803/08H	Closed Star Single Instrument Holder 3 - 8mm diameter	38 x 25	Most Popular
STAS3808/15H	Closed Star Single Instrument Holder 8 - 15mm diameter	38 x 30	Most Popular



Product Code	Product Description	Dimensions L x H (mm)
STASD3803/08H	Small Star Shape Single Instrument Holder 3 - 8mm 2 places	72 x 30
STASD120308	Star Shape Single Instrument Holder 3 - 8mm with Keyhole 12mm diameter	72 x 30
STASD150308	Star Shape Single Instrument Holder 3 - 8mm with Keyhole 15mm diameter	72 x 30
STASD120815	Star Shape Single Instrument Holder 8 - 15mm with Keyhole 12mm diameter	72 x 30
STASD150815	Star Shape Single Instrument Holder 8 - 15mm with Keyhole 15mm diameter	72 x 30




Product Code	Product Description	Dimensions L x H (mm)	
STASH3812K	Endoscope Holder High Open End 12mm diameter	38 x 56	Most Popular
STASH3815K	Endoscope Holder High Open End 15mm diameter	38 x 56	Most Popular
STASH3803/08H	Endoscope Holder High Hinged 3 - 8mm diameter	38 x 50	Most Popular
STASH3808/15H	Endoscope Holder High Hinged 8 - 15mm diameter	38 x 54	Most Popular
STAS3802V	Endoscope Holder V Shape 2mm diameter	38 x 30	



Inserts & Accessories

Product Code	Product Description	Dimensions L x H (mm)
STAD5030-PERF	Stainless Steel Tray Divider	50 x 30
STAD5035-PERF	Stainless Steel Tray Divider	50 x 35
STAD12035-PERF	Stainless Steel Tray Divider	120 x 35
STAD20035-PERF	Stainless Steel Tray Divider	200 x 35
STAD45035-PERF	Stainless Steel Tray Divider	450 x 35




Stainless Steel Tray Divider
STAD12035-PERF
Includes fixing bolts

Product Code	Product Description	Dimensions L x H (mm)
STAT6035B-M	Stainless Steel Tag (Blank) Modified with Shorter Neck	60 x 35
STAT6030B-M	Stainless Steel Tag (Blank) Modified with Shorter Neck	60 x 30
STAT6025-CB	Stainless Steel Pronged Clip on Label (Blank)	60 x 25
STAT9025-CB	Stainless Steel Pronged Clip on Label (Blank)	90 x 25




Stainless Steel Tag
STAT6035B-M



STAT6025-CB

Silicone Mats

Product Code	Product Description	Dimensions L x H (mm)
STASM165130MW	Silicone Mat Micro Weave	165 x 130
STASM190165MW	Silicone Mat Micro Weave	190 x 165
STASM195135MW	Silicone Mat Micro Weave	195 x 135
STASM220160MW	Silicone Mat Micro Weave	220 x 160
STASM220220MW	Silicone Mat Micro Weave	220 x 220
STASM260160MW	Silicone Mat Micro Weave	260 x 160
STASM260220MW	Silicone Mat Micro Weave	260 x 220
STASM280120MW	Silicone Mat Micro Weave	280 x 120
STASM280220MW	Silicone Mat Micro Weave	280 x 220
STASM320220MW	Silicone Mat Micro Weave	320 x 220 Most Popular
STASM360220MW	Silicone Mat Micro Weave	360 x 220
STASM400220MW	Silicone Mat Micro Weave	400 x 220
STASM460140MW	Silicone Mat Micro Weave	460 x 140
STASM460220MW	Silicone Mat Micro Weave	460 x 220
STASM500210MW	Silicone Mat Micro Weave	500 x 210
STASM520220MW	Silicone Mat Micro Weave	520 x 220 Most Popular




Soft silicone with multiple perforations that allows effective cleaning and steam penetration for sterilization.

Designed to protect your instruments inside a sterilization basket throughout both the wash cycle and the sterilization cycle.

They are available in many different sizes to meet your needs but can easily be cut to the exact size to meet your requirements.


Examples of Instrument Set Layouts

Small Joint Arthroscopy Set 1




Product Code	Product Description	Quantity
ST54050	Perforated Sheet Basket without Lid	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	1
STAS7205S	Key Shape Instrument Holder 5 places 5mm diameter	1
STAS7203/08H	Star Shape Small Instrument Holder 3 places: 3 - 8mm diameter	1
STAS3819K	Endoscope Holder Open End 19mm diameter	1
STAS3806C	Endoscope Holder Closed End 6mm diameter	1
STAFB032	Fixing Bolts for Ring Handled Instruments	5

Small Joint Arthroscopy Set 2




Product Code	Product Description	Quantity
ST25050	Perforated Sheet Basket without Lid	1
STAS72FS	Comb Shape Instrument Holder 5 places	1
STAS7203/08H	Star Shape Small Instrument Holder 3 places: 3 - 8mm diameter	1
STAS3819K	Endoscope Holder Open End 19mm diameter	2
STAS3806C	Endoscope Holder Closed End 6mm diameter	2
STAFB032	Fixing Bolts for Ring Handled Instruments	3

Resectoscope Instrument Set 1




Product Code	Product Description	Quantity
ST40012045RC	Perforated Sheet Basket Round Corners without Lid and Handles	1
STAS8903/08SS	Star Shape Instrument Holder Small 5 places	2

Resectoscope Instrument Set 2




Product Code	Product Description	Quantity
ST40012045RSB	Perforated Sheet Basket Round Corners without Lid and Handles	1
STAS8903/08SS	Star Shape Instrument Holder Small 5 places	2

Plastics Instrument Set 1




Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAS232FS	Comb Shape Instrument Holder 18 places	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	1
STAS7231LR SB	Endoscope Holder Open End 31mm diameter	1
STAS100FS	Comb Shape Instrument Holder 7 places	1
STAFB032	Fixing Bolts for Ring Handled Instruments	6

Plastics Instrument Set 2




Product Code	Product Description	Quantity
ST25050	Perforated Sheet Basket without Lid	1
STAFB032	Fixing Bolts for Ring Handled Instruments	9
STAS100FS	Comb Shape Instrument Holder 7 places	2

Micro Plastics Instrument Set



Product Code	Product Description	Quantity
ST42035	Perforated Sheet Basket without Lid	1
STAS232FS	Comb Shape Instrument Holder 18 places	3
STAS178FS	Comb Shape Instrument Holder 14 places	2
STAS100FS	Comb Shape Instrument Holder 7 places	1
STC606020MINI	Mini Perforated Sheet Basket with Lid	1
STASB667-1	Fixing Bolt with Rounded Head to secure Mini Basket w/lid	1

D & C Instrument Set



Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAFB032	Fixing Bolts for Ring Handled Instruments	9
STAS8903/08SS	Star Shape Instrument Holder Small 5 places	1

General Laparoscopic Instrument Set




Product Code	Product Description	Quantity
ST54050	Perforated Sheet Basket without Lid	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	2
STAD20035-PERF	Stainless Steel Tray Divider	1

Minor General Instrument Set




Product Code	Product Description	Quantity
ST54070	Perforated Sheet Basket without Lid	1
STAFB050	Fixing Bolts for Ring Handled Instruments	15
STAS165FS	Comb Shape Instrument Holder 13 places	1

Phaco Extras Instrument Set




Product Code	Product Description	Quantity
ST27017035RC	Perforated Sheet Basket without Lid	1
STAS100FS	Comb Shape Instrument Holder 7 places	1
STASD150815	Star Shap Single Instrument Holder 8-15mm with Keyhole 15mm diameter	1
STAFB020	Fixing Bolts for Ring Handled Instruments	3

Cataract Instrument Set




Product Code	Product Description	Quantity
ST25050	Perforated Sheet Basket without Lid	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	1
STAS232FS	Comb Shape Instrument Holder 18 places	1
STAS7203/08SS	Star Shape Small Instrument Holder 6 places: 3 - 8mm diameter	2
STAS7206S	Key Shaped Instrument Holder XL 5 places 6mm diameter	1
STAFB032	Fixing Bolts for Ring Handled Instruments	2

Myringoplasty Instrument Set




Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAS178FS	Comb Shape Instrument Holder 14 places	1
STAS232FS	Comb Shape Instrument Holder 18 places	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	1
STAS3812K	Endoscope Holder Open End 12mm Diameter	1
STAFB032	Fixing Bolts for Ring Handled Instruments	14
STASSLL350	Movable Soft Lip Silicone Bar for underside of Lid	1

Micro Ear Instrument Set




Product Code	Product Description	Quantity
ST42035	Perforated Sheet Basket without Lid	1
STL420255-R1	Lid with Side Cut-Outs for Handles	1
STC1705025	Perforated Sheet Basket with Lid	1
STAS232FS	Comb Shape Instrument Holder 18 places	1
STAS72FS	Comb Shape Instrument Holder 6 places	1
STAS23203/08SS	Star Shape Small Instrument Holder 18 places: 3 - 8mm diameter	1
STAS7203/08SS	Star Shape Small Instrument Holder 6 places: 3 - 8mm diameter	1
STASSLL350	Soft Lip Silicone Bar for Lid	2
STAFB667-1	M4 Fixing Screw + Nut	2
STAFB032	Fixing Bolts for Ring Handled Instruments	6

FESS Instrument Set




Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAS17803/08SS	Star Shape Small Instrument Holder 10 places: 3 - 8mm diameter	1
STAS178FS	Comb Shape Instrument Holder 14 places	1
STAFB032	Fixing Bolts for Ring Handled Instruments	12

Trans Canal Instrument Set




Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAS232FS	Comb Shape Instrument Holder 18 places	1
STAS178FS	Comb Shape Instrument Holder 14 places	1
STAFB032	Fixing Bolts for Ring Handled Instruments	21

Probe Tray Set




Product Code	Product Description	Quantity
ST42070	Perforated Sheet Basket without Lid	1
STL420255-R1	Lid with Side Cut-Outs for Handles	1
STAS3812K	Endoscope Holder Open End 12mm diameter	1
STAS3803/08H	Closed Star Single Instrument Holder 3 - 8mm diameter	1
STAFB050	Fixing Bolts for Ring Handled Instruments	3

Cystoscopy Instrument Set




Product Code	Product Description	Quantity
ST42050	Mini Perforated Sheet Basket with Lid	1
STAS7205S	Key Shape Instrument Holder 5 places 5mm diameter	4

Olympus Videoscope Tray




Product Code	Product Description	Quantity
ST54070	Perforated Sheet Basket without Lid	1
STL540255-R1	Lid with Side Cut-Outs for Handles	1
STAS3825RSB	Endoscope Holder Open End 25mm diameter	1
STAS3815K	Endoscope Holder Open End 15mm diameter	2
STAS3808C	Endoscope Holder Closed End 8mm diameter	1
STAS3804C	Endoscope Holder Closed End 4mm diameter	1
STAS7231HRSB	Slotted Silicone Holder Large Hole	1

Olympus Fibre Scope, Light Lead & Accessories Tray Set




Product Code	Product Description	Quantity
ST54070	Perforated Sheet Basket without Lid	1
STASM165130MW	Silicone Mat Micro Weave	1
STAS3808/15H	Closed Star Singe Instrument Holder 8 - 15mm diameter	1
STAS3819K	Endoscope Holder Open End 15mm diameter	1
STAS5027RSB	Endoscope Holder Open End 27mm diameter	1

Single Instruments




Product Code	Product Description	Quantity
STC1705030MINI	Mini Perforated Sheet Basket with Lid	1
STAS3803/08H	Closed Star Single Instrument Holder 3-8mm diameter	1

Dental Implant Instrument Set




Product Code	Product Description	Quantity
STC1008030MINI	Mini Perforated Sheet Basket with Lid	1
STAFB020	Fixing Bolts for Ring Handled Instruments	3

Neuro Micro Instrument Set



Product Code	Product Description	Quantity
ST42050	Perforated Sheet Basket without Lid	1
STAS232FS	Comb Shape Instrument Holder 18 places	2
STAS17803/08SS	Star Shape Small Instrument Holder 10 places: 3 - 8mm diameter	2

Ziess Lens Tray Set



	Product Description	Quantity
STC25012060	Perforated Sheet Basket DIN	1
STLC120250	Lid for Basket	1
STAS72CC	Silicone Insert Separator Solid	1
STASD03/08	Star Shape Instrument Holder 2 places	1
STAT6025-CB	Stainless Steel Clip on Label	1

Hard to Find Instruments?

Let us do the looking for you



Ortho Medics
586 Great South Road
Ellerslie Auckland 1051
New Zealand
admin@orthomedics.co.nz
TOLL FREE: 0800 333 909
GST No. 19 079 279
2024

