


Container Seal - 2,7 cm long pin



Secures your sterile container...

In order to secure the container for opening by un-authorized staff the container seals are offering an important shackle in the quality management chain. Prevents opening by un-authorized staff. Only by breaking the seal the container can be opened. A check at the OR for broken seals is guaranteeing the sterility of materials from the moment of sterilization till the actual use.

REF: 902.111.1000	Container seal, blue 2,7 cm	Unit: 1000 pieces
REF: 902.112.1000	Container seal, blue 2,7 cm with steam ind.	Unit: 1000 pieces
REF: 902.113.1000	Container seal, red 2,7 cm	Unit: 1000 pieces
REF: 902.114.1000	Container seal, red 2,7 cm with steam ind.	Unit: 1000 pieces
REF: 902.115.1000	Container seal, orange 2,7 cm	Unit: 1000 pieces
REF: 902.116.1000	Container seal, orange 2,7 cm with steam ind.	Unit: 1000 pieces
REF: 902.117.1000	Container seal, green 2,7 cm	Unit: 1000 pieces
REF: 902.118.1000	Container seal, green 2,7 cm with steam ind.	Unit: 1000 pieces

Characteristics:

- ◆ Produced in polypropylene , softening above 150°C.
- ◆ Suitable for steam sterilization processes.
- ◆ Fulfilling EN 868 part 8 section 4.2.3
- ◆ Steam indicator conform to ISO 11140 - Type 1
- ◆ Clear color transition of the steam indicator



902.111.1000



902.112.1000



902.113.1000



902.114.1000



902.115.1000



902.116.1000