

## Plasma Tape - 19 mm x 60 meter



### Sticky, Indicating with reference color ....

Especially selected for the Plasma sterilization processes the Sterintech™ Plasma tapes provide a good stickiness to most surfaces. Now with a reference color as a post-sterilization reference

REF: 504.019.0001

Plasma tape 19 mm x 60 meter

Unit: 1 piece

Providing a quick recognition of sterile goods by it's indicator while securing the sterility of wrapped items our plasma tape performs in Hospitals anywhere in the world.

#### Characteristics:

- ◆ Single packed to guarantee shelf-life specifications up to 2 years after production
- ◆ Easy unwinding from the roll and easy tearing by the side
- ◆ Plasma indicator conform to ISO 11140 - Type 1
- ◆ Clear red - yellow color transition of the indicator
- ◆ Latex free glue
- ◆ Non-Toxic / Heavy metal free

Other types and dimensions available upon request.