


Steam Tapes



Whatever purpose they stick

Especially selected for the various steam sterilization processes the Steritech™ steam tapes provide a good stickiness to most surfaces.

REF: 502.019.0048

Steam tape 19 mm x 50 mtr

Unit: 48 pieces

REF: 502.025.0036

Steam tape 25 mm x 50 mtr

Unit: 36 pieces

Providing a quick recognition of sterile goods by it's indicator (yellow-dark) while securing the sterility of wrapped items our steam tape performs in Hospitals anywhere in the world.

Characteristics:

- ◆ Single packed to guarantee shelf-life specifications up to 3 years after production
- ◆ Easy unwinding from the roll
- ◆ Steam indicator conform to ISO 11140 - Type 1
- ◆ Clear yellow - dark color transition of the indicator
- ◆ Latex free glue
- ◆ Non-Toxic / Heavy metal free
- ◆ Adhesive works well on Soft Crepe and Non-Woven materials
- ◆ Suitable for steam cycles 121 and 134 Degrees Celsius

