

Sterilization Rolls - 60 gsm - green film



Sterile barrier ... be sure about it.

Today Infection control is getting increasingly important since more and more bacterias starts to get resistant to the latest anti-biotics. Infection Prevention is the key word. Sterilization packaging materials have been developed over many years to increase the stability of the sterile barrier over time.

Sterintech™ Sterilization rolls are used to wrap medical devices. The flexibility of length cutting is making a dedicated packaging of individual instruments possible and offer a number of advantages to the end-user

The Sterintech™ Sterilization rolls are characterized by the following:


- ◆ Easy recognition of content by the multi-layer transparent film
- ◆ Strong , multi-line seals
- ◆ Indicators and text printed outside the packing area
- ◆ Non-tearing film avoids tearing at opening the pack
- ◆ Clean and fibre-free opening of pack
- ◆ Compatible to STEAM, EO and FORM sterilization by three imprinted indicators
- ◆ Compliant to standard ISO 11607-1, EN 868 - 5 and ISO 11140-1

Please refer to the next side for all standard dimensions.



Sterilization Rolls

Packaging



Sterile barrier ... be sure about it.

REF: 301.050.0001	Flat reel, 60 gr/sqm, 50 mm x 200 mtr	Unit: 1 pcs
REF: 301.075.0001	Flat reel, 60 gr/sqm, 75 mm x 200 mtr	Unit: 1 pcs
REF: 301.100.0001	Flat reel, 60 gr/sqm, 100 mm x 200 mtr	Unit: 1 pcs
REF: 301.125.0001	Flat reel, 60 gr/sqm, 125 mm x 200 mtr	Unit: 1 pcs
REF: 301.150.0001	Flat reel, 60 gr/sqm, 150 mm x 200 mtr	Unit: 1 pcs
REF: 301.200.0001	Flat reel, 60 gr/sqm, 200 mm x 200 mtr	Unit: 1 pcs
REF: 301.250.0001	Flat reel, 60 gr/sqm, 250 mm x 200 mtr	Unit: 1 pcs
REF: 301.300.0001	Flat reel, 60 gr/sqm, 300 mm x 200 mtr	Unit: 1 pcs
REF: 301.350.0001	Flat reel, 60 gr/sqm, 350 mm x 200 mtr	Unit: 1 pcs
REF: 301.400.0001	Flat reel, 60 gr/sqm, 400 mm x 200 mtr	Unit: 1 pcs

REF: 302.075.0001	Gusseted reel, 60 gr/sqm, 75 mm x 100 mtr	Unit: 1 pcs
REF: 302.100.0001	Gusseted reel, 60 gr/sqm, 100 mm x 100 mtr	Unit: 1 pcs
REF: 302.150.0001	Gusseted reel, 60 gr/sqm, 150 mm x 100 mtr	Unit: 1 pcs
REF: 302.200.0001	Gusseted reel, 60 gr/sqm, 200 mm x 100 mtr	Unit: 1 pcs
REF: 302.250.0001	Gusseted reel, 60 gr/sqm, 250 mm x 100 mtr	Unit: 1 pcs
REF: 302.300.0001	Gusseted reel, 60 gr/sqm, 300 mm x 100 mtr	Unit: 1 pcs
REF: 302.350.0001	Gusseted reel, 60 gr/sqm, 350 mm x 100 mtr	Unit: 1 pcs
REF: 302.400.0001	Gusseted reel, 60 gr/sqm, 400 mm x 100 mtr	Unit: 1 pcs