


Sterilization Rolls (Tyvek 1059B)



Sterile barrier ... be sure about it.

Today Infection control is getting increasingly important since more and more bacterias starts to get resistant to the latest anti-biotics. Infection Prevention is the key word. Sterilization packaging materials have been developed over many years to increase the stability of the sterile barrier over time.

Sterintech™ Sterilization rolls (Tyvek) are used to pack medical devices for Plasma Sterilization. The flexibility of length cutting is making a dedicated packaging of individual instruments possible and offer a number of advantages to the end-user


The Sterintech™ Sterilization rolls (Tyvek) are characterized by the following:

- ◆ Easy recognition of content by the multi-layer transparent film
- ◆ Strong seal
- ◆ Indicators and text printed outside the packing area
- ◆ Non-tearing film avoids tearing at opening the pack
- ◆ Clean and fibre-free opening of pack
- ◆ Compatible to Plasma sterilization with imprinted indicators
- ◆ Compliant to standard ISO 11607-1, EN 868 - 5 and ISO 11140-1

Please refer to the next side for all standard dimensions.



Sterilization Rolls (Tyvek)



Sterile barrier ... be sure about it.

Standard dimensions. Bespoke products available upon request.

REF: 303.075.0001	Flat reel, 1059B,	75 mm x 70 mtr	Unit: 1 pcs
REF: 303.100.0001	Flat reel, 1059B,	100 mm x 70 mtr	Unit: 1 pcs
REF: 303.150.0001	Flat reel, 1059B,	150 mm x 70 mtr	Unit: 1 pcs
REF: 303.200.0001	Flat reel, 1059B,	200 mm x 70 mtr	Unit: 1 pcs
REF: 303.250.0001	Flat reel, 1059B,	250 mm x 70 mtr	Unit: 1 pcs
REF: 303.300.0001	Flat reel, 1059B,	300 mm x 70 mtr	Unit: 1 pcs
REF: 303.350.0001	Flat reel, 1059B,	350 mm x 70 mtr	Unit: 1 pcs
REF: 303.400.0001	Flat reel, 1059B,	400 mm x 70 mtr	Unit: 1 pcs
REF: 303.500.0001	Flat reel, 1059B,	500 mm x 70 mtr	Unit: 1 pcs

Material: Tyvek 1059B (64 gr/sqm) and 52um PET/PPP film