


Plasma Indicator - Type 4

SAP / Monitoring



when you want to check inside....

These indicator strips can be used as in-pack indicator in order to detect sterilization conditions have been met inside the pack.

REF: 104.022.0200


Plasma Indicator

Unit: 200 pieces

The Sterintech™ Plasma Indicator has been tested by means of numerous cycles of Plasma sterilization in order to secure a proper working when used in various Plasma sterilizers.

Characteristics:

- ◆ Shelf life : 5 years after production
- ◆ Clear color change from red to yellow
- ◆ Conforms to ISO 11140 - 1:2014 Type 4
- ◆ Very easy interpretation of result
- ◆ Based upon a PP Material
- ◆ Non-Toxic



Dimensions: 75 x 20 mm (per strip)

