


Emulator - Type 6



The Emulator with a big plus...

Emulators can be used in addition to the batch control Helix which will allow you to release load. Especially when OR staff wants to be re-affirmed about the sterilization conditions.

REF: 106.103.0500 Emulator 134°C-3,5 min/121°C-15 min Unit: 500 pieces
REF: 106.103.0500SA Emulator 134°C-3,5 min/121°C-15 min self-adhesive Unit: 500 pieces

The Sterintech™ Emulators are offering you more information about the steam sterilization process by means of the very accurate indicator fields of the Emulator. These Emulators can detect many disturbances of the sterilization parameters inside packs or containers.

Characteristics:

- ◆ Steam emulator conform to ISO 11140 - Type 6
- ◆ Shelf life : 5 years after production
- ◆ Clear color transition from Blue into Pink
- ◆ Chemical emulator covered with lacquer
- ◆ Non-Toxic / Lead Free