

a1 medical



O P H T H A L M O L O G Y

COMPANY

a1 medical is a manufacturer of highest quality micro-surgical instruments used in Ophthalmology, offering customers the freedom of choice. These hallmarks have proven to address a gap in the market, since incorporation of the company in 2012. Consequently more and more surgeons and nurses worldwide experience and value the consistent high quality level of a1 medical products and services.

The a1 medical team combines decades of experience in medical technology, with passion for instruments, skills, team spirit and mastery of latest manufacturing technologies. We deliver on time, are responsive and committed to medical progress. Driving improved approaches for medical routines by realizing novel instrument designs together with our partners and customers is therefore our continued effort.

We are independent, rooted in the Tuttlingen region, manufacture and distribute highest quality level instruments Made in Germany, enabling medical professionals to help patients as best as humanly possible.

a1 medical provides you with the reliability patients deserve.

NOTE

This catalogue has been made for you with great care. However, we would like to point out that we assume no liability for the correctness of every detail in the catalogue.

The manufacturer reserves the right to change technical details. No part of this catalogue may be reproduced (printed, photocopied or any other means) or used without the written permission of a1 medical GmbH.

Photos may vary slightly from actual size.

We would like to thank you in advance for your understanding.

TABLE OF CONTENTS

	Numeric Prefix	CHAPTER	Page
1.		Ceramic Coating Information	6-7
2.	ES-	Chisels, Eye Speculums, Fixation Hooks, Fixation Rings, Lid Plates, Muscle Hooks, Needles, Retractors, Scleral Depressors, Strabismus Hooks, Instruments for MIG-Surgery	9-34
3.	C-	Blade Handles, Burrs, Dissectors, Diamond Knives, Foreign Body Instruments, IOL Cutter, Knives, Punches, Spuds, Trephines, Needles, Corneal Scarifiers	35-50
4.	HS-	Fixation Rings, Hooks, Keratometer, Loops, Marker, Manipulators, Marking Pen, Phaco Choppers, Preassure Plate, Spatulas & Curettes, Spoon, Titanium Hooks, a1 medical Single- and Double Hooks	51-88
5.	F-	Forceps, MIC Forceps, Titanium Forceps, Bipolar Forceps, ErgoGrip Forceps, TC Forceps, SIGNUM MICS Forceps, SMILE Grip Forceps	89-154
6.	N-	Blade Breakers, Needle Holders & Titanium Needle Holders, ErgoGrip Needle Holders, TC Needle Holders	155-170
7.	S-	Micro Scissors, Super Cut Scissors, Ring Scissors	171-190
8.	VR-	Vitrectomy Instruments	191-196
9.	L-	Bimanual I/A Instruments, Cannulas, Lacrimal Instruments, Detachable I/A Instruments, Semi-Detachable Punches	197-238
10.	M-	Marking Instruments	239-242
11.	MIS-	Boles, Eye Shields, Miscellaneous, Trays & Cassettes	243-257
12.		Content	258-267
13.		Alphabetical Index	268-275

M A D E I N G E R M A N Y

a1 medical

CERAMIC COATING

- ... LONGER PRODUCT LIFE
- ... EASIER CLEANING
- ... REDUCED EYE-FATIGUE
- ... 100 % BIOCOMPATIBLE
- ... ANTI ALLERGIC



CERAMIC COATING

ORDER INFORMATION

Please add „C“ to the end of the article number. e.g. F-0010C

... improves instruments:

- cost reduction and longer product life through
 - higher surface hardness** – up to 10 times harder than stainless steel
 - lower friction** – improves wear and tear of scissor blades
 - more resistance** against corrosion
- reduced eye-fatigue through less glare due to light absorbing black surface
- easier cleaning and reprocessing due to lower surface energy and non-stick effect
- anti-allergic due to absence of Chromium and Nickel in surface layer

... established in cardiology and neurology

... consists of TiAlN, applied in a magnetron sputtering process

... is 100% quality controlled and comes with lifetime warranty

EYE SPECULUMS

RETRACTORS

FIXATION RINGS

FIXATION HOOKS

MUSCLE HOOKS

CHISELS

LID PLATES

SCLERAL DEPRESSORS

FASCIA NEEDLE

STRABISMUS HOOKS

INSTRUMENTS FOR

MIG-SURGERY

FOR MIG – SURGERY

Specially designed products for the XEN® Gel Stent Implantation

XEN® Gel Stent (Allergan) is used to lower intraocular pressure in patients with primary open angle glaucoma, where previous treatments have failed.



Vera Fixation Hook

ES-0725

- ① Provides easy eye fixation and manipulation through a standard cornea side port incision during XEN® Gel Stent Implantation

The small cone in the front of the Vera Hook enables an easy and safe entry into the side port incision and provides holding resistance to inhibit the instrument from slipping out. A 1.8mm ball towards the back end prevents the full instrument from entering the eye through a 0.9mm side port incision while allowing rotation to any angle to accommodate various holding positions during XEN® Gel Stent Implantation

For Ceramic Coating please add "C"



Mod. Krumeich Eye Speculum for adults,
length 7.0cm, adjustable angulation,
with lock and 15.0mm open blades

ES-0287

- With superior at the internal and inferior at the external part removed valve-bridges for an optimal use during the XEN® Gel Stent Implantation

The XEN® Gel Stent Applicator can not stick to the inferior valve-bridge during application, especially when inserting the applicator

The speculum can be removed shortly after the application, as the superior valve-bridge can not exert any pressure on the bulbus when removing the speculum and thereby not delocalize the stent out of his position

The speculum can, by the adjustable arms, be used both temporal and nasal, on the right and the left eye, and is also perfectly adaptable to the facial shape of the patient, so that the speculum does not protrude and the surgeon has a free temporal access.

By the adjustability of the speculum, it can even be used to lift the eye lid, which is resulting in relief of the pressure on the eye, already during surgery

Through the safe locking mechanism, the speculum can be fixed for an optimal and secure spreading of the eyelids during surgery




Barrett Eye Speculum adjustable,
14.5mm open blades

ES-0311

- The Barrett Type adjustable Speculum is a alternative speculum for use during XEN® Gel Stent Implantation. It comes with 14.5mm open blades, is angled to rest only nasally and adjustable with a thumb screw

For Ceramic Coating please add "C"



Barraquer Wire Speculum (standard),
10.0mm blades

ES-0010

Barraquer Wire Speculum (standard),
15.0mm blades

ES-0020

Barraquer Wire Speculum (heavy wire),
15.0mm blades

ES-0030

For Ceramic Coating please add "C"




Temporal Approach, Wire Speculum,
with spring on nasal side,
15.0mm open blades

ES-0050

Temporal Approach, Wire Speculum,
with spring on nasal side,
15.0mm closed blades

ES-0060




Barraquer Wire Speculum,
10.0mm solid blades

ES-0070

Barraquer Wire Speculum,
15.0mm solid blades

ES-0080

For Ceramic Coating please add "C"



Kratz-Barraquer Wire Speculum,
15.0mm, standard wire

ES-0090

Kratz-Barraquer Wire Speculum,
11.0mm, solid blades

ES-0100



Barraquer Wire Eye Speculum,
14.0mm open blades,
length 4.0cm

ES-0122

For Ceramic Coating please add "C"



Lieberman Aspirating Speculum complete,
15.0mm blades

ES-0130

Replacement silicone tubing

ES-0140



Lieberman Eye Speculum, with Aspiration and Irrigation ports,
15.0mm blades

ES-0150

Replacement silicone tubing

ES-0160

For Ceramic Coating please add "C"



Sauer Lid Speculum, for Infants,
11.0mm blades

ES-0170

Sauer Lid Speculum, for Infants,
6.0mm blades


ES-0180

Sauer Lid Speculum, for premature Infants,
4.0mm blades

ES-0190

Sauer Lid Speculum, for premature Infants,
2.0mm blades


ES-0200



Infant Lid Speculum,
4.0mm blades

ES-0210

For Ceramic Coating please add "C"



Barraquer Infant Wire Speculum,
9.0mm blades

ES-0220



Cook Eye Speculum #1,
8.0mm blades

ES-0230

Cook Eye Speculum #2,
11.0mm blades

ES-0240

Cook Eye Speculum #3,
15.0mm blades

ES-0250

For Ceramic Coating please add "C"




Murdoch Eye Speculum #2,
12.0mm blades

ES-0260

Murdoch Eye Speculum #3,
15.0mm blades

ES-0270



Mellinger Eye Speculum,
open wire, 16.0mm blades

ES-0280

Mellinger Eye Speculum, for children
open wire, 11.0mm blades

ES-0290



Krumeich Eye Speculum for adults,
15.0mm open blades, length 7.0cm,
adjustable angulation, with lock

ES-0285

- can be used temporal and nasal
- perfectly adaptable to the facial shape of the patient
- free temporal access without protruding parts
- safe locking mechanism for an optimal and secure spreading during surgery

For Ceramic Coating please add "C"



Williams Eye Speculum with lock screw,
11.0mm blades

ES-0300

Williams Eye Speculum with lock screw,
15.0mm blades

ES-0310



Lancaster Eye Speculum with locking nut,
16.0mm solid blades

ES-0320



Lieberman Eye Speculum, adjustable mechanism,
10.0mm open wire blades

ES-0330

Lieberman Eye Speculum, adjustable mechanism,
15.0mm open wire blades

ES-0340

Lieberman Eye Speculum, adjustable mechanism,
7.0mm open wire blades

ES-0341

For Ceramic Coating please add "C"



Lieberman Eye Speculum, adjustable mechanism,
15.0mm V-shaped wire blades

ES-0350



Lieberman Eye Speculum, adjustable mechanism,
10.0mm solid blades

ES-0360

Lieberman Eye Speculum, adjustable mechanism,
15.0mm solid blades

ES-0370

Kraff Eye Seculum temporal, adjustable,
14.0mm solid blades

ES-0371



Ginsberg Eye Speculum,
15.0mm blades

ES-0380

For Ceramic Coating please add "C"



Lieberman Reversible Eye Speculum, adjustable mechanism,
for nasal or temporal placement,
15.0mm solid blades

ES-0390



Lieberman Eye Speculum, Brown Triple,
12.0mm blades

ES-0400



Lieberman Eye Speculum,
15.0mm round blades

ES-0410

Lieberman Eye Speculum,
8.0mm round blades,
length 6.0cm

ES-0420

For Ceramic Coating please add "C"



Castroviejo Eye Speculum, adjustable mechanism,
14.0mm blades

ES-0430

Castroviejo Eye Speculum, adjustable mechanism,
16.0mm blades

ES-0440



Fishkind Eye Speculum, adjustable mechanism,
15.0mm blade


ES-0450



Stevenson Lacrimal Sac Retractor,
with 3x3 blunt prongs,
15.0mm blades

ES-0460

For Ceramic Coating please add "C"



Flieringa Scleral Fixation Ring, set of 8 sizes
Ø 15.0mm - 22.0mm (1.0mm increments)

ES-0500

To receive individual rings please add the requested 'mm' at the end of the article number, e.g. ES-0500-15MM



McNeil-Goldman Scleral Fixation Ring & Blepharostat,
medium 16.0mm

ES-0520

For Ceramic Coating please add "C"



Schepens Scleral Depressor 7.0mm wide tip,
small thimble, #8

ES-0540

Schepens Scleral Depressor 7.0mm wide tip,
medium thimble, #10

ES-0550

Schepens Scleral Depressor 7.0mm wide tip,
large thimble, #12

ES-0560



Gass Scleral Marker,
Ø 1.5mm marking tip

ES-0570



Schocket Scleral Depressor,
with 6.0mm wide bar & 2.5mm teardrop shaped ends,
length 13.5cm

ES-0580

For Ceramic Coating please add "C"



Jaffe Wire Lid Retractor,
(pair) 15.0mm blades

ES-0590

Jaffe Wire Lid Retractor,
(pair) 13.0mm blades

ES-0600



Desmarres Lid Retractor, Fig. 0,
11.0mm

ES-0610

Desmarres Lid Retractor, Fig.1,
14.0mm

ES-0620

Desmarres Lid Retractor, Fig. 2,
15.0mm

ES-0630

For Ceramic Coating please add "C"



Jaeger Lid Plate,
length 11.0cm

ES-0640



Wright Fascia Needle,
1.0 mm x 6.0 mm oval hole,
length 13.0cm

ES-0650



Needle & Gouge,
foreign body spud

ES-0660



Helvoston Tissue Retractor,
7.0mm

ES-0670

Helvoston Tissue Retractor,
9.0mm

ES-0680

Helvoston Tissue Retractor,
11.0mm

ES-0690

For Ceramic Coating please add "C"



Agricola Lacrimal Sac Retractor,
3x3 semi-sharp prongs,
length 4.5cm

ES-0700



Guthrie Fixation Hook, double, sharp,
small, length 13.0cm

ES-0710

Guthrie Fixation Hook, double, sharp,
large, length 13.0cm

ES-0720



Rollet Lacrimal Sac Retractor, 4 prongs,
sharp, length 13.5cm

ES-0730

Rollet Lacrimal Sac Retractor, 4 prongs
blunt, length 13.5cm

ES-0740



Axenfeld Hook, 2 prongs,
blunt, length 14.0cm

ES-0760

Axenfeld Hook, 3 prongs,
sharp, length 14.0cm

ES-0770

Axenfeld Hook, 3 prongs,
blunt, length 14.0cm

ES-0780



Graefe Muscle Hook, small,
Size 1

ES-0790



Graefe Muscle Hook, medium,
Size 2

ES-0800



Graefe Muscle Hook, large,
Size 3

ES-0810



Helveston Muscle Hook extra delicate, small,
with blunt tip, 8.0mm

ES-0820



Helveston Muscle Hook, extra delicate, medium,
with blunt tip, 10.0mm

ES-0830



Helveston Muscle Hook, extra delicate, large,
with blunt tip, 12.0mm

ES-0840

For Ceramic Coating please add "C"



Jameson Muscle Hook,
small

ES-0850



Jameson Muscle Hook,
medium

ES-0860



Jameson Muscle Hook,
large

ES-0870



Green Muscle Hook

ES-0880

Stevens Tenotomy Hook,
blunt small

ES-0890



Gass Retinal Detachment Hook,
with oval opening

ES-0900

Helveston „Teaser“ Hook, extra delicate,
blunt

ES-0910

For Ceramic Coating please add "C"



Jaeger Strabismus Hook,
Fig. 1

ES-0920

Jaeger Strabismus Hook,
Fig. 2

ES-0930



Manson Double Ended Strabismus Hook

ES-0940

For Ceramic Coating please add "C"



Urretts-Zavalía Depressor, for fixation & pick

ES-0960



Iris Retractor 10.0mm blue micro hooks,
6 retractors in autoclavable case, reusable

ES-1050

Iris Retractor 10.0mm blue micro hooks,
1 set of 5 retractors, disposable, sterile

ES-1060

ES-1050:



Howard Instruments II, Inc.
Alabama, USA

ES-1060:



SIDAPHARM
Thessaloniki, Greece

For Ceramic Coating please add "C"

TREPHINES
SPUDS & GOUGES
NEEDLES
CORNEAL SCARIFIERS
DISSECTORS
BURRS
PUNCHES
IOL CUTTERS
KNIVES
BLADE HANDLES
DIAMOND KNIVES
FOREIGN BODY-
INSTRUMENTS



Tooke Corneal Knife,
straight blade & curved cutting edge,
3.0x18.0mm, length 11.5cm

C-0010



Gill Corneal Knife,
slightly curved, curved cutting edge

C-0030



Fukasaku „Hokey“ Knife,
length 11.5cm

C-0040

For Ceramic Coating please add "C"



Bonn Scleral Knife,
4.0mm blade,
length 12.5cm

C-0050



Paton Corneal Dissector,
length 11.5cm

C-0060



Kuhnt Corneal Scarifier,
Fig. 1, length 11.5cm

C-0070



Kuhnt Corneal Scarifier,
Fig. 2, length 11.5cm

C-0080

For Ceramic Coating please add "C"



Foreign Body Needle, straight,
length 12.5cm

C-0210

Foreign Body Needle, curved,
length 12.5cm

C-0220



Plange Foreign Body Needle,
length 12.0cm

C-0240



Francis Foreign Body Spud,
length 12.0cm

C-0250



La Forces Foreign Body Spud,
length 12.5cm

C-0260

For Ceramic Coating please add "C"



Algerbrush II Corneal Rust Ring Remover, complete C-0280

Burr 0.5mm C-0290

Burr 1.0mm C-0300

Pterygium Drill, round, 3.5mm burr, battery powered, complete C-0310

Replacement burr with chuck C-0320

All products on this page:

 Alger Inc.
Lago Vista, USA



MICS Diamond Knife,
for 2.2mm incision

C-0362



Diamond Phaco Knife, straight, keratome shaped blade,
2.5mm

C-0370



Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 2.5mm C-0430

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, 2.8mm C-0440

Phaco Diamond Knife, angled, keratome shaped blade & 4 cutting edges, trapezoid from 2.7mm to 3.2mm C-0470




Diamond Knife with 30 ° cutting edge, 1.0mm wide C-0500

Diamond Knife with 45 ° cutting edge, 1.0mm wide C-0510

Diamond Knife with trie-faced, 1.0mm wide C-0520

Diamond Knife Spear, 1.0mm wide C-0530



Universal Trephine Handle

C-0550



Trephine Blade „short model“,
6.0mm - 9.5mm Ø (0.25mm increments)

C-0560 + size



Bard-Parker Scalpel Handle #3

C-0670



Bard-Parker Scalpel Handle #9 C-0680



Chuck Scalpel Handle round, serrated, 10.0cm C-0690

Bard-ParkerScalpel blade, sterile, disposable PK/100 #11 C-0700

Bard-Parker Scalpel blade, sterile, disposable PK/100 #15 C-0710

Bard-Parker Scalpel blade, sterile, disposable PK/10 #65 C-0711

Bard-Parker Scalpel blade, sterile, disposable PK/10 #64 C-0712

Bard-Parker Scalpel blade, sterile, disposable PK/10 #69 C-0714

C-0700, C-0710, C-0711, C-0712, C-0714:



Luntz-Dodick Trabeculectomy Punch,
1.0 x 0.5mm

C-0720



Kelly Descemet Membrane Punch,
1.0 x 0.75mm

C-0730



Soft IOL Cutter,
6.0mm long serrated blades,
for cutting soft IOL's

C-0750

a1 medical HOOKS,
HOOKS, SPATULAS,
MANIPULATORS,
FIXATION RINGS,
KERATOMES, LOOPS, MARKERS,
MARKING PEN,
PHACO CHOPPERS,
PREASSURE PLATES, SPOONS





SINGLE HOOK



Choose a favorite color:

Get in touch with our newly designed, patented handle, which has been developed to enhance the performance of surgeons! Feel the difference, with our new ergonomic handle. The titanium base material is the best there is. It is half the weight of stainless steel, causes less trauma and is antimagnetic. It helps to reduce fatigue during extended operations thanks to its ultralight weight. Our patented handle is ergonomically round and has perfect balance during surgery. Its projecting grooves not only provide an excellent grip, but the handle is also easy to maintain and sterilize. And finally, to help you save on costs, we have designed the handle so that the instrument tips can be replaced.

Total length 11.0cm



Titanium Handle only, for Single Hook System, used with tip variations HS-0001TSX*

Kuglen Iris Hook and Lens Manipulator, push-pull model, angled HS-0000TSX*

Replacement tip HS-0000TS-TIP A-1



Rosen Phaco Splitter, blunt tips, 10.0mm angled HS-0010TSX*

Replacement tip HS-0010TS-TIP A-2



Rosen Phaco Splitter, blunt tips, 10.0mm angled 60° HS-0020TSX*

Replacement tip HS-0020TS-TIP A-3



Haefliger Phaco cleaver, 45° angled, right HS-0030TSX*

Replacement tip HS-0030TS-TIP A-4



Nagahara Phaco-Chopper, 1.25mm tip, offset blade, angled 60°, right HS-0040TSX*

Replacement tip HS-0040TS-TIP A-5

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)



Nagahara Phaco-Chopper, 1.5mm tip, offset blade,
angled 60°, right

HS-0050TSX*

Replacement tip HS-0050TS-TIP

A-6




Pisacano Nucleus Rotator,
vertical flat disc

HS-0060TSX*

Replacement tip HS-0060TS-TIP

A-7



Pisacano Nucleus Rotator,
horizontal flat disc

HS-0070TSX*

Replacement tip HS-0070TS-TIP

A-8



Bechert Nucleus Rotator, blunt forked tip, Y-horizontal,
7.0mm angled

HS-0080TSX*

Replacement tip HS-0080TS-TIP

A-9



Drysdale Nucleus Manipulator

HS-0090TSX*

Replacement tip HS-0090TS-TIP

A-10

Drysdale Nucleus Manipulator, extra delicate

HS-0091TSX*

Replacement tip HS-0091TS-TIP

A-39



Lester IOL-Manipulator,
straight

HS-0110TSX*

Replacement tip HS-0110TS-TIP

A-11

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)



Lester IOL-Manipulator,
angled, 10.0mm

HS-0120TSX*

Replacement tip HS-0120TS-TIP

A-12



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
straight

HS-0140TSX*

Replacement tip HS-0140TS-TIP

A-14



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
angled

HS-0150TSX*

Replacement tip HS-0150TS-TIP

A-15



Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip,
straight

HS-0160TSX*

Replacement tip HS-0160TS-TIP

A-16



Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip,
angled

HS-0170TSX*

Replacement tip HS-0170TS-TIP

A-17

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)



Barraquer Iris Spatula, smooth Ø 0.5mm,
angled, 12.0mm

HS-0220TSX*

Replacement tip HS-0220TS-TIP

A-22



Barraquer Iris Spatula, smooth Ø 0.25mm,
angled, 12.0mm

HS-0230TSX*

Replacement tip HS-0230TS-TIP

A-23



Osher Y-Hook Lens Manipulator,
straight

HS-0240TSX*

Replacement tip HS-0240TS-TIP

A-24



Osher Y-Hook Lens Manipulator,
angled

HS-0250TSX*

Replacement tip HS-0250TS-TIP

A-25



Ligabue Nucleus Manipulator & Splitter,
12.0mm curved arm flat sandplasted oive shaped
90° rotated

HS-0256TSX*

Replacement tip HS-0256TS-TIP

A-45

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)



DOUBLE HOOK

You can either choose a tip from our selection below or choose a tip according to your needs..

Total length 12.0cm

		Titanium Handle only, for Double Hook System, used with tip variations	HS-0002TSX*
1.	2.		
		Rowen Phaco Spatula & Bechert Nucleus Rotator	HS-0260TSX*
		1. Replacement tip HS-0260TS-Tip Rowen	A-26
		2. Replacement tip HS-0260TS-Tip Bechert	A-27
		Nucleus Rotator (push-pull) & Dividing Hook	HS-0280TSX*
		1. Replacement tip HS-0280TS-Tip Push-Pull	A-30
		2. Replacement tip HS-0280TS-Tip Dividing	A-31
		Chang Phaco Chopper, for left hand	HS-0310TSX*
		1. Replacement tip HS-0310TS-Tip Microfinger	A-32
		2. Replacement tip HS-0310TS-Tip Quick-Chop	A-33
		Chang Phaco Chopper, for right hand	HS-0320TSX*
		1. Replacement tip HS-0320TS-Tip Microfinger	A-34
		2. Replacement tip HS-0320TS-Tip Quick-Chop	A-35
		Masket Phaco Spatula	HS-0330TSX*
		1. Replacement tip HS-0330TS-Tip Splitter	A-36
		2. Replacement tip HS-0330TS-Tip Spatula	A-3

* replace "X" with either B (for blue), M (for magenta), Y (for yellow), G (for green), N (for neutral)



Haefliger Phaco cleaver, 45° angled,
right

HS-0340



Phaco Manipulator, finger shaped,
curved tip

HS-0350



Nagahara Phaco-Chopper, 1.25mm tip, offset blade,
angled 60°, right

HS-0360



Nagahara Phaco-Chopper, 1.5mm tip, offset blade,
angled 60°, right

HS-0370



Lund Iris Hook,
14.0mm tip, 45° angled

HS-0380



Dardenne Iris Hook,
10.0mm, angled shaft

HS-0390




Miyohshi Nucleus Divider,
straight

HS-0400



Miyoshi Nucleus Divider,
left

HS-0410



Miyoshi Nucleus Divider,
right

HS-0420



Neuhann Nucleus Rotator,
10.0mm angled tip, blunt forked tip

HS-0430



Neuhann Nucleus Spatula / Phaco Instrument,
curved spatulated tip with 45° shaft

HS-0440



Manipulator, Spatula 0.50mm x 0.25mm, round, Ø 0.25mm,
angled

HS-0450

For Ceramic Coating please add "C"



Barraquer Iris Spatula, round smooth, Ø 0.5mm,
angled, 12.0mm

HS-0460



Barraquer Iris Spatula, round smooth, Ø 0.25mm,
angled, 12.0mm

HS-0470



Hirschman Iris Hook, blunt,
angled

HS-0480



Osher Y-Hook Lens Manipulator,
straight

HS-0490



Osher Y-Hook Lens Manipulator,
angled

HS-0500



Deitz ICL Manipulator, blunt forked tip Y-Type,
angled

HS-0510



Koch Stop & Chop Manipulator, 1.5mm tip,
for „Stop & Chop technique“.

HS-0521

For Ceramic Coating please add "C"



Manipulator, Spatula 1.0 x 0.25mm,
angled

HS-0530



Knolle Nucleus Spatula, 0.5mm wide,
straight

HS-0540



Hunkeler IOL Manipulator, with Ø 0.3mm semi-ball tip,
angled

HS-0550



Koch Nucleus Spatula, sharp & smooth tip,
gently curved

HS-0570



Koch Nucleus Spatula, duckbill tip,
gently curved

HS-0571



Kuglen Iris Hook and Lens Manipulator, push-pull model, angled

HS-0580

Kuglen Iris Hook, delicate tip

HS-0581



Rosen Phaco Splitter, wedge shaped blunt tip, inferior edges, angled 10.0mm

HS-0590




Rosen Phaco Splitter, wedge shaped blunt tip, inferior edges, angled 10.0mm, 60°

HS-0600



Pisacano Nucleus Rotator, vertical flat disc

HS-0610



Pisacano Nucleus Rotator, horizontal flat disc

HS-0620



Pisacano Nucleus Rotator, 11.0mm angled, olive round tip

HS-0621



Bechert Nucleus Rotator, blunt forked tip, Y-horizontal, 7.0mm angled

HS-0630



Drysdale Nucleus Manipulator, paddle shaped tip

HS-0640

For Ceramic Coating please add "C"



Lester IOL-Manipulator, Ø 0.25mm,
straight

HS-0650



Lester IOL-Manipulator, Ø 0.25mm,
10.0mm angled

HS-0660



Rentsch IOL Manipulator, blunt,
straight

HS-0670



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
straight

HS-0680



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip,
angled

HS-0690



Sinskey Lens Manipulating Hook, Ø 0.2mm, blunt tip,
angled, 180° reversed

HS-0691



Sinskey Lens Manipulating Hook, Ø 0.25mm blunt tip,
straight

HS-0700



Sinskey Lens Manipulating Hook, Ø 0.25mm blunt tip,
angled

HS-0710

For Ceramic Coating please add "C"



Fenzl Lens Manipulating Hook, Ø 0.15mm tip,
angled

HS-0730



Maltzman-Fenzl Lens Manipulator Hook, Ø 0.15mm,
V-shaped tip, straight

HS-0740



Maltzman-Fenzl Lens Manipulator Hook, Ø 0.15mm,
V-shaped tip, angled

HS-0750




Kagawa Phaco-Chopper, offset blade,
angled 90°, right

HS-0760



Graether Collar Button Micro Iris Hook,
straight

HS-0770



Graether Collar Button Micro Iris Hook,
angled

HS-0780



Barraquer Iris Spatula, with flat disc,
angled 8.0mm

HS-0790



Barraquer Iris Spatula, with flat disc,
angled 10.0mm

HS-0800



Barraquer Iris Spatula, with flat disc,
angled 12.0mm

HS-0810



Barraquer Iris Spatula, with flat disc,
angled 15.0mm

HS-0820

For Ceramic Coating please add "C"



Bonn Micro Iris Hook, extra delicate, blunt tip, straight

HS-0830



Lieberman Micro Finger Nucleus Manipulator, for left hand

HS-0850



Gimbel Nucleus Spatula, round, smooth, Ø 0.25mm, angled 30 °

HS-0860



Fukasaku Snapper Hook, with blunt tip & front notch

HS-0870



Fukasaku Snapper Hook, angled

HS-0880



LASIK Spatula rounded anterior & flat posterior surface with dissecting tip

HS-0890



Shepard Iris Hook, blunt, 0.5mm hook,
straight

HS-0900



Bonn Micro Iris Hook, blunt, 0.3mm hook,
extra delicate

HS-0910



Kuglen Lens Manipulator push-pull, "Clover Leaf" model,
straight

HS-0930



Kuglen Iris Hook & Lens Manipulator, push-pull model,
angled 10.0mm

HS-0940



Jaffe Lens Manipulating Hook, Ø 0.15mm blunt tip,
angled

HS-0960



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, straight

HS-0970



Sinskey II Lens Manipulating Hook, Ø 0.2mm, blunt tip, angled

HS-0980



Sinskey Lens Manipulating Hook, Ø 0.25mm, blunt tip, straight

HS-0990



Jaffe Lens Spatula, smooth, 0.5mm wide blade, angled 10.0mm

HS-1000



Obstbaum Synechiae Spatula, 0.25mm wide blade, 10.0mm angled

HS-1010



DALK Corneal Dissector,
12.0mm long spatula, flat posterior surface,
blunt, beveled tip

HS-1022

For Ceramic Coating please add "C"



Rowen Phaco Spatula & Bechert Nucleus Rotator

HS-1030



Nucleus Rotator & Dividing Hook

HS-1040



Nucleus Rotator (push-pull) & Dividing Hook

HS-1050



Pfaffl Femto Flaplifter

HS-1051

For Ceramic Coating please add "C"



Castroviejo Spatula, double ended,
10.0mm & 15.0mm long

HS-1060



Chang Phaco Chopper, a combination of a modified MicroFinger,
and a pointed „Quick Chop“ tip,
for left hand

HS-1090



Chang Phaco Chopper, a combination of a modified MicroFinger,
and a pointed „Quick Chop“ tip,
for right hand

HS-1100



Masket Phaco Spatula, a combination of the Rosen Phaco Splitter
and a 12.0mm long round spatula with disc shaped tip

HS-1110

For Ceramic Coating please add "C"



Culler Iris Spatula,
1.0mm wide

HS-1120

Culler Iris Spatula,
2.0mm wide

HS-1130



Bangerter Iris Spatula,
1.0 x 0.25mm

HS-1140

Bangerter Iris Spatula,
1,5 x 0.25mm

HS-1150

Bangerter Iris Spatula,
2,0 x 0.25mm

HS-1160

For Ceramic Coating please add "C"



Giessen Orbital Spatula

HS-1170



Schepens Orbital Spatula

HS-1250

For Ceramic Coating please add "C"



Arruga Globe Retractor, cup shaped,
length 14.5cm

HS-1260



Wells Enucleation Spoon

HS-1270



Paton Spatula & Spoon

HS-1290

Lewis Lens Loop,
small

HS-1300

Lewis Lens Loop,
large

HS-1310



New Orleans Lens Loop

HS-1320

Wilder Lens Loop,
15.0mm x 4.0mm smooth loop

HS-1330

Snellen Lens Loop,
15.0mm x 4.0mm, serrated

HS-1340

Arlt Lens Loop,
serrated, 3.0mm wide

HS-1350

Keener-Arlt Lens Loop, serrated, 3.0mm wide,
strong curved

HS-1360



Kansas Nucleus Loop, serrated, 3.0mm wide

HS-1370

For Ceramic Coating please add "C"



Evisceration Spoon, bunge small,
6.0mm

HS-1460

Evisceration Spoon, bunge small,
8.0mm

HS-1470

Evisceration Spoon, bunge medium,
10.0mm

HS-1480

Evisceration Spoon, bunge large,
12.0mm

HS-1490

For Ceramic Coating please add "C"



Meyerhoefer Chalazions Curette,
size 00, Ø 1.0mm

HS-1500

Meyerhoefer Chalazions Curette,
size 0, Ø 1.5mm

HS-1510

Meyerhoefer Chalazions Curette,
size 1, Ø 1.8mm

HS-1520

Meyerhoefer Chalazions Curette,
size 2, Ø 2.0mm

HS-1530

Meyerhoefer Chalazions Curette,
size 3, Ø 2.5mm

HS-1540

Meyerhoefer Chalazions Curette,
size 3.5, Ø 3.0mm

HS-1545


Meyerhoefer Chalazions Curette,
size 4, Ø 3.5mm

HS-1550

Meyerhoefer Chalazions Curette,
size 5, Ø 4.0mm

HS-1560

For Ceramic Coating please add "C"




MARKING
PATTERN

Bores Axis Marker,
used with Mendez Degree Gauge

HS-1570

Toric IOL Marker
two 2.55mm radial marks extending to 11.0mm

HS-1571




MARKING
PATTERN

Mendez Degree Gauge 12.0mm,
0 - 180° in 10° increments

*HS-1580

For Ceramic Coating please add "C"
*Not available in Ceramic Coating



MARKING
PATTERN

Fine-Thornton Fixation Ring,
Ø 16.0mm, 20 blunt teeth, 15°

HS-1600

Hofmann-Thornton Fixation Ring,
Ø 13.0mm, 8 blunt teeth, 15°

HS-1610

For Ceramic Coating please add "C"
*Not available in Ceramic Coating




Globe Fixation Ring

HS-1620

For Ceramic Coating please add "C"

* Not available in Ceramic Coating



MARKING
PATTERN

Anis Radial Marker,
8 blades

HS-1640

Anis Radial Marker,
12 blades

HS-1650

Anis Radial Marker,
16 blades

HS-1660



MARKING
PATTERN

Radial Marker

12 blades

HS-1663

16 blades

HS-1664



MARKING PATTERN

Mendez LASIK Corneal Marker

*HS-1700

For Ceramic Coating please add "C"

* Not available in Ceramic Coating



Hoffer Zone Marker, 3.0mm	HS-1710
------------------------------	---------

Hoffer Zone Marker, 3.25mm	HS-1720
-------------------------------	---------

Hoffer Zone Marker, 3.5mm	HS-1730
------------------------------	---------

Hoffer Zone Marker, 4.0mm	HS-1750
------------------------------	---------

For Ceramic Coating please add "C"



Hoffer Zone Marker,
5.0mm

HS-1790

Hoffer Zone Marker,
5.5mm

HS-1800

Hoffer Zone Marker,
6.0mm

HS-1810

Hoffer Zone Marker,
8.0mm

HS-1830

Hoffer Zone Marker,
9.0mm

HS-1840

For Ceramic Coating please add "C"



Bores Optic Zone Marker,
4.0mm

HS-1880

Bores Optic Zone Marker,
6.0mm

HS-1930

Bores Optic Zone Marker,
7.0mm

HS-1940



Bore's Optic Zone Marker,
8.0mm

HS-1950

Bore's Optic Zone Marker,
9.0mm

HS-1960



Lasek Epithelial Trephine,
10.0mm

HS-1990

For Ceramic Coating please add "C"



Camellin Style, LASEK Alcohol Well,
Ø 8.5mm

HS-2000

Camellin Style, LASEK Alcohol Well,
Ø 9.5mm

HS-2010



Buratto LASIK Flap Protector

HS-2020




Fechtner Trabeculectomy Marker

HS-2021



MARKING
PATTERN

For Ceramic Coating please add "C"



Membran-Peeler, blunt,
angled 110°,
length 13.0cm

HS-2040



Membrane Hook,
length 13.0cm

HS-2060



Membrane Scratcher,
length 13.0cm

HS-2080



Eckardt Pressure Plate, 23G/ 0.6mm,
for Vitrectomy System, left hand, used with F-1540,
length 9.5cm

HS-2100

Eckardt Pressure Plate, 23G/ 0.6mm,
for Vitrectomy System, right hand, used with F-1540,
length 9.5cm

HS-2110

For Ceramic Coating please add "C"

* Not available in Ceramic Coating



Maloney Intra-Operative Keratometer,
titanium

HS-2120



Marking Pen, sterile,
„single use“ (Box of 5 psc.)

HS-2130

HS-2130:


 Distributed by a1 medical

For Ceramic Coating please add "C"

**FORCEPS
MIC FORCEPS
TITANIUM FORCEPES
BIPOLAR FORCEPS
ERGO GRIP FORCEPS
TC FORCEPS
SIGNUM MICS FORCEPS
SMILE GRIP FORCEPS**




GENERAL INFORMATION




Pplatform Forceps

The usable area of the platform with guaranteed parallel closure. In this way a secure grip on the sutures can be ensured.



Pierse Type Forceps

These forceps are ideal for tissue fixation, manipulation of limbal area, corneal lip and minimize trauma to the tissue.




Ring Forceps

These ring forceps have been designed for perfect grasping soft tissue & handling of the conjunctiva. Additionally it provides an optimal counter-pressor function.



Tothed Forceps

These delicate suture forceps are ideal for fixation & manipulating of conjunctiva, limbal area or other delicate tissue & grasping fine suture tying.



Capsulorhexis Forceps

These forceps with triangular jaws, small platform and ultra fine grasping tips are optimal to perform the capsulorhexis.

...sharp tip is used as a cystotome.

...blunt tip is used to grip, lift and pull the iris.

ULTRA LIGHT ERGO GRIP

...Perfect balance for the surgeon's use

Every surgical operation requires good skill and sensitivity. Operations with longer instruments under the microscope, in particular require thoroughness and precision work. Long operations are especially challenging for the surgeon. In order to facilitate the work of the surgeon, a1 medical has developed an ergonomically shaped and light-weight line of instruments.

The "Ultra-Light Ergo Grip". At a1 medical, our micro instruments are manufactured only from the finest stainless steel; their protruding ribs give you a perfect grip, and they can be easily cleaned. Their ergonomic design is not only pleasing to the eye but also offers perfect balance in the hand. Particularly noteworthy is the price-performance ratio as compared to titanium. Instruments made of titanium usually cost almost twice as much as stainless steel instruments.

Due to their excellent design, "Ultra-Light Ergo Grip" instruments weigh hardly any more than those made of titanium and are thus a good cost alternative.

Of course, our "Ultra-Light Ergoline" is also available with features such as TC coating and ceramic coating.

ADVANTAGES AT A GLANCE

- provides perfect balance during surgery
- ergonomic & attractive design
- excellent grip




- high quality stainless steels
- easy to maintain & sterilize
- saves costs

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.5cm




Ergo Grip Tennant Tying Forceps, straight,
with 6.0mm platform, for 9.0-11.0 sutures

F-4800




Ergo Grip Tennant Tying Forceps, curved,
with 6.0mm tying platform, for 9.0-11.0 sutures,

F-4805



Ergo Grip Kelman-McPherson Forceps, angled,
7.5mm long smooth jaws
(ideal for implanting IOL)

F-4810




Ergo Grip McPherson Suturing Forceps, angled,
with 5.0mm tying platform

F-4820

For Ceramic Coating please add "C"




With guide pin & stop pin, Ø 8.0mm, length 11.5cm



Ergo Grip Steinert Direct Action Forceps, angled,
with 7.0mm long cone-shaped jaws

F-4880



Ergo Grip Utrata Capsulorhexis Forceps, angled,
triangular shaped tips

F-4910



Ergo Grip Lehner Capsulorhexis Forceps, curved,
triangular shaped tips

F-4920



Ergo Grip Nevyas Capsulorhexis Forceps, angled,
cystotome shaped tips

F-4940



Ergo Grip Bonn Forceps, straight, with 6.0mm tying platform,
1x2 teeth, 0.12mm 90° angle

F-4950

SIGNUM MICS FORCEPS



Workingends made of titanium are too flexible, sensitive and deform more easily.

> Therefore the workingend of the SIGNUM MICS FORCEPS are made of stainless steel.

A handle made of stainless steel would be too heavy.

> Therefore the handle of the SIGNUM MICS FORCEPS are made of titanium.

STAINLESS STEEL

TITANIUM



ULTRA LIGHT ERGO GRIP-HANDLE, IN TITANIUM

ADVANTAGES AT A GLANCE

- provides perfect balance during surgery
- ergonomic & attractive design
- excellent grip & haptic
- fatigue-free work
- easy to maintain and to sterilize

For 1.5mm incisions



SIGNUM LEHNER MICS Forceps, angled, curved,
Cross Action, 1.5 mm, graduated at 2.5 and 5.0mm,
9.5mm shanks from pivot to tip

F-5000TS



SIGNUM GIANETTI MICS Forceps, angled, curved,
Cross Action, 1.5mm, 9.5mm shanks from pivot to tip

F-5010TS



SIGNUM INAMURA MICS Forceps, angled, curved,
Cross Action, 1.5 mm, 90°, sharp Tips, graduated at 2.5 and 5.0mm,
9.5mm shanks from pivot to tip

F-5020TS

TC COATING


TC-Coating



When coated with tungsten carbide, tips of instruments last significantly longer since tungsten carbide is much harder than stainless steel. This allows an atraumatic and save handling of delicate tissue, vessels and fine needles.



With guide pin & stop pin, Ø 8.0mm, length 11.5cm



Ergo GripTennant Tying Forceps, straight, with 6.0mm tying platform, for 9.0-11.0 sutures

F-4800TC

SMILE GRIP FLAT HANDLE


Introducing the a1 medical - flat handle design



SMILE GRIP IOL Loading Forceps for the Monarch® II and III cartridges, 11.0cm


F-6000

For Ceramic Coating please add "C"




Harms Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm

F-0010




Barraquer Corneal Utility Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0020




Barraquer Corneal Utility Forceps, with tying platform,
1x2 teeth, 0.3mm, 90° angle,
length 7.5cm

F-0030




Barraquer Corneal Utility Forceps, with tying platform,
1x2 teeth, 0.4mm, 90° angle,
length 7.5cm

F-0040




Barraquer Corneal Utility Forceps,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0050



Micro Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm

F-0070




Micro Colibri Forceps, with tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 7.5cm

F-0080




Troutman Corneal Fixation Forceps,
1x2 teeth, 0.12mm, 45° angle,
length 7,5.cm

F-0100



Barraquer Corneal Utility Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0120




Barraquer Corneal Utility Forceps,
1x2 teeth, 0.12m, 90° angle,
length 7.5cm

F-0140




Suture Forceps (Zurich model) straight, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm

F-0160




Corneal Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0230




Barraquer Forceps for removing sutures, with tying platform,
length 7.5cm

F-0250




Polack Double Corneal Forceps,
with 2x 1x2 teeth, 0.12mm,
length 7.5cm

F-0260



Fixation Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 7.5cm


F-0270



Corneal Colibri Utility Forceps, with tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 7.5cm

F-0280

For Ceramic Coating please add "C"




Gill Iris Forceps,
criss-cross serrated tips,
length 7.5cm

F-0410



Micro Suture Forceps (Bonn model), straight,
5.0mm tying platform
length 7.5cm

F-0440



Iris Forceps (Bonn model), straight, fine serrated,
length 7.5cm

F-0450

For Ceramic Coating please add "C"




Bonn Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0510



Iris Forceps Bonn model, straight,
with 1x2 teeth, 0.12mm, 90° angle,
length 7.5cm

F-0520



Corneal Tying Forceps, straight,
with extra fine tying platform,
length 8.5cm

F-0560

Same in titanium

F-0560T




Corneal Tying Forceps, curved ,
with extra fine tying platform,
length 8.5cm

F-0570

Same in titanium

F-0570T



Corneal Tying Forceps, angled,
with extra fine 7.0mm tying platform,
length 8.5cm

F-0580

Same in titanium

F-0580T

For Ceramic Coating please add "C"



Corneal Fixation Forceps, straight, with tying platform,
1x1 Pierse teeth, 0.1mm,
length 8.5cm

F-0590

Same in titanium

F-0590T



Corneal Fixation Forceps, straight, with tying platform,
1x1 Pierse teeth, 0.2mm,
length 8.5cm

F-0600



Corneal Fixation Forceps, straight,
1x1 Pierse teeth, 0.2mm,
length 8.5cm

F-0610



Corneal Fixation Forceps, straight,
1x1 Pierse teeth, 0.3mm,
length 8.5cm

F-0620



Corneal Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-0650

Same in titanium


F-0650T



Corneal Suturing Forceps, curved, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 8.5cm

F-0680

For Ceramic Coating please add "C"



Corneal Fixation Forceps, straight,
1x1 Pierse teeth, 0.2mm,
length 8.5cm


F-0740



Corneal Fixation Forceps, angled, with
1x1 Pierse teeth, 0.2mm,
length 8.5cm


F-0770

For Ceramic Coating please add "C"



Paufigue Suture Forceps, with 5.0mm tying platform,
standard jaws,
length 8.5cm

F-0830



Paufigue Suture Forceps, with 5.0mm tying platform,
very delicate jaws,
length 8.5cm

F-0840



Suture Forceps (St. Martin model), with 5.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-0850



Suture Forceps (St. Martin model), with 5.0mm tying platform,
1x2 teeth, 0.3mm, 90° angle,
length 8.5cm

F-0860



Suture Forceps (Paufigue model), with 5.0mm tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 8.5cm

F-0870



Paufigue Suture Forceps, with 5.0mm tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 8.5cm

F-0880



Manhattan Eye & Ear Suturing Forceps, with 5.0mm tying platform,
1x2 teeth, 0.5mm, 45° angle,
length 8.5cm

F-0890

For Ceramic Coating please add "C"




Paufigue Suture Forceps, straight, with 5.0mm tying platform, 1x2 teeth, 0.7mm, 45° angle, length 8.5cm

F-0900



Sauer Suturing Forceps, straight, 1x2 teeth, 0.5mm, 45° angle, length 8.5cm

F-0930



Suture Forceps (Tuebingen model), with 5.0mm tying platform very delicate, length 8.5cm

F-0940

For Ceramic Coating please add "C"




McPherson Tying Forceps straight, with 5.0mm tying platform, length 8.5cm

F-0950




McPherson Tying Forceps angled, with 5.0mm tying platform, length 8.5cm

F-0960




Kelman-McPherson Tying Forceps, 7.5mm, angled, with 5.0mm platform, length 8.5cm

F-0970




Bechert-McPherson Tying Forceps, 10.0mm, angled, with 5.0mm tying platform, length 8.5cm

F-0980



Sheets McPherson IOL Forceps, 10.0mm, angled, with serrated tips, length 8.5cm

F-0990




Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 8.5cm

F-1000



Lehner Capsulorhexis Forceps, angled, triangular shaped tips, length 8.5cm

F-1010




For Ceramic Coating please add "C"



Gimbel Capsulorhexis Forceps, angled,
only 3.0mm opening, triangular shaped tips,
length 8.5cm

F-1020




Utrata Capsulorhexis Forceps, angled,
criss-cross serrated, triangular shaped tips,
length 8.5cm

F-1030



McPherson Corneal Forceps, straight, with 5.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-1040




Bonn Tying Forceps, straight, with 5.0mm tying platform,
for 9.0-11.0 sutures,
length 9.5cm

F-1050



Bonn Tying Forceps, curved, with 5.0mm tying platform,
for 9.0-11.0 sutures,
length 9.5cm

F-1060



Bonn Tying Forceps, angled, with 5.0mm tying platform,
for 9.0-11.0 sutures,
length 9.5cm


F-1070



Bonn Forceps straight, long model, with 5.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 9.5cm


F-1090

For Ceramic Coating please add "C"





Thorpe Conjunctiva Fixation Forceps, straight, with tying platform, with 2x3 teeth, 0.6mm, 90° angle, length 9.5cm

F-1120




Lester Fixation Forceps, straight, with 2x3 teeth, 0.6mm, 90° angle, length 9.5cm

F-1130



Castroviejo Tying Forceps, straight, with 5.0mm tying platform, length 11.0cm

F-1140



Kelman McPherson Tying Forceps, 7.5mm angled, with tying platform, length 11.0cm

F-1150




Utrata Capsulorhexis Forceps, angled, triangular shaped tips, length 11.0cm


F-1160

For Ceramic Coating please add "C"




 	Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.12, 45° angle, length 11.0cm	F-1190
 	Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm	F-1200
 	Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.5mm, 45° angle, length 11.0cm	F-1210
 	Castroviejo Suturing Forceps, with 6.0mm tying platform, 1x2 teeth, 0.9mm, 45° angle, length 11.0cm	F-1220
 	Snider-Castroviejo Sutures Forceps, with 6.0mm tying platform, 1x2 teeth, 0.12mm, 45° angle, with lock, length 11.0cm	F-1240
 	Castroviejo Suturing Forceps, 1x2 teeth, 0.3mm, 45° angle, length 11.0cm	F-1250
 	Castroviejo Suturing Forceps, 1x2 teeth, 0.5mm, 45° angle, length 11.0cm	F-1260
 	McPherson Corneal Forceps, with 5.0mm tying platform, 1x2 teeth, 0.4mm, 90° angle, length 11.0cm	F-1290

For Ceramic Coating please add "C"




Castroviejo Suturing Forceps angled, with 6.0mm tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 11.0cm

F-1310



Castroviejo Suturing Forceps angled, with 6.0mm tying platform,
1x2 teeth, 0.5mm, 45° angle,
length 11.0cm

F-1320




Colibri Forceps, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 11.0cm

F-1380




Stern-Castroviejo Fixation Forceps, with 6.0mm tying platform,
1x2 teeth, 0.5mm, 45° angle, with lock,
length 11.0cm

F-1400



Moody Fixation Forceps, right,
1x2 teeth, 0.5mm, 45° angle, with lock,
length 11.0cm


F-1410



Moody Fixation Forceps, left,
1x2 teeth, 0.5mm, angle, with lock,
length 11.0cm

F-1420

For Ceramic Coating please add "C"




Harms Tying Forceps, straight, with 6.0mm tying platform,
length 11.0cm

F-1510




Harms Tying Forceps, curved, with 6.0mm tying platform,
length 11.0cm

F-1530



Removal Cannula Forceps for 23G,
length 11.0cm

F-1540



Dodick Nucleus Cracker, cross-action,
longitudinally serrated paddle shaped jaws,
length 11.0cm

F-1550



Castroviejo Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 11.0cm

F-1560

Same in titanium

F-1560T



Castroviejo Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 11.0cm

F-1570



Castroviejo Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.5 mm, 45° angle,
length 11.0cm

F-1580




Castroviejo Suturing Forceps, with 6.0mm tying platform,
1x2 teeth, 0.9mm, 45° angle,
length 11.0cm

F-1590



Stern-Castroviejo Fixation Forceps,
1x2 teeth, 0.5mm, 45° angle, with lock,
length 11.0cm


F-1610



Moody Fixation Forceps, curved, with
1x2 teeth, 0.5mm, 45° angle, with lock, left,
length 11.0cm


F-1620

For Ceramic Coating please add "C"




Moody Fixation Forceps, curved,
1x2 teeth, 0.5mm 45°angle, with lock, right,
length 11.0cm

F-1630



Colibri Forceps,
1x2 teeth, 0.12mm, 90° angle,
length 11.0cm

F-1670




McPherson Suture Forceps straight, with 5.0mm tying platform,
length 10.5cm

F-1680



Helveston Tying Forceps,
standard smooth jaws,
length 10.5cm

F-1690




Kelman-McPherson Forceps, angled,
very delicate 7.5mm long smooth jaws,
length 10.5cm

F-1710


Same in titanium

F-1710T




McPherson Suturing Forceps, angled,
with 5.0mm tying platform,
length 10.5cm

F-1720




Same in titanium

F-1720T




Gaskin Fragment Forceps, 11.0mm angled,
smooth jaws,
length 10.5cm

F-1740




Bechert Forceps, 10.0mm angled,
length 10.5cm

F-1750



Faulkner Lens Holding Forceps, angled,
11.0mm cone shaped jaws with, rounded surfaces & edges,
length 10.5cm


F-1760



Livernois Lens Forceps, angled,
5.0mm cone-shaped jaws with flat inside surfaces & rounded edges,
length 10.5cm

F-1770

For Ceramic Coating please add "C"




Steinert Direct Action Forceps, angled,
with 7.0mm long cone-shaped jaws,
length 10.5cm

F-1780



Haberle-McPherson Lens Holding Forceps, 7.5mm angled,
smooth jaws, diagonal groove inside tips,
length 10.5cm

F-1790



Spaleck Acrylic IOL Implantation Forceps, angled, 15.0mm,
length 10.5cm

F-1800



Buratto Lens Folding Forceps, convex jaws,
length 10.5cm

F-1810



Utrata Capsulorhexis Forceps, angled,
triangular shaped tips,
length 10.5cm

F-1820

Same in titanium

F-1820T




Lehner Capsulorhexis Forceps, curved,
triangular shaped tip.
length 10.5cm

F-1830

Same in titanium

F-1830T



Jaffe Capsulorhexis Forceps ,angled,
blunt, olive shaped tips,
length 10.5cm

F-1840

For Ceramic Coating please add "C"



LASIK Flap Forceps, curved disc shaped & serrated jaws,
length 10.5cm

F-1860




Bonn Forceps, straight, with tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 10.5cm

F-1910




Bonn Forceps, straight, with tying platform,
1x2 teeth, 0.3mm, 90° angle,
length 10.5cm

F-1920




McPherson Corneal Forceps, straight, with 5.0mm tying platform,
1x2 teeth, 0.4mm 90° angle,
length 10.5cm

F-1930



Kelman-McPherson Forceps, angled 7.5mm, with tying platform,
1x2 teeth, 0.3mm 90° angle,
length 10.5cm

F-1940



Fechtner Conjunctiva Forceps Ø 0.8mm ring,
length 10.5cm

F-2421




Arita Meibomian Gland Compressor,
length 10.5cm

F-1971



Endothelial / Descemetorhexis Forceps, extra-delicate,
reversed triangular tips, angled 75°,
length 10.5cm

F-1975



Jaffe Forceps, straight, delicate jaw,
with 6.0mm tying platform,
length 11.0cm

F-1990




Jaffe Forceps, curved, delicate jaw,
with 6.0mm tying platform,
length 11.0cm

F-2000

For Ceramic Coating please add "C"



MIC Forceps for small incisions



MIC Utrata Capsulorhexis Forceps, triangular grasping tips, extremely thin shanks, for 2.2mm incision, controlled opening, length 10.5cm

F-1960



MIC Corydon Capsulorhexis Forceps, cystotome shaped tips, extremely thin shanks, for 2.2mm incision, controlled opening, length 10.5cm

F-1970



MIC GIANNETTI Capsulorhexis Forceps, curved, for 2.0mm incision, controlled opening, length 10.5cm

F-1976



MIC Utrata Capsulorhexis Forceps, triangular shaped tips, 14.0mm, angled, very delicate shanks, length 8.5cm

F-1003



MIC Jaffe Capsulorhexis Forceps, olive shaped blunt tips, 14.0mm, angled, very delicate shanks, length 8.5cm

F-1004

For Ceramic Coating please add "C"



MIC Forceps for small incisions



MIC Corydon Capsulorhexis Forceps, curved, 12.0mm, very delicate shanks, cystotome shaped tips, for 2.0mm incision, controlled opening, Ø 6.0mm, length 11.0cm

F-3190



MIC GIANNETTI Capsulorhexis Forceps, curved, for 2.0mm incision, controlled opening, Ø 6.0mm, length 11.0cm

F-3200




MIC Utrata Capsulorhexis Forceps, 13.0mm, angled, extra delicate shanks, triangular shaped tips, Ø 6.0mm, length 10.5cm

F-3151

For Ceramic Coating please add "C"



Capsulorhexis Forceps with markings on jaw



Corydon Capsulorhexis Forceps 12.0mm angled, cystotome tips, marking 2.5mm & 5.0mm, Ø 8.0mm, length 11.0cm

F-2981



MIC Forceps for small incisions



MIC Lehner Capsulorhexis Forceps, 12.0mm, angled, for 2.2mm incision, triangular tips, delicate shaft, controlled opening, Ø 8.0mm, length 11.0cm

F-2982



GIANNETTI MIC Capsulorhexis Forceps, curved, for 2.0mm incision, controlled opening, Ø 8.0mm, length 11.0cm

F-2983

For Ceramic Coating please add "C"



Harms Tying Forceps (Tuebingen model), straight,
with 6.0mm tying platform,
length 10.0cm

F-2460



McPherson Suture Forceps, straight, with 5.0mm tying platform,
length 10.0cm

F-2470



McGregor Conjunctiva Fixation Forceps, with 5.0mm tying platform,
cross-serrated tips,
length 10.0cm

F-2480



Harms Tying Forceps (Tuebingen model), curved,
with 6.0mm tying platform,
length 10.0cm

F-2490



Harms Tying Forceps (Tuebingen model), angled,
with 6.0mm tying platform,
length 10.0cm

F-2500



Suture Forceps Kelman-McPherson, 7.5mm angled,
length 10.0cm

F-2510



Utrata Capsulorhexis Forceps, 12.0mm angled,
triangular shaped tips,
length 10.0cm

F-2520



For Ceramic Coating please add "C"



Iris Forceps Bonn model, straight, with 6.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 10.0cm

F-2540



Bonn Iris Forceps,
1x2 teeth, 0.12mm, 90° angle,
length 10.0 cm

F-2550



Castroviejo Suturing Forceps, straight, with 6.0mm tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 10.0cm

F-2560



Castroviejo Suturing Forceps, straight, with 6.0mm tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 10.0cm

F-2570



Castroviejo Suturing Forceps, straight, with 6.0mm tying platform,
1x2 teeth, 0.5mm, 45° angle,
length 10.0cm


F-2580



Castroviejo Suturing Forceps, straight, with tying platform,
1x2 teeth, 0.9mm, 45° angle,
length 10.0cm


F-2590

For Ceramic Coating please add "C"



Troutman Tying Forceps, straight, with extra fine 6.0mm tying platform, for 8-0 to 11-0 sutures, length 10.0cm

F-2800




Troutman Tying Forceps, curved, with extra fine 6.0mm tying platform, for 8-0 to 11-0 sutures, length 10.0cm

F-2810



Troutman Suturing Forceps, straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle, length 10.0cm

F-2820



Akahoshi Phaco Pre-Chopper, with 4.5mm delicate jaws & sharp tips, length 10.0cm

F-2850



Ernest-McDonald II IOL Inserting Forceps, 5.0mm straight jaws, length 10.0cm


F-2860

For Ceramic Coating please add "C"




Calf III IOL Forceps,
folder for foldable IOL's,
length 10.0cm

F-2830



Ernest-Calf Nucleus Cracking Forceps (cross-action),
extremely thin paddle,
length 10.0cm

F-2840



Livernois Pickup & Folding Forceps, angled,
L-shaped jaws with concave inside surfaces,
length 10.5cm

F-2880

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.0cm



Tennant Tying Forceps, straight, with 6.0mm platform, for 9.0-11.0 sutures

F-2890




Same in titanium

F-2890T




Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9.0-11.0 sutures

F-2900



Kelman-McPherson Tying Forceps, angled, with 8.0mm tying platform for 9.0-11.0 sutures

F-2910




Gaskin Fragment Forceps, angled, 11.0mm, ultra thin smooth jaws

F-2920

Same in titanium

F-2920T



Faulkner Lens Holding Forceps, 11.0mm, angled

F-2930




Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular tips

F-2940

Same in titanium

F-2940T



Jaffe Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips,

F-2960

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 8.0mm, length 11.0cm




Lehner Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips

F-2970



CORYDON Capsulorhexis Forceps, angled, 12.0mm, curved, cystotome shaped tips

F-2980



Bonn Forceps straight, with 6.0mm tying platform, 1x2 teeth, 0.12mm

F-2990

Same in titanium

F-2990T




Castroviejo Suture Forceps, straight, with tying platform, 1x2 teeth, 0.12mm 45° angle

F-3000

For Ceramic Coating please add "C"



With guide pin & stop pin, Ø 6.0mm, length 11.0cm




Tennant Tying Forceps straight, with 6.0mm tying platform, for 9.0-11.0 sutures

F-3100




Tennant Tying Forceps curved, with 6.0mm tying platform, for 9.0-11.0 sutures

F-3110



Kelman-McPherson Tying Forceps, angled, with 8.0mm tying platform, for 9.0-11.0 sutures

F-3120



Gaskin Fragment Forceps, angled, 11.0mm, with ultra thin smooth jaws

F-3130



Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips

F-3140



Jaffe Utrata Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips

F-3160



Bonn Forceps straight, with tying platform, 1x2 teeth, 0.12mm, 90° angle

F-3170

For Ceramic Coating please add "C"



Tennant Tying Forceps, straight, 6.0mm tying platform
for 9-0 to 1-0 sutures,
length 10.0cm

F-3210



Tennant Tying Forceps, curved, with 6.0mm tying platform,
for 9-0 to 1-0 sutures,
length 10.0cm

F-3220




Utrata Capsulorhexis Forceps, angled, 12.0mm,
triangular tips,
length 10.0cm

F-3250




Lehner Capsulorhexis Forceps, curved, 12.0mm,
triangular shaped tips,
length 10.0cm


F-3260



Colibri Forceps, straight, with tying platform,
1x2 teeth, 0.12mm, 45° angle,
length 10.0cm

F-3300






Tennant Tying Forceps, straight, with 6.0mm tying platform, for 9.0-1.0 sutures, length 10.5cm

F-3310




Tennant Tying Forceps, curved, with 6.0mm tying platform, for 9.0-1.0 sutures, length 10.5cm

F-3320




Utrata Capsulorhexis Forceps, angled, 12.0mm, triangular shaped tips, length 10.5cm

F-3330




Jaffe Capsulorhexis Forceps, angled, 12.0mm, blunt shaped tips, length 10.5cm

F-3331




O'Gawa Tying Forceps, straight, with 5.5mm tying platform, length 10.5cm

F-3340



O'Gawa Tying Forceps, curved, with 5.5mm tying platform, length 10.5cm

F-3350



Utrata Capsulorhexis Forceps, angled, 13.0mm, triangular shaped tips, length 10.5cm

F-3360

For Ceramic Coating please add "C"



Designed for the „Phaco Chopper Technique“



Salvitti Akahoshi Combo Pre-Chopper,
length 12.5cm

F-3390



Akahoshi Combo Pre-Chopper,
length 12.5cm

F-3400



Akahoshi Combo Pre-Chopper, Universal,
length 12.5cm


F-3410



Akahoshi Combo Pre-Chopper, cross-action,
length 12.5cm


F-3420

For Ceramic Coating please add "C"




Bishop Harmon Fixation Forceps, serrated, 0.5mm tips,
length 8.5cm

F-3530



Bishop Harmon Fixation Forceps, serrated, 0.8mm tips,
length 8.5cm

F-3540



Bishop Harmon Fixation Forceps, serrated, 1.0mm tips,
length 8.5cm

F-3550



Bishop-Harmon Dressing Forceps, criss-cross serrated, 0.5mm tips,
length 8.5cm

F-3560



Bishop-Harmon Dressing Forceps, criss-cross serrated, 0.8mm tips,
length 8.5cm

F-3570



Bishop-Harmon Forceps(Bonn model), with 5.0mm tying platform,
1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-3590



Bishop-Harmon Forceps (Bonn model), with 5.0mm tying platform,
1x2 teeth, 0.3mm, 90° angle,
length 8.5cm

F-3600

For Ceramic Coating please add "C"




Bishop-Harmon Forceps (Bonn model),
1x2 teeth, 0.4mm, 90° angle,
length 8.5cm

F-3610




Bishop-Harmon Forceps (Bonn model),
1x2 teeth, 0.7mm, 90° angle,
length 8.5cm

F-3620



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.12mm, 90° angle,
length 8.5cm

F-3630



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.3mm, 90° angle,
length 8.5cm

F-3640



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.5mm, 90° angle,
length 8.5cm

F-3650



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.6mm, 90° angle,
length 8.5cm

F-3660



Bishop-Harmon Tissue Forceps,
with 1x2 teeth, 0.8mm, 90°, angle,
length 8.5cm


F-3670

For Ceramic Coating please add "C"




Lester Conjunctiva Fixation Forceps, straight,
with 1x2 teeth, 0.6mm, 90° angle,
length 9.5cm

F-3690




Lester Conjunctiva Fixation Forceps, straight,
with 2x3 teeth, 0.6mm, 90° angle,
length 9.5cm

F-3700



Bonaccolto Utility Forceps, 1.2mm wide,
length 11.0cm

F-3720



Bonaccolto Utility Forceps, 1.7mm wide,
length 11.0cm

F-3730



Elschnig Skleral Fixation Forceps, straight,
with 1x2 teeth, 0.6mm, 45° angle,
length 11.0cm


F-3750



Elschnig Skleral Fixation Forceps, straight,
1x2 teeth, 0.8mm, 45° angle,
length 11.0cm

F-3760


For Ceramic Coating please add "C"



Remky Superior Rectus Forceps,
length 11.0cm F-3780



Lister Conjunctiva Forceps,
length 11.0cm F-3810



Bonaccolto Forceps (French model), straight,
serrated jaws, 0.8mm,
length 11.0cm F-3820




Jayle Forceps, straight, with 7.0mm tying platform,
1x2 teeth, 0.3mm, 45° angle,
length 11.0cm F-3840



Jayle Forceps, straight, with 7.0mm tying platform,
1x2 teeth, 0.8mm, 45° angle,
length 11.0cm F-3860

For Ceramic Coating please add "C"




Dressing Forceps, serrated, straight, 0.5mm,
length 10.0cm

F-3870




Dressing Forceps, serrated, straight, 0.8mm,
length 10.0cm

F-3880




Dressing Forceps, serrated, straight, 1.0mm,
length 10.0cm

F-3890




Dressing Forceps, serrated, straight, 1.2mm,
length 10.0cm

F-3900




Dressing Forceps, serrated, curved, 0.5mm,
length 10.0cm

F-3910



Dressing Forceps, serrated, strong curved, 0.6mm,
length 10.0cm


F-3920



Iris Forceps serrated, curved, 0.8mm,
length 10.0cm


F-3930

For Ceramic Coating please add "C"




Dressing Forceps, serrated, strong curved, 0.8mm, length 10.0cm

F-3940




Dressing Forceps, serrated, curved, 1.0mm, length 10.0cm

F-3950




Dressing Forceps, serrated, strong curved, 1.0mm, length 10.0cm

F-3960



Dressing Forceps, serrated, curved, 1.2 mm, length 10.0cm

F-3970



Nugent Utility Forceps, curved shafts, 45° angle, with smooth 10.0mm jaws, length 10.0cm

F-3980




Iris Forceps, straight, with 1x2 teeth, 0.7mm, 90° angle, length 10.0cm

F-4000



Iris Forceps, straight, with 1x2 teeth, 1.0mm, 90° angle, length 10.0cm


F-4010



Iris Forceps, curved, with 1x2 teeth, 0.8mm, 90° angle, length 10.0cm


F-4020

For Ceramic Coating please add "C"




Barraquer Fixation Forceps,
length 10.0cm

F-4040



Graefe Fixation Forceps, without lock,
length 11.0cm

F-4050



Graefe Fixation Forceps, with lock,
length 11.0cm

F-4060



Barraquer Cilia Forceps,
length 10.5cm

F-4070

For Ceramic Coating please add "C"



Barraquer Cilia Forceps, slightly angled with 6.0mm tying platform,
length 11.0cm

F-1980



Cilia Forceps, 2.0mm,
length 9.0cm

F-4080

Cilia Forceps, 3.0mm,
length 9.0cm

F-4090




Gradle Cilia Forceps,
smooth jaws,
length 9.5cm

F-4100

Gradle Cilia Forceps,
with positive gripping jaws,
length 9.5cm


F-4110

For Ceramic Coating please add "C"



Jeweler's Forceps #5, straight,
extra delicate

F-4120



Jeweler's Forceps #4, straight,
extra delicate pointed tips


F-4130



Jeweler's Forceps #7, curved,
extra delicate, pointed,
length 11.5 cm


F-4140

For Ceramic Coating please add "C"



Moorfields Conjunctiva Forceps,
standard model,
length 11.0cm

F-4170



Moorfield Conjunctiva Forceps,
fine model,
length 11.0cm

F-4180



Watzke Sleeve Spreader,
length 11.0cm

F-4200

For Ceramic Coating please add "C"



Shaaf Cilia Foreign Body Forceps,
length 9.5cm

F-4190



Ayer Chalazion Forceps,
round, Ø 8.0mm,
length 9.2cm

F-4210




Lambert Chalazion Forceps,
round, tip width: 11.0mm,
length 9.2cm

F-4220

Lambert Chalazion Forceps,
round, tip width: 15.0mm,
length 9.2cm

F-4230

For Ceramic Coating please add "C"



Desmarres Chalazion Forceps, 20.0mm,
length 9.2cm

F-4240

Desmarres Chalazion Forceps, 26.0mm,
length 9.2cm

F-4250

Desmarres Chalazion Forceps, 31.0mm,
length 9.2cm

F-4260



Hunt Chalazion Forceps,
round, Ø 12.0mm,
length 9.7cm


F-4270

For Ceramic Coating please add "C"



Baird Chalazion Forceps, oval,
11.0 x 8.0mm, small,
length 9.5cm

F-4280



Francis Chalazion Forceps oval,
14.0 x 11.0mm I.D., medium,
length 9.5cm

F-4290



Wies Chalazion Forceps,
length 9.5cm

F-4300

For Ceramic Coating please add "C"



Desmarres cross-action Forceps, oval,
14.0mm x 11.0mm,
Length 9.5cm

F-4310



Lordan Chalazion Forceps,
length 9.5cm

F-4320



Berke Ptosis Forceps,
with slide lock, longitudinally grooved jaws,
20.0mm, length 10.5cm

F-4330

Berke Ptosis Forceps,
with slide lock, longitudinally grooved jaws,
27.0mm, length 10.5cm

F-4340

For Ceramic Coating please add "C"



Snellen Entropium Forceps, right,
with small solid blades, 28.0mm,
length 9.5cm

F-4350

Snellen Entropium Forceps, left,
with small solid blades, 28.0mm,
length 9.5cm

F-4360

Snellen Entropium Forceps, right,
with solid blades, 35.0mm,
length 9.5cm

F-4370

Snellen Entropium Forceps, left,
with solid blades, 35.0mm,
length 9.5cm

F-4380



Downs Lid Forceps,
length 9.8cm

F-4395

For Ceramic Coating please add "C"



Putterman Forceps, Lid Clamp,
with 6 pins, serrated jaws, sliding lock,
length 10.0cm

F-4390



Jameson Muscle Forceps, right,
with slide lock, 4 teeth & holes,
length 10.0cm

F-4400

Jameson Muscle Forceps, left,
with slide lock, 4 teeth & holes,
length 10.0cm

F-4405

Jameson Muscle Forceps, right,
with slide lock, 6 teeth & holes,
length 10.0cm

F-4410

Jameson Muscle Forceps, left,
with slide lock, 6 teeth & holes,
length 10.0cm

F-4415

For Ceramic Coating please add "C"



Ocular Landers Lens Forceps,
length 10,3cm

F-4420



Scleral Plug Forceps,
length 10.3cm

F-4430

Designed to hold scleral plugs during insertion and removal,
cross-action handle



Bangerter Muscle Clamp,
length 10.5cm

F-4440

For Ceramic Coating please add "C"



Dieffenbach Bulldog Clamp, serrifine straight,
length 3.5cm

F-4450

Dieffenbach Bulldog Clamp, serrifine curved,
length 3.5cm

F-4460



Dieffenbach Bulldog Clamp, serrifine straight,
length 5.0cm

F-4470

Dieffenbach Bulldog Clamp, serrifine straight, with hollow rivet,
length 3.5cm


F-4485

For Ceramic Coating please add "C"




Blaydes Implantation Forceps, without lock,
very delicate narrow jaws angled 45°,
length 11.0cm

F-4490



Shepard Lens Forceps, without lock,
length 11.0cm



F-4540



Inamura Capsulorhexis cross-action Forceps, titanium,
13.0mm curved, triangular shaped tips,
length 10.5cm


F-4575T

For Ceramic Coating please add "C"



Blaydes Implantations Forceps,
without lock, very delicate narrow jaws, angled 45°,
length 11.0cm

F-4600



Shepard Lens Holding Forceps,
without lock, gently curved,
length 11.0cm

F-4610



Hartmann Mosquito Forceps, straight,
length 9.5cm

F-4640

Hartmann Mosquito Forceps, curved,
length 9.5cm

F-4650

Halsted Mosquito Forceps, straight,
length 14.5cm

F-4660

Halsted Mosquito Forceps, curved,
length 14.5cm

F-4670



Baby Jones Towel Clamp,
length 5.5cm

F-4680

For Ceramic Coating please add "C"



Backhaus towel clamp,
length 9.0cm

F-4690

Backhaus towel clamp,
length 13.5cm

F-4700



Bipolar Forceps straight,
slightly blunted tips, 0.25mm,
length 11.0cm

F-4710

Bipolar Forceps angled,
slightly blunted tips, 0.25mm,
length 11.0cm

F-4720

Bipolar Cord, autoclavable,
5.0m

F-4730

F-4710, F-4720, F-4730:



Günter Bissinger Medizintechnik GmbH
Teningen, Germany

For Ceramic Coating please add "C"



McPherson Bipolar Forceps, straight,
slightly blunted tips, 0.25mm, isolated,
length 9.5cm

F-4740

McPherson Bipolar Forceps, angled,
slightly blunted tips, uninsulated,
length 9.5cm

F-4750

Bipolar Cord, autoclavable,
length 5.0m

F-4730

F-4730, F-4740, F-4750:



Günter Bissinger Medizintechnik GmbH
Teningen, Germany

For Ceramic Coating please add "C"

NEEDLE HOLDERS ERGOGRIP NEEDLE HOLDERS BLADE BREAKERS TITANIUM NEEDLE HOLDERS TC NEEDLE HOLDERS



ULTRA LIGHT ERGO GRIP

...Perfect balance for doctors' use

Every surgical operation requires skill and good sensitivity. Operations with longer instruments under the microscope in particular require thoroughness and precision work. Long operations are especially challenging for the surgeon. In order to facilitate the work of the surgeon, a1 medical has developed an ergonomically shaped and light-weight line of instruments.

The "Ultra-Light Ergo Grip". At a1 medical, our micro instruments are manufactured only from the finest stainless steel; their protruding ribs give you a perfect grip, and they can be easily cleaned. Their ergonomic design is not only pleasing to the eye but also offers perfect balance in the hand. Particularly noteworthy is the price-performance ratio as compared to titanium. Instruments made of titanium usually cost almost twice as much as stainless steel instruments.

Due to their excellent design, "Ultra-Light Ergo Grip" instruments weigh hardly any more than those made of titanium and are thus a good cost alternative.

Of course, our "Ultra-Light Ergoline" is also available with features such as TC coating and ceramic coating.

ADVANTAGES AT A GLANCE

- provides perfect balance during surgery
- ergonomic & attractive design
- excellent grip
- high quality stainless steel
- easy to maintain & sterilize
- saves costs



For Ceramic Coating please add "C"



Ergo Grip Barraquer Needle Holder, straight, 10.0 x 1.0mm, N-0800
smooth jaw, with lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, straight, 10.0 x 1.0mm, N-0810
smooth jaw, without lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 1.0mm, N-0820
smooth jaw, with lock, length 13.5cm

Ergo Grip Barraquer Needle Holder, curved, 10.0 x 1.0mm, N-0830
smooth jaw, without lock, length 13.5cm



Ergo Grip Barraquer Needle Holder, curved, 10.0 x 0.6mm,
smooth jaw, without lock, length 13.5cm

N-0870

For Ceramic Coating please add "C"

TC-COATING

TC-Coating



When coated with tungsten carbide, tips of instruments last significantly longer since tungsten carbide is much harder than stainless steel. This allows a save handling of fine needles.



Ergo Grip Barraquer Needle Holder with TC-Coating, straight,
10.0 x 1.0mm, with lock,
length 13.5cm

N-0800TC




Ergo Grip Barraquer Needle Holder, with TC-Coating, straight,
10.0 x 1.0mm, without lock,
length 13.5cm

N-0810TC



Ergo Grip Barraquer Needle Holder with TC-Coating, curved,
10.0 x 1.0mm, without lock,
length 13.5cm

N-0830TC




McPherson Needle Holder, straight, 10.0mm x 0.5mm, N-0010
smooth jaw, with lock, 11.0cm

McPherson Needle Holder, straight, 10.0mm x 0.5mm, N-0020
smooth jaw, without lock, 11.0cm

McPherson Needle Holder, curved, 10.0mm x 0.5mm, N-0030
smooth jaw, with lock, 11.0cm

McPherson Needle Holder, curved, 10.0mm x 0.5mm, N-0040
smooth jaw, without lock, 11.0cm



McPherson Needle Holder, straight, 10.0mm x 0.5mm, N-0050
smooth jaw, with lock, 10.0cm

McPherson Needle Holder, straight, 10.0mm x 0.5mm, N-0060
smooth jaw, without lock, 10.0cm

McPherson Needle Holder, curved, 10.0mm x 0.5mm, N-0070
smooth jaw, with lock, 10.0cm

McPherson Needle Holder, curved, 10.0mm x 0.5mm, N-0080
smooth jaw, without lock, 10.0cm

For Ceramic Coating please add "C"



Castroviejo Needle Holder, straight, 10.0mm x 1.0mm,
smooth jaw, with lock, length 14.0cm N-0090

Castroviejo Needle Holder, straight, 10.0mm x 1.0mm,
smooth jaw, without lock, length 14.0cm N-0100

Castroviejo Needle Holder, curved, 10.0mm x 1.0mm,
smooth jaw, with lock, length 14.0cm N-0110

Castroviejo Needle Holder, curved, 10.0mm x 1.0mm,
smooth jaw, without lock, length 14.0cm N-0120

Castroviejo Needle Holder, straight, 10.0mm x 0.75mm,
smooth jaw, with lock, length 14.0cm N-0130

Castroviejo Needle Holder, straight, 10.0mm x 0.75mm,
smooth jaw, without lock, length 14.0cm N-0140

Castroviejo Needle Holder, curved, 10.0mm x 0.75mm,
smooth jaw, with lock, length 14.0cm N-0150

Castroviejo Needle Holder, curved, 10.0mm x 0.75mm,
smooth jaw, without lock, length 14.0cm N-0160

For Ceramic Coating please add "C"



Kalt Needle Holder, straight, 10.0 x 1.0mm,
length 14.0cm

N-0250



Barraquer Needle Holder (Baby Model), straight, 9.0 x 0.6mm,
smooth jaw, with lock, length 10.0cm N-0280

Barraquer Needle Holder (Baby Model), straight, 9.0 x 0.6mm,
smooth jaw, without lock, length 10.0cm N-0290


Barraquer Needle Holder (Baby Model), curved, 9.0 x 0.6mm,
smooth jaw, with lock, length 10.0cm N-0300

Barraquer Needle Holder (Baby Model), curved, 9.0 x 0.6mm,
smooth jaw, without lock, length 10.0cm N-0310

Barraquer Needle Holder, curved, 9.0 x 0.5mm,
smooth jaw, with lock, length 12.0cm N-0320

Barraquer Needle Holder, curved, 9.0 x 0.5mm,
smooth jaw, without lock, length 12.0cm N-0330

For Ceramic Coating please add "C"



Barraquer Needle Holder, curved, 9.0 x 0.75mm, smooth jaw, with lock, length 12.0cm	N-0340
--	--------


Barraquer Needle Holder, curved, 9.0 x 0.75mm, smooth jaw, without lock, length 12.0cm	N-0350
---	--------

Barraquer Needle Holder, straight, 10.0 x 1.0mm, smooth jaw, with lock, length 13.5cm	N-0360
--	--------

Barraquer Needle Holder, straight, 10.0 x 1.0mm, smooth jaw, without lock, length 13.5cm	N-0370
---	--------

Barraquer Needle Holder, curved, 10.0 x 1.0mm, smooth jaw, with lock, length 13.5cm	N-0380
--	--------

Barraquer Needle Holder curved, 10.0 x 1.0mm, smooth jaw, without lock, length 13.5cm	N-0390
--	--------



Barraquer Needle Holder, straight, 10.0 x 0.6mm,
with lock, length 13.5cm

N-0400

Barraquer Needle Holder, straight, 10.0 x 0.6mm,
without lock, length 13.5cm

N-0410

Barraquer Needle Holder, curved, 10.0 x 0.6mm,
smooth jaw, with lock, length 13.5cm

N-0420

Barraquer Needle Holder, curved, 10.0 x 0.6mm,
smooth jaw, without lock, length 13.5cm

N-0430

Barraquer Needle Holder, curved, 10.0 x 0.5mm,
smooth jaw, with lock, length 13.5cm


N-0440

Barraquer Needle Holder, curved, 10.0 x 0.5mm,
smooth jaw, without lock, length 13.5cm

N-0450

Needle Holder in titanium upon request!

For Ceramic Coating please add "C"



Anis Needle Holder, straight, 7.0 x 0.4mm,
smooth jaw, with lock, length 12.0cm N-0460

Anis Needle Holder, straight, 7.0 x 0.4mm,
smooth jaw, without lock, length 12.0cm N-0470

Anis Needle Holder, curved, 7.0 x 0.4mm,
smooth jaw, with lock, length 12.0cm N-0480

Anis Needle Holder, curved, 7.0 x 0.4mm,
smooth jaw, without lock, length 12.0cm N-0490



Jakobi Needle Holder, straight, 7.0 x 0.5mm,
smooth jaw, with lock, length 12.0cm

N-0500

Jakobi Needle Holder, straight, 7.0 x 0.5mm,
smooth jaw, without lock, length 12.0cm

N-0510

Jakobi Needle Holder, curved, 7.0 x 0.5mm,
smooth jaw, with lock, length 12.0cm

N-0520

Jakobi Needle Holder, curved, 7.0 x 0.5mm,
smooth jaw, without lock, length 12.0cm

N-0530



Barraquer Needle Holder, French style, straight, 10.0 x 0.6mm
smooth jaw, without lock, length 12.0cm N-0540

Barraquer Needle Holder, French style, straight, 10.0 x 0.6mm
smooth jaw, with lock, length 12.0cm N-0550

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth jaw, without lock, length 12.0cm N-0560

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth jaw, with lock, length 12.0cm N-0570

Barraquer Needle Holder, French style, curved, 10.0 x 0.6mm
smooth, without lock, length 13.5cm N-0630



Castroviejo Blade Breaker and Holder,
with sliding lock,
length 13.5cm

N-0660



Troutman-Chris Blade Breaker & Holder,
with concave/ convex inside surface,
length 9.5cm

*N-0690



Derf Needle Holder straight,
serrated jaws with groove,
length 12.5cm

N-0700

For Ceramic Coating please add "C"

* Not available in Ceramic Coating

SCISSORS

SUPER-CUT SCISSORS

RING SCISSORS

MICRO SCISSORS



GENERAL INFORMATION




Sharp scissors (Adventitia)

For fine trimming and secure suture cutting.



Blunt scissors (Dissection)

Blunt blades avoid possible damages during dissection



Super Cut Scissors

A razor sharp cutting edge on one blade and a micro-serration on the opposite blade assure a perfectly precise and non traumatic cut with minimal force.



Vannas blade

The vannas scissors are also applied in vascular surgery because due to their more delicate, smaller and shorter blades.



Castroviejo Universal Corneal Scissors, curved,
blunt tips & small blades, cutting length 5.0mm,
length 11.5cm

S-0010

Castroviejo Universal Corneal Scissors, curved,
blunt tips & medium blades, cutting length 6.5mm,
length 11.5cm

S-0020

Castroviejo Universal Corneal Scissors, curved,
blunt tips & large blades, cutting length 8.0mm,
length 11.5cm

S-0030



Westcott Tenotomy Suture Scissors, curved,
blunt tips & medium blades, cutting length 9.0mm,
length 11.5cm

S-0040

McPherson-Westcott Suture Scissors, curved,
sharp tips & small blades, cutting length 6.0mm,
length 11.5cm

S-0050

McPherson-Westcott Suture Scissors, curved,
blunt tips & small blades, cutting length 6.0mm,
length 11.5cm

S-0060

For Ceramic Coating please add "C"



Jaffe Stitch Scissors, curved,
sharp tips & medium blades, cutting length 8.0mm,
length 11.5cm

S-0070



Fine Stitch Scissors, straight,
sharp tips & medium blades, cutting length 11.0mm,
length 11.5cm

S-0080

Fine Stitch Scissors, curved,
sharp tips & medium blades, cutting length 11.0mm,
length 11.5cm

S-0090

For Ceramic Coating please add "C"



Vannas Scissors, straight,
sharp tips & 10.0mm blades, cutting length 8.0mm,
length 11.5cm

S-0100

Vannas Scissors, curved,
sharp tips & 10.0mm blades, cutting length 8.0mm,
length 11.5cm

S-0110



Castroviejo Corneoscleral Scissors, right,
small blades, cutting length 5.0mm,
length 10.0cm

S-0120

Castroviejo Corneoscleral Scissors, left,
small blades, cutting length 5.0mm,
length 10.0cm

S-0130

Castroviejo Corneoscleral Scissors, right,
medium blades, cutting length 7.0mm,
length 10.0cm

S-0140

Castroviejo Corneoscleral Scissors, left,
medium blades, cutting length 7.0mm,
length 10.0cm

S-0150

For Ceramic Coating please add "C"



Castroviejo Universal Corneal Scissors, curved,
blunt tips & small blades, cutting length 6.0mm,
length 10.0cm

S-0160

Castroviejo Universal Corneal Scissors, curved,
blunt tips & medium blades, cutting length 10.0mm,
length 10.0cm

S-0170

Castroviejo Universal Corneal Scissors, curved,
blunt tips & large blades, cutting length 10.0mm,
length 10.0cm

S-0180



Castroviejo Keratoplasty Scissors, angled to side,
blunt tips & medium blades, cutting length 10.0mm,
length 10.0cm

S-0190



McPherson-Westcott Conjunctival Scissors, curved, blunt tips & small blades, cutting length 5.0mm, length 10.5cm	S-0200
Westcott Tenotomy Scissors, curved, blunt tips & medium blades, cutting length 8.0mm, length 11.5cm	S-0210
Westcott Tenotomy Scissors, curved right, blunt tips & standard blades, cutting length 11.0mm, length 11.5cm	S-0382
Westcott Tenotomy Scissors, curved right, blunt tips & standard blades, cutting length 14.0mm, length 11.5cm	S-0220
Westcott Tenotomy Scissors, curved left, blunt tips & standard blades, cutting length 9.0mm, length 11.5cm	S-0230
Westcott Tenotomy Scissors, curved, blunt tips & small blades, cutting length 5.0mm, length 10.0cm	S-0240
Westcott Tenotomy Scissors, curved, pointed tips & standard blades, cutting length 14.0mm, length 11.5cm	S-0250

For Ceramic Coating please add "C"



Iris Scissors, straight, sharp tips & small blades
cutting length 6.0mm,
length 10.0cm

S-0260

Iris Scissors, curved, sharp tips & small blades,
cutting length 6.0mm,
length 10.0cm

S-0270



Noyes Iris Scissors, straight, sharp tips & small blades,
cutting length 10.0mm,
length 11.0cm

S-0280

Noyes Iris Scissors, lightly curved, sharp tips & small blades,
cutting length 10.0mm,
length 11.0cm

S-0290

Noyes Iris Scissors, straight, sharp tips & standard blades,
cutting length 12.5mm,
length 11.0cm

S-0300

Noyes Iris Scissors, lightly curved, sharp tips & standard blades,
cutting length 12.5mm,
length 11.0cm

S-0310

For Ceramic Coating please add "C"



Westcott Tenotomy Scissors, curved,
blunt tips & 13.0mm blades, cutting length 7.5mm,
length 14.0cm

S-0380

Westcott Tenotomy Scissors, curved,
blunt tips & 22.0mm blades, cutting length 11.0mm,
length 14.0cm

S-0390



Troutman-Katzin Corneal Transplant Scissors, strong curved, right,
pointed tips & small blades, cutting length 6.0mm,
length 10.0cm

S-0400

Troutman-Katzin Corneal Transplant Scissors, strong curved, left,
pointed tips & small blades, cutting length 6.0mm
length 10.0cm

S-0410

For Ceramic Coating please add "C"



Vannas Capsulotomy Scissors, straight,
sharp tips, 5.0mm blade,
length 8.5cm

S-0420

Vannas Capsulotomy Scissors, straight,
blunt tips, 5.0mm blade,
length 8.5cm

S-0430

Vannas Capsulotomy Scissors, curved,
sharp tips, 5.0mm blade,
length 8.5cm

S-0440

Vannas Capsulotomy Scissors, curved,
blunt tips, 5.0mm blade,
length 8.5cm

S-0450

Gills-Vannas Capsulotomy Scissors, straight,
sharp tips, 7.0mm blade,
length 8.5cm

S-0460

Gills-Vannas Capsulotomy Scissors, curved,
sharp tips, 7.0mm blade,
length 8.5cm

S-0470

Vannas Capsulotomy Scissors, angled forward,
sharp tips, 5.0mm blade,
length 8.5cm

S-0480

For Ceramic Coating please add "C"



Gills-Vannas Capsulotomy Scissors, angled forward,
sharp tips, 7.0mm blade,
length 8.5cm

S-0490

Gills-Welsh Capsulotomy Scissors, angled forward,
sharps tips, 11.0mm blade,
length 8.5cm

S-0501

Uribe-Stern Capsulotomy Scissors, angled forward,
concave curved blade,
length 8.5cm

S-0505



DeWecker Iridectomy Scissors, angled,
sharp tips & 7.0mm blade,
length 11.0cm

S-0510

For Ceramic Coating please add "C"



Vannas Scissors, sharp, curved, cutting length 6.0mm,
Ø 6.0mm handle, length 11.5cm

S-0415

Vannas Scissors, blunt, curved, cutting length 6.0mm,
Ø 6.0mm handle, length 11.5cm

S-0416



Barraquer Iris Scissors,
blunt tips & 7.0mm blades

S-0530

Barraquer Iris Scissors,
pointed tips & 7.0mm blades

S-0540



Eye Scissors, straight, pointed tips, small blades,
14.0mm blades,
length 9.0cm

S-0570

Eye Scissors, curved, pointed tips, small blades,
14.0mm blades,
length 9.0cm

S-0580

Eye Scissors, curved, blunt tips, small blades,
14.0mm blades,
length 9.0cm

S-0581

Eye Scissors, straight, pointed tips, large blades,
23.0mm blades,
length 9.0cm

S-0590

Eye Scissors, curved, pointed tips, large blades,
24.0mm blades,
length 9.0cm

S-0600

Eye Scissors, straight, pointed tips, standard model,
23.0mm blades,
length 10.5cm

S-0610

Eye Scissors, straight, blunt tips, standard model,
23.0mm blades,
length 10.5cm

S-0611

Eye Scissors, curved, pointed tips, standard model,
24.0mm blades,
length 10.5cm

S-0620

For Ceramic Coating please add "C"



Eye Scissors, curved, blunt tips, standard model,
24.0mm blades,
length 10.5cm

S-0621

Eye Scissors, straight, pointed tips, standard model,
23.5mm blades,
length 9.5cm

S-0630

Eye Scissors, curved, pointed tips, big rings,
22.5mm blades,
length 9.5cm

S-0640

Eye Scissors, curved, blunt tips, big rings,
22,5mm blades,
length 9.5cm

S-0641



Stitch Scissors, ribbon type, curved,
pointed tips, 12.50mm blades,
length 11.0cm

S-0650

For Ceramic Coating please add "C"



Stitch Scissors, lightly curved,
pointed tips, 12.5mm blades,
length 9.5cm

S-0660




Stevens Tenotomy Scissors, straight,
blunt tips, 14.0mm blades,
length 11.5cm

S-0670

Stevens Tenotomy Scissors, curved,
blunt tips, 13.0mm blades,
length 11.5cm

S-0680

For Ceramic Coating please add "C"



Stevens Tenotomy Scissors, ribbon type, straight,
blunt tips, 13.0mm blades,
length 9.5cm

S-0690

Stevens Tenotomy Scissors, ribbon type, curved,
blunt tips, 12.0mm blades,
length 9.5cm

S-0700



Strabismus Scissors, ribbon type, straight,
blunt tips, 30.0mm blades,
length 10.5cm

S-0710

Strabismus Scissors, ribbon type, curved,
blunt tips, 29.0mm blades,
length 10.5cm

S-0720

For Ceramic Coating please add "C"



Knapp Strabismus Scissors, straight,
blunt-blunt tips, 31.0mm blades,
length 10.3cm

S-0730

Knapp Strabismus Scissors, curved,
blunt-blunt tips, 30.0mm blades,
length 10.3cm

S-0740



Enucleation Scissors, medium curved,
blunt tips, 43.5mm blades,
length 13.5cm

S-0750

Enucleation Scissors, strong curved,
blunt tips, 40.2mm blades,
length 13.5cm

S-0760

For Ceramic Coating please add "C"



Utility Scissors, straight,
sharp/ blunt tips,
length 13.5cm

S-0770



Metzenbaum Scissors, curved,
blunt tips, 48.0mm blades,
length 14.5cm

S-0780



Spencer Ligature Scissors, straight,
with one notched tip, 21.0mm blades,
length 9.0cm

S-0790

For Ceramic Coating please add "C"



Super-Cut Scissors

Perfectly precise and non-traumatic cutting through any kind of tissue with minimal force requires a SuperCut Scissor. SuperCut Scissors have

- a razor sharp cutting edge on one blade
- and a standard edge combined with a micro-serration on the opposite blade.

The micro-serration prevents tissue from slipping away during the cutting process.

Super-Cut Stevens Tentotomy Scissors, straight,
blunt tips, 14.0mm blades,
length 11.5cm

S-0670SC

Super-Cut Stevens Tentotomy Scissors, curved,
blunt tips, 13.0mm blades,
length 11.5cm

S-0680SC



Super-Cut Metzenbaum Scissors, curved,
blunt tips, 48.0mm, length 14.5cm

S-0780SC

For Ceramic Coating please add "C"



Super-Cut Scissors

Perfectly precise and non-traumatic cutting through any kind of tissue with minimal force requires a SuperCut Scissor. SuperCut Scissors have

- a razor sharp cutting edge on one blade
- and a standard edge combined with a micro-serration on the opposite blade.

The micro-serration prevents tissue from slipping away during the cutting process.

Super-Cut Westcott Tenotomy Scissors, curved,
blunt tips, standard blade, cutting length 14.0mm,
length 11.5cm

S-0220SC

Super-Cut Westcott Tenotomy Scissors, curved,
pointed tips, standard blades, cutting length 14.0mm,
length 11.5cm

S-0250SC

For Ceramic Coating please add "C"

VITRECTOMY INSTRUMENTS



Titanium handle & stainless steel tip



Vitrectomy Forceps, smooth,
with grip straight

23G

VR-0010-23T



Vitrectomy Forceps, serrated tip,
straight

23G

VR-0030-23T



Vitrectomy Forceps, serrated tip,
straight

23G

VR-0050-23T



Vitrectomy Forceps,
with notch, 0.3mm, straight

23G

VR-0090-23T



Vitrectomy Forceps,
grasping tip

23G

VR-0110-23T



Endgripping Forceps, straight,
extra fine, 0.2mm x 0.5mm

23G

VR-0205-23T



Titanium handle & stainless steel tip



Vitrectomy scissors,
straight

23G

VR-0220 -23T



Vitrectomy Scissors,
angled 60°

23G


VR-0240-23T



Vitrectomy Scissors,
with scoop tips, 60°


23G

VR-0250-23T



Vitrectomy Foreign Body Forceps, 20G,
three 1.9mm long arms, slightly angled inwards,
titanium handle & stainless steel tip

VR-0270



Vitrectomy Foreign Body Forceps, 20G,
titanium handle & stainless steel tip

VR-0280



Lens cutting scissors with 4mm stainless steel blades to cut and explant silicone or acrylic lenses through a 1.8mm incision. Both blades have small serrations inside to prevent slippage of lens during cutting and removal procedure.



Lens cutting Scissors, straight blades, 15mm shaft, 18G

VR-1600

Lens cutting Scissors, curved blades, 15mm shaft, 18G

VR-1610



Lens holding forceps for removal of acrylic or silicone IOL's, with multiple grooved jaws, for securely gripping lens during explantation procedure



Grasping Forceps, 20G

VR-1620



Detachable System

Choose the instrument head you require and attach to your handle. Each head will be supplied in a protection cap. This guarantees an easier cleaning & storage.

All instruments are suitable to use with 23G trocar system.

Titanium handle & stainless steel tip

Titanium Handle only,
used with tip variations.



Micro Gripping Forceps, serrated,
1.5mm tip, 30.0mm shaft. Excellent
for peeling and manipulating thicker
membrane.

Tip for VR-0300-23TDT



ILM Forceps, 2.0mm tip length, the
jaw surface maximizes the gripping
ability during ILM peeling.

Tip for VR-0310-23TDT



Straight Scissors, sharp, tip length
2.5mm. Ideal for cutting in the iris
plane and peripheral membranes.

Tip for VR-0400-23TDT



Curved Scissors, sharp, tip length
1.25mm. Ideal for cutting parallel to
the retina and cutting membranes
adherent to the retina surface.

Tip for VR-0410-23TDT

Handle only

VR-5000T

23G

Complete instrument

VR-0300-23TDT*

Tip only

VR-0300-23-TIP

23G

Complete instrument

VR-0310-23TDT*

Tip only

VR-0310-23-TIP

23G

Complete instrument

VR-0400-23TDT*

Tip only

VR-0400-23-TIP

23G

Complete instrument

VR-0410-23TDT*

Tip only

VR-0410-23-TIP

*DT= Detachable

CANNULAS

LACRIMAL INSTRUMENTS

BIMANUAL I/A INSTRUMENTS

DETACHABLE I/A INSTRUMENTS

SEMI-DETACHABLE PUNCHES

SEMI-DETACHABLE PUNCH

**a1 medical presents the shortest semi-detachable Punch,
on the market, for the conventional lacrimal sac surgery**

The Punches have a working length of 11.0cm, which is specially designed and adapted to the conventional lacrimal sac surgery, a punching angle of 90° and an ejector built into the shaft.



Ferris-Smith-Kerrison Punch, with semi-detachable shaft,
1.0mm, 90°, with ejector, length: 11.0cm

L-0001

Ferris-Smith-Kerrison Punch, with semi-detachable shaft,
2.0mm, 90°, with ejector, length: 11.0cm

L-0002

For Ceramic Coating please add "C"

- ① Semi-Detachable shaft for a better and optimal maintenance of the instrument. All angles and lumen are easily reachable

Semi-detachable means, that the punch does not and can not be disassembled completely. No parts could get lost or mixed up during and after maintenance process and when reassembling the punch

Easy open and close system by pushing the button of the release mechanism. Very fast training- and learning-curve

Specially sharpened Punching tips. Less force for more efficiency



Scan this QR-Code for watching the video

about disassembling and reassembling the semi-detachable Punch



Kerrison Bone Nibbling Rongeur,
size # 0, 3.0mm

L-0010

Kerrison Bone Nibbling Rongeur,
size # 1, 3.5mm

L-0020

Kerrison Bone Nibbling Rongeur,
size # 2, 4.0mm

L-0030

Kerrison Bone Nibbling Rongeur,
size # 3, 4.5mm

L-0040



Citelli Bone Nibbling Rongeur,
1.5mm wide

L-0050



Beyer Rongeur, single action, curved

L-0080



Nasal Speculum, infant L-0090

Nasal Speculum, child L-0100

Nasal Speculum, adult L-0110



Castroviejo Lacrimal Dilator, double-ended, fine & medium taper, length 14.0cm L-0120

For Ceramic Coating please add "C"



Wilder Lacrimal Dilator,
size # 1, fine taper

L-0130

Wilder Lacrimal Dilator,
size # 2, medium taper

L-0140

Wilder Lacrimal Dilator,
size # 3, heavy taper

L-0145

Ruedemann Infant Dilator,
length 8.0cm

L-0150



Freer Periosteal Elevator,
length 18.0cm

L-0160

For Ceramic Coating please add "C"



Partsch Mallet,
stainless steel head,
length 18.0cm

L-0170



Joseph Periosteal Elevator,
gently curved tip

L-0180

For Ceramic Coating please add "C"



Partsch Lacrimal Sac Chisel,
length 13.5cm

3.0mm wide

L-0200

5.0mm wide

L-0220

For Ceramic Coating please add "C"



West Lacrimal Chisel,
straight

L-0290

For Ceramic Coating please add "C"



Freer Lacrimal Chisel,
curved, 6.0mm,
length 16.0cm

L-0300



Harms Trabeculotomy Probe,
9.0mm long pointed tips, right

L-0310

Harms Trabeculotomy Probe,
9.0mm long pointed tips, left

L-0320

For Ceramic Coating please add "C"



Bowman Lacrimal Probe, Stainless steel, set of 6		L-0330
	ø mm	
size 0000 & 000	0.4 / 0.6	L-0340
size 00 & 0	0.7 / 0.8	L-0350
size 1 & 2	0.9 / 1.1	L-0360
size 3 & 4	1.3 / 1.4	L-0370
size 5 & 6	1.5 / 1.6	L-0380
size 7 & 8	1.8 / 1.9	L-0390

For Ceramic Coating please add "C"



Kellnar Retrograde Probe,
right

*L-0400



Catalano Lacrimal Intubations Set,
0.7mm x 14.0cm

*L-0420

Silicone tubing

*L-0430

* Not available in Ceramic Coating



Crawford Lacrimal Intubation Set,
0.6mm x 14.0cm

*L-0440

Silicone tubing

*L-0450



Crawford Retrieving Hook,
length 15.0cm

L-0460

For Ceramic Coating please add "C"
* Not available in Ceramic Coating



Pigtail Probe with suture holes,
length 14.8cm

L-0470



Simcoe Cortex Extractor, 0.3mm side port,
23G

L-0480



Anis Cortex Aspirating Cannula, straight,
0.4mm, 23G, textured tip

L-0520



Anis Cortex Aspirating Cannula, angled,
0.4mm, 23G, textured tip

L-0530



Gills Aspirating Cannula,
5.0mm tip, 25G

L-0560



Gills Aspirating Cannula,
7.0mm tip, 25G

L-0570




Gills Aspirating Cannula,
8.0mm tip, 25G

L-0580



Gills Aspirating Cannula,
10.0mm tip, 25G

L-0590



Lacrimal Cannula, straight, tapered,
23G

L-0600



Lacrimal Cannula, curved, tapered,
23G

L-0610



Lacrimal Cannula, angled, 10.0mm conical

L-0620



Lacrimal Cannula, strong curved, conical

L-0630



McIntyre Lacrimal Cannula, straight, 0.3mm side port, 23G

L-0640



McIntyre Lacrimal Cannula, curved, 0.3mm side port, 23G

L-0650



Anel Lacrimal Cannula, straight,
23G

L-0660



Anel Lacrimal Cannula, curved,
23G

L-0670




Bangerter Lacrimal Probe Cannula,
40.0mm, 21G

L-0680



Bangerter Lacrimal Probe Cannula,
60.0mm, 21G

L-0690



Rycroft Air Injection Cannula, angled, 5.0mm, 27G L-0700



Rycroft Air Injection Cannula, angled, 5.0mm, 30G L-0710




Air Injection Cannula, angled, 7.0mm, 27G L-0720



Air Injection Cannula, angled, 7.0mm, 30G L-0730



Rainin Injection Cannula, angled, 7.0mm, 27G L-0740



Hydrodissection Cannula, angled, 8.0mm, flat, front opening 25G L-0750



Seeley Hydrodissection Cannula, angled, 12.0mm, flat, front opening 26G

L-0760



Pearce Hydrodissection Cannula, J-shaped, straight, front opening 25G & 22.0mm tube

L-0780




Pearce Hydrodissection Cannula, J-shaped, 8.0mm, angled, front opening 25G, 22.0mm tube

L-0790



Frankfurt Model Hydrodissections Cannula, 27G

L-0800



Sauter Hydrodissection Cannula, vertically & flat tip,
front opening 27G

L-0810




McIntyre-Binkhorst Cannula, right, blunt tip,
26G

L-0820



McIntyre-Binkhorst Cannula, left, blunt tip,
26G

L-0840



Bishop-Harmon Anterior Chamber Irrigator, complete,
(Cannula, adaptor and silicone bulb),
19G

L-0850

silicone bulb & adaptor

L-0860

adaptor

L-0870

Cannula,
19G

L-0880

Cannula,
21G

L-0890



McIntyre Anterior Chamber Irrigating Cannula, straight,
15.0mm long tube, blunt tip,
30G

L-0900



McIntyre Anterior Chamber Irrigating Cannula, straight,
15.0mm long tube, blunt tip,
26G

L-0910



McIntyre Anterior Chamber Irrigating Cannula, straight,
15.0mm long tube, blunt tip,
23G

L-0920



McIntyre Anterior Chamber Irrigating Cannula, angled,
15.0mm long tube, blunt tip,
30G

L-0930



McIntyre Anterior Chamber Irrigating Cannula, angled,
15.0mm long tube, blunt tip,
26G

L-0940



McIntyre Anterior Chamber Irrigating Cannula, angled,
15.0mm long tube, blunt tip,
23G

L-0950



Knolle Anterior Chamber Irrigating Cannula, angled, 2.0mm, front opening, 27G

L-0960



Knolle Anterior Chamber Irrigating Cannula, angled, 2.0mm, front opening, 23G

L-0970




Heslin Anterior Chamber Maintaining Cannula, 25G

L-0980




Heslin Anterior Chamber Maintaining Cannula, 23G

L-0990



Jensen Capsule Polisher, front opening, sandblasted, 23G

L-1010



Jensen Capsule Polisher, front opening, sandblasted, 25G

L-1020



Jensen Capsule Polisher, front opening, polished,, 25G

L-1021



Jensen Posterior Capsule Polisher, olive tip, front opening sandblasted, 23G

L-1040



Gills Welsh Olive Tip Cannula , 25G L-1041




Kratz Capsule Polisher, 22G L-1050



Rubman-Katzin Capsule Polisher, 22G L-1060



Capsule Polishing Cannula with bulbous tip, 18G L-1070




Knolle-Pearce Irrigating Vectis Cannula,
with three front opening ports,
23G Loop

L-1080



Sheets Irrigating Vectis Cannula,
with three front opening ports, serrated,
20G Loop

L-1090



Fukasaku Anaesthesia Cannula, gently curved,
24G

L-1120



Buratto LASIK Cannula, on front and two side ports,
25G

L-1130



Buratto LASIK Cannula, on front and two side ports,
23G

L-1140



Viscoelastic Aspirating Cannula, angled, 9.0mm, 23G L-1150



Viscoelastic Aspirating Cannula, angled, 9.0mm, 21G L-1160



Viscoelastic Aspirating Cannula, angled, 10.0mm, 23G L-1170



Viscoelastic Aspirating Cannula, angled, 10.0mm, 21G L-1180



Viscoelastic Aspirating Cannula, angled, 10.0mm, 22G L-1200



Viscoelastic Aspirating Cannula, angled, 10.0mm, 19G L-1210




McIntyre Cystotome Cannula, 10.0mm, angled tube,
23G

L-1260



Kelman Cannula, 30.0mm, straight tube,
with sharp point & sharp cutting edge,
22G

L-1270



Randolph Cyclodialysis Cannula, angled, 13.0mm
19G

L-1280





Detachable I/A Instrument complete

Instruments complete = titanium handle + sealing ring + stainless steel tip



Irrigation Handpiece,
2.0 x 0.5mm, 21G, polished

L-1311



Irrigation Handpiece,
2.0 x 0.3mm, 23G, rough

L-1321




Aspiration Handpiece,
0.3mm, 21G, polished

L-1331



Aspiration Handpiece,
0.5mm, 21G, rough

L-1341



Brauweiler Aspiration Handpiece,
0.5mm, 21G, rough

L-1351



Brauweiler Irrigation Handpiece,
0.8mm, 21G, rough

L-1361



Tip only, for detachable I/A system



Tip for Irrigation Handpiece L-1311
2.0 x 0.5mm, 21G, polished

L-1311-TIP



Tip for Irrigation Handpiece, L-1321,
2.0 x 0.3mm, 23G, rough

L-1321-TIP



Tip for Aspiration Handpiece, L-1331,
0.3mm, 21G, polished

L-1331-TIP



Tip for Aspiration Handpiece, L-1341,
0.5mm, 21G, rough

L-1341-TIP



Tip for Aspiration Handpiece, L-1351,
0.5mm, 21G, rough

L-1351-TIP



Tip for Brauweiler Irrigation Handpiece, L-1361,
0.8mm, 21G, rough

L-1361-TIP



Titanium handles only,
for detachable I/A system


Aspiration Handpiece only

L-1300-A



Irrigation Handpiece only

L-1300-I



Gills Aspirating Cannula,
23G, with silicone tubing

L-1400



Simcoe Irrigating/ Aspirating Cannula,
23G

L-1410



Simcoe Irrigating/ Aspirating Cannula,
reversed, 23G

L-1420



Gess Irrigating-Aspirating Cannula, curved right,
with end opening aspiration port,
23G

L-1430



Gess Irrigating-Aspirating Cannula, curved left,
with end opening aspiration port,
23G

L-1440



Infusions Plug, 20G,
45° bevelled tip, 2.5mm long tip

L-1450




Monomanual Phaco Safe-Lock I/A Handpiece,
with male/female luer slip fittings, compatible with
all I/A and phaco systems

L-1700



Safe - Lock I/A tip, 16G, straight,
for standard aspiration procedures with 0.35mm aspirating port

L-1710



Safe - Lock I/A tip 16G, straight,
for removal of viscoelastic fluids with 0.5mm aspirating port

L-1720



Safe - Lock I/A tip, 16G, angled 45°,
for peripheral cortex removal

L-1730




Safe - Lock I/A tip, 16G, angled 90°,
for peripheral and superior cortex removal

L-1740



Safe - Lock I/A tip, 16G, curved,
for initial cortex removal

L-1750




O-Ring Replacement Set,
for I/A tips, 5 each rings 2x1, 3x1
and 5x1mm, (pkg/15)

L-1760

All products on this page:



Bürki inno med AG
Widnau, Switzerland



Safe - Lock I/A SI-Sleeve tip, 16G, straight,
without soft style si-sleeve

L-1770



Safe - Lock I/A SI-Sleeve tip, 20G, straight,
without soft style si-sleeve

L-1780



Safe - Lock I/A SI-Sleeve tip, 45°, 16G, straight,
without soft style si-sleeve

L-1790



Silicone Irrigation Sleeve,
non-sterile
used for L-1770 - L-1790

L-1800

All products on this page:



Bürki inno med AG
Widnau, Switzerland



Luer-Lock Adaptor,
2x Luer male

L-1820

MARKING INSTRUMENTS



Castroviejo Caliper, 20.0mm,
straight


*M-0010

Castroviejo Caliper, 20.0mm,
curved

*M-0020

Castroviejo Caliper, 20.0mm,
with 0.5mm scale, straight

*M-0030



Jameson Caliper, adjustable,
0-80mm scale in 1.0mm increments


*M-0040

* Not available in Ceramic Coating



Ruler, stainless steel *M-0060

*M-0060:
 Bayerische Maßindustrie A. Keller GmbH
Hersbruck, Germany




Braunstein Measuring Instrument,
3.0mm/ 3.5mm *M-0070

Braunstein Measuring Instrument,
3.0mm/ 4.0mm *M-0080

Braunstein Measuring Instrument,
3.5mm/ 4.0mm *M-0090

* Not available in Ceramic Coating



Schiotz Tonometer, autoclavable, straight scale,
incl. certificate

MIS-0011

Schiotz Tonometer, autoclavable, curved scale,
with metrological controll report "add M to the article number"


MIS-0020
MIS-0020M

All products on this page:



Biro Medizintechnik
Burladingen-Hausen im Killertal, Germany

MISCELLANEOUS



Carter Sphere Introducer,
length 12.5cm

MIS-0010



Steri Oil Pen for lubrication of instruments with articulations,
12ml, 0,4 fl. oz., 7.0mm tip

MIS-0770

MIS-0770:



Ahrens-Innovation-Design
Wernigerode, Germany



Scleral Eye Shield (pair),
large, stainless steel

MIS-0142


Scleral Eye Shield (pair),
small, stainless steel

MIS-0143



Metal Bole, Ø 58.0mm,
height 40.0mm

MIS-0150




Sterilizing Case, autoclavable plastic, base, lid & silicone mat, outside dimension 38x68x19mm	MIS-0160
Sterilizing Case, autoclavable plastic, base deep, lid & silicone mat, outside dimension 38x68x25mm	MIS-0170
Sterilizing Case, autoclavable plastic, base deep, lid & silicone mat, outside dimension 38x68x32mm	MIS-0180
Sterilizing Case, autoclavable plastic, base & silicone mat, outside dimension 200x110x25mm	MIS-0190
Sterilizing Case, autoclavable plastic, base, silicone mat, outside dimension 265x160x25mm	MIS-0200
Sterilizing Case, autoclavable plastic, base, double & silicone mat, outside dimension 265x160x50mm	MIS-0210
Sterilizing Case, autoclavable plastic, base lid & silicone mat, outside dimension 64x150x19mm	MIS-0220

All products on this page:



PST Corporation
Pittsfield, USA



Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat ,
outside dimension 64 x150x32mm

MIS-0230

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 64x190x19mm

MIS-0240

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 102x62x19mm

MIS-0250

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 102x190x38mm

MIS-0260

Sterilization Case, autoclavable plastic,
medium base, lid & silicone mat,
outside dimension 152x254x38mm

MIS-0270

Phaco Tray, autoclavable plastic,
base, lid & 2 bars / 1 Slot,
outside dimension 100x19x38mm


MIS-0280

Phaco Tray, autoclavable plastic,
base, lid & 2 bars / 1 Slot
outside dimension 150x25x38mm

MIS-0290

All products on this page:

 PST Corporation
Pittsfield, USA



Phaco Trays, Base, autoclavable plastic,
lid, two silicone bars / 3 Slot,
outside dimension 150x25x38mm

MIS-0300

Knife Tray, autoclavable plastic,
1 slot-base, lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0310

Knife Tray, autoclavable plastic,
2 slot-base,lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0320

Knife Tray, autoclavable plastic,
3 slot-base, lid, two silicone bars,
outside dimension 64x152x19mm

MIS-0330

Knife Tray, autoclavable plastic,
5 slot-base, lid, 2 silicone bars,
outside dimension 102x162x19mm

MIS-0340

All products on this page:

 PST Corporation
Pittsfield, USA



Double Sterilisation Case KSI 1, stainless steel,
with lid & bottom perforated, 2 Silicone mat, paper filter,
outside dimension 300x220x75mm

MIS-0400

Sterilisation Case KSI 1, stainless steel,
with lid & bottom perforated, double, 2 Silicone mat & filter,
outside dimension 400x290x75mm

MIS-0420



The aluminium sterilisation containers are specially designed to accommodate 1/1, 1/2 & 1/3 Trays assuring complete sterilisation, as well as safe in-house transportation and storage.

Sterilisation Case aluminium, complete with quick-release, MIS-0430
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x40

Sterilisation Case aluminium, complete with quick-release, MIS-0440
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x65

Sterilisation Case aluminium, complete with quick-release, MIS-0450
with positive action spring mechanism on the lid, bottom perforated,
outside dimension 310x190x80



Paper Filters,with chemical indicator,
disposable (pack of 100) for single use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0460

Textile Filters, reusable (set of 2) for frequent use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0470

PTFE Teflon Filters, long-life (set of 2) for permanent use.
For MIS-0430; MIS-0440; MIS-0450.

MIS-0480

All filters are suitable for all steam sterilisation methods.

All products on this page:

 ERMIS MedTech GmbH
Tuttlingen, Germany




Identification Plates ALUMINIUM (50.0 x 18.0mm),
can be engraved to your specification;
max. 16 letters

MIS-0490



Indicator label,
disposable (pack of 100)

MIS-0500



Plastic Safety Seals,
disposable (pack of 100)

MIS-0510

All products on this page:



ERMIS MedTech GmbH
Tuttlingen, Germany



Wash Tray stainless steel, complete,
with 2 silicone instrument supports
and 1 instrument fixation bar,
1/1 size: 275x178x24mm
for 16 hand instruments

MIS-0520

Wash Tray stainless steel, complete,
with 2 silicone instrument supports
and 1 instrument fixation bar,
1/2 size: 178x140x24mm
for 8 hand instruments

MIS-0530

Wash Tray stainless steel, 24 mm height,
complete with 2 silicone instrument supports
and 1 instrument fixation bar,
1/3 size: 178x90x24mm
for 5 hand instruments

MIS-0540




Silicone Finger Mats, extra soft, light blue material,
with narrow spaced fine nubbins,
275x178x24mm

MIS-0560

Silicone Finger Mats, extra soft, light blue material,
with narrow spaced fine nubbins,
178x140x24mm

MIS-0570




Sterilisation tray with lid, stainless steel, with silicone mats,
storage of up to 16 instruments,
outer dimensions: 261x219x26mm

MIS-0614



Double-level sterilization tray, stainless steel, with silicone mats,
storage of up to 50 instruments,
outer dimensions: 325x265x 42mm

MIS-0615



Sterilization tray with lid, stainless steel, electropolished,
270x175x30 mm; 4 lower and 2 upper instrument supports,
for use in autoclave

MIS-0611



Upper Instrument Support element,
for 10 Instruments

MIS-0612



Lower instrument support element,
for 10 Instruments,

MIS-0613

CONTENT

A

A-1	52
A-2	52
A-3	52
A-3	56
A-4	52
A-5	52
A-6	53
A-7	53
A-8	53
A-9	53
A-10	53
A-11	53
A-12	54
A-14	54
A-15	54
A-16	54
A-17	54
A-22	55
A-23	55
A-24	55
A-25	55
A-26	56
A-27	56
A-30	56
A-31	56
A-32	56
A-33	56
A-34	56
A-35	56
A-36	56
A-39	53
A-45	55

C

C-0010	36
C-0030	36
C-0040	36
C-0050	37
C-0060	37
C-0070	37
C-0080	37
C-0210	39
C-0220	39
C-0240	40
C-0250	40
C-0260	40
C-0280	41
C-0290	41
C-0300	41
C-0310	41
C-0320	41
C-0362	43
C-0370	44
C-0430	45
C-0440	45
C-0470	45
C-0500	46
C-0510	46
C-0520	46
C-0530	46
C-0550	47
C-0560 + size	47
C-0670	47
C-0680	48
C-0690	48
C-0700	48
C-0710	48
C-0711	48
C-0712	48

C-0714	48
C-0720	49
C-0730	49
C-0750	50

ES

ES-0010	12
ES-0020	12
ES-0030	12
ES-0050	13
ES-0060	13
ES-0070	13
ES-0080	13
ES-0090	14
ES-0100	14
ES-0122	14
ES-0130	15
ES-0140	15
ES-0150	15
ES-0160	15
ES-0170	16
ES-0180	16
ES-0190	16
ES-0200	16
ES-0210	16
ES-0220	17
ES-0230	17
ES-0240	17
ES-0250	17
ES-0260	18
ES-0270	18
ES-0280	18
ES-0285	18
ES-0287	11
ES-0290	18

ES-0300	19	ES-0760	30	F-0440	100
ES-0310	19	ES-0770	30	F-0450	100
ES-0311	11	ES-0780	30	F-0510	101
ES-0320	19	ES-0790	31	F-0520	101
ES-0330	19	ES-0800	31	F-0560	101
ES-0340	19	ES-0810	31	F-0560T	101
ES-0341	19	ES-0820	31	F-0570	101
ES-0350	20	ES-0830	31	F-0570T	101
ES-0360	20	ES-0840	31	F-0580	101
ES-0370	20	ES-0850	32	F-0580T	101
ES-0371	20	ES-0860	32	F-0590	102
ES-0380	20	ES-0870	32	F-0590T	102
ES-0390	21	ES-0880	32	F-0600	102
ES-0400	21	ES-0890	32	F-0610	102
ES-0410	21	ES-0900	32	F-0620	102
ES-0420	21	ES-0910	32	F-0650	102
ES-0430	22	ES-0920	33	F-0650T	102
ES-0440	22	ES-0930	33	F-0680	102
ES-0450	22	ES-0940	33	F-0740	103
ES-0460	22	ES-0960	34	F-0770	103
ES-0500	23	ES-1050	34	F-0830	104
ES-0520	23	ES-1060	34	F-0840	104
ES-0540	24	F		F-0850	104
ES-0550	24			F-0860	104
ES-0560	24			F-0870	104
ES-0570	24			F-0880	104
ES-0580	24			F-0890	104
ES-0590	25	F-0010	97	F-0900	105
ES-0600	25	F-0020	97	F-0930	105
ES-0610	25	F-0030	97	F-0940	105
ES-0620	25	F-0040	97	F-0950	106
ES-0630	25	F-0050	97	F-0960	106
ES-0640	26	F-0070	98	F-0970	106
ES-0650	26	F-0080	98	F-0980	106
ES-0660	27	F-0100	98	F-0990	106
ES-0670	27	F-0120	98	F-1000	106
ES-0680	27	F-0140	98	F-1003	118
ES-0690	27	F-0160	98	F-1004	118
ES-0700	28	F-0230	99	F-1010	106
ES-0710	29	F-0250	99	F-1020	107
ES-0720	29	F-0260	99	F-1030	107
ES-0725	10	F-0270	99	F-1040	107
ES-0730	30	F-0280	99	F-1050	107
ES-0740	30	F-0410	100	F-1060	107

F-1070	107	F-1780	115	F-2890	125
F-1090	107	F-1790	115	F-2890T	125
F-1120	108	F-1800	115	F-2900	125
F-1130	108	F-1810	115	F-2910	125
F-1140	108	F-1820	115	F-2920	125
F-1150	108	F-1820T	115	F-2920T	125
F-1160	108	F-1830	115	F-2930	125
F-1190	109	F-1830T	115	F-2940	125
F-1200	109	F-1840	115	F-2940T	125
F-1210	109	F-1860	116	F-2960	125
F-1220	109	F-1910	116	F-2970	126
F-1240	109	F-1920	116	F-2980	126
F-1250	109	F-1930	116	F-2981	120
F-1260	109	F-1940	116	F-2982	120
F-1290	109	F-1960	118	F-2983	120
F-1310	110	F-1970	118	F-2990	126
F-1320	110	F-1971	117	F-2990T	126
F-1380	110	F-1975	117	F-3000	126
F-1400	110	F-1976	118	F-3100	127
F-1410	110	F-1980	138	F-3110	127
F-1420	110	F-1990	117	F-3120	127
F-1510	111	F-2000	117	F-3130	127
F-1530	111	F-2421	116	F-3140	127
F-1540	111	F-2460	121	F-3151	119
F-1550	111	F-2470	121	F-3160	127
F-1560	112	F-2480	121	F-3170	127
F-1560T	112	F-2490	121	F-3190	119
F-1570	112	F-2500	121	F-3200	119
F-1580	112	F-2510	121	F-3210	128
F-1590	112	F-2520	121	F-3220	128
F-1610	112	F-2540	122	F-3250	128
F-1620	112	F-2550	122	F-3260	128
F-1630	113	F-2560	122	F-3300	128
F-1670	113	F-2570	122	F-3310	129
F-1680	113	F-2580	122	F-3320	129
F-1690	113	F-2590	122	F-3330	129
F-1710	113	F-2800	123	F-3331	129
F-1710T	113	F-2810	123	F-3340	129
F-1720	114	F-2820	123	F-3350	129
F-1720T	114	F-2830	124	F-3360	129
F-1740	114	F-2840	124	F-3390	130
F-1750	114	F-2850	123	F-3400	130
F-1760	114	F-2860	123	F-3410	130
F-1770	114	F-2880	124	F-3420	130

F-3530	131	F-4080	138	F-4575T	150
F-3540	131	F-4090	138	F-4600	151
F-3550	131	F-4100	138	F-4610	151
F-3560	131	F-4110	138	F-4640	152
F-3570	131	F-4120	139	F-4650	152
F-3590	131	F-4130	139	F-4660	152
F-3600	131	F-4140	139	F-4670	152
F-3610	132	F-4170	140	F-4680	152
F-3620	132	F-4180	140	F-4690	153
F-3630	132	F-4190	141	F-4700	153
F-3640	132	F-4200	140	F-4710	153
F-3650	132	F-4210	141	F-4720	153
F-3660	132	F-4220	141	F-4730	153
F-3670	132	F-4230	141	F-4730	154
F-3690	133	F-4240	142	F-4740	154
F-3700	133	F-4250	142	F-4750	154
F-3720	133	F-4260	142	F-4800	92
F-3730	133	F-4270	142	F-4800TC	96
F-3750	133	F-4280	143	F-4805	92
F-3760	133	F-4290	143	F-4810	92
F-3780	134	F-4300	143	F-4820	92
F-3810	134	F-4310	144	F-4880	93
F-3820	134	F-4320	144	F-4910	93
F-3840	134	F-4330	144	F-4940	93
F-3860	134	F-4340	144	F-4950	93
F-3870	135	F-4350	145	F-5000TS	95
F-3880	135	F-4360	145	F-5010TS	95
F-3890	135	F-4370	145	F-5020TS	95
F-3900	135	F-4380	145	F-6000	96
F-3910	135	F-4390	146		
F-3920	135	F-4395	145		
F-3930	135	F-4400	146		
F-3940	136	F-4405	146		
F-3950	136	F-4410	146		
F-3960	136	F-4415	146	HS-0000TS	52
F-3970	136	F-4420	147	HS-0001TS	52
F-3980	136	F-4430	147	HS-0002TS	56
F-4000	136	F-4440	147	HS-0010TS	52
F-4010	136	F-4450	148	HS-0020TS	52
F-4020	136	F-4460	148	HS-0030TS	52
F-4040	137	F-4470	148	HS-0040TS	52
F-4050	137	F-4485	148	HS-0050TS	53
F-4060	137	F-4490	149	HS-0060TS	53
F-4070	137	F-4540	149	HS-0070TS	53

HS

HS-0080TS	53	HS-0581	61	HS-1040	68
HS-0090TS	53	HS-0590	61	HS-1050	68
HS-0091TS	53	HS-0600	61	HS-1051	68
HS-0110TS	53	HS-0610	61	HS-1060	69
HS-0120TS	54	HS-0620	61	HS-1090	69
HS-0140TS	54	HS-0621	61	HS-1100	69
HS-0150TS	54	HS-0630	61	HS-1110	69
HS-0160TS	54	HS-0640	61	HS-1120	70
HS-0170TS	54	HS-0650	62	HS-1130	70
HS-0220TS	55	HS-0660	62	HS-1140	70
HS-0230TS	55	HS-0670	62	HS-1150	70
HS-0240TS	55	HS-0680	62	HS-1160	70
HS-0250TS	55	HS-0690	62	HS-1170	71
HS-0256TS	55	HS-0691	62	HS-1250	71
HS-0260TS	56	HS-0700	62	HS-1260	72
HS-0280TS	56	HS-0710	62	HS-1270	72
HS-0310TS	56	HS-0730	63	HS-1290	73
HS-0320TS	56	HS-0740	63	HS-1300	73
HS-0330TS	56	HS-0750	63	HS-1310	73
HS-0340	57	HS-0760	63	HS-1320	73
HS-0350	57	HS-0770	64	HS-1330	73
HS-0360	57	HS-0780	64	HS-1340	73
HS-0370	57	HS-0790	64	HS-1350	73
HS-0380	58	HS-0800	64	HS-1360	73
HS-0390	58	HS-0810	64	HS-1370	73
HS-0400	58	HS-0820	64	HS-1460	74
HS-0410	58	HS-0830	65	HS-1470	74
HS-0420	58	HS-0850	65	HS-1480	74
HS-0430	58	HS-0860	65	HS-1490	74
HS-0440	58	HS-0870	65	HS-1500	75
HS-0450	58	HS-0880	65	HS-1510	75
HS-0460	59	HS-0890	65	HS-1520	75
HS-0470	59	HS-0900	66	HS-1530	75
HS-0480	59	HS-0910	66	HS-1540	75
HS-0490	59	HS-0930	66	HS-1545	75
HS-0500	59	HS-0940	66	HS-1550	75
HS-0510	59	HS-0960	66	HS-1560	75
HS-0521	59	HS-0970	67	HS-1570	76
HS-0530	60	HS-0980	67	HS-1571	76
HS-0540	60	HS-0990	67	HS-1580	76
HS-0550	60	HS-1000	67	HS-1600	77
HS-0570	60	HS-1010	67	HS-1610	77
HS-0571	60	HS-1022	67	HS-1620	78
HS-0580	61	HS-1030	68	HS-1640	79

HS-1650	79	L-0090	202	L-0660	215
HS-1660	79	L-0100	202	L-0670	215
HS-1663	79	L-0110	202	L-0680	215
HS-1664	79	L-0120	202	L-0690	215
HS-1700	80	L-0130	203	L-0700	216
HS-1710	81	L-0140	203	L-0710	216
HS-1720	81	L-0145	203	L-0720	216
HS-1730	81	L-0150	203	L-0730	216
HS-1750	81	L-0160	203	L-0740	216
HS-1790	82	L-0170	204	L-0750	216
HS-1800	82	L-0180	204	L-0760	217
HS-1810	82	L-0200	205	L-0780	217
HS-1830	82	L-0220	205	L-0790	217
HS-1840	82	L-0290	206	L-0800	217
HS-1880	83	L-0300	207	L-0810	218
HS-1930	83	L-0310	207	L-0820	218
HS-1940	83	L-0320	207	L-0840	218
HS-1950	84	L-0330	208	L-0850	219
HS-1960	84	L-0340	208	L-0860	219
HS-1990	84	L-0350	208	L-0870	219
HS-2000	85	L-0360	208	L-0880	219
HS-2010	85	L-0370	208	L-0890	219
HS-2020	85	L-0380	208	L-0900	220
HS-2021	85	L-0390	208	L-0910	220
HS-2040	86	L-0400	209	L-0920	220
HS-2060	86	L-0420	209	L-0930	220
HS-2080	86	L-0430	209	L-0940	220
HS-2100	87	L-0440	210	L-0950	220
HS-2110	87	L-0450	210	L-0960	221
HS-2120	88	L-0460	210	L-0970	221
HS-2130	88	L-0470	211	L-0980	221
L		L-0480	211	L-0990	221
		L-0520	212	L-1010	222
		L-0530	212	L-1020	222
		L-0560	213	L-1021	222
		L-0570	213	L-1040	222
L-0001	198	L-0580	213	L-1041	223
L-0002	198	L-0590	213	L-1050	223
L-0010	200	L-0600	213	L-1060	223
L-0020	200	L-0610	213	L-1070	223
L-0030	200	L-0620	214	L-1080	224
L-0040	200	L-0630	214	L-1090	224
L-0050	201	L-0640	214	L-1120	225
L-0080	201	L-0650	214	L-1130	225

L-1140	225	M		MIS-0420	250	
L-1150	226			MIS-0430	251	
L-1160	226			MIS-0440	251	
L-1170	226		M-0010	240	MIS-0450	251
L-1180	226		M-0020	240	MIS-0460	252
L-1200	226		M-0030	240	MIS-0470	252
L-1210	226		M-0040	240	MIS-0480	252
L-1260	227		M-0060	241	MIS-0490	253
L-1270	228		M-0070	241	MIS-0500	253
L-1280	228		M-0080	241	MIS-0510	253
L-1300-A	232	M-0090	241	MIS-0520	254	
L-1300-I	232			MIS-0530	254	
L-1311	230			MIS-0540	254	
L-1311-TIP	231	MIS		MIS-0560	255	
L-1321	230			MIS-0570	255	
L-1321-TIP	231			MIS-0611	257	
L-1331	230		MIS-0010	244	MIS-0612	257
L-1331-TIP	231		MIS-0011	242	MIS-0613	257
L-1341	230		MIS-0020	242	MIS-0614	256
L-1341-TIP	231		MIS-0020M	242	MIS-0615	256
L-1351	230		MIS-0142	245	MIS-0770	244
L-1351-TIP	231		MIS-0143	245		
L-1361	230		MIS-0150	245	N	
L-1361-TIP	231	MIS-0160	247			
L-1400	233	MIS-0170	247			
L-1410	233	MIS-0180	247			
L-1420	233	MIS-0190	247	N-0010		160
L-1430	234	MIS-0200	247	N-0020		160
L-1440	234	MIS-0210	247	N-0030		160
L-1450	234	MIS-0220	247	N-0040		160
L-1700	235	MIS-0230	248	N-0050		160
L-1710	235	MIS-0240	248	N-0060		160
L-1720	235	MIS-0250	248	N-0070	160	
L-1730	235	MIS-0260	248	N-0080	160	
L-1740	235	MIS-0270	248	N-0090	161	
L-1750	235	MIS-0280	248	N-0100	161	
L-1760	235	MIS-0290	248	N-0110	161	
L-1770	236	MIS-0300	249	N-0120	161	
L-1780	236	MIS-0310	249	N-0130	161	
L-1790	236	MIS-0320	249	N-0140	161	
L-1800	236	MIS-0330	249	N-0150	161	
L-1820	237	MIS-0340	249	N-0160	161	
		MIS-0400	250	N-0250	162	

N-0280	163	S		S-0420	180
N-0290	163			S-0430	180
N-0300	163			S-0440	180
N-0310	163	S-0010	173	S-0450	180
N-0320	163	S-0020	173	S-0460	180
N-0330	163	S-0030	173	S-0470	180
N-0340	164	S-0040	173	S-0480	180
N-0350	164	S-0050	173	S-0490	181
N-0360	164	S-0060	173	S-0501	181
N-0370	164	S-0070	174	S-0505	181
N-0380	164	S-0080	174	S-0510	181
N-0390	164	S-0090	174	S-0530	182
N-0400	165	S-0100	175	S-0540	182
N-0410	165	S-0110	175	S-0570	183
N-0420	165	S-0120	175	S-0580	183
N-0430	165	S-0130	175	S-0581	183
N-0440	165	S-0140	175	S-0590	183
N-0450	165	S-0150	175	S-0600	183
N-0460	166	S-0160	176	S-0610	183
N-0470	166	S-0170	176	S-0611	183
N-0480	166	S-0180	176	S-0620	183
N-0490	166	S-0190	176	S-0621	184
N-0500	167	S-0200	177	S-0630	184
N-0510	167	S-0210	177	S-0640	184
N-0520	167	S-0220	177	S-0641	184
N-0530	167	S-0220SC	190	S-0650	184
N-0540	168	S-0230	177	S-0660	185
N-0550	168	S-0240	177	S-0670	185
N-0560	168	S-0250	177	S-0670SC	189
N-0570	168	S-0250SC	190	S-0680	185
N-0630	168	S-0260	178	S-0680SC	189
N-0660	169	S-0270	178	S-0690	186
N-0690	170	S-0280	178	S-0700	186
N-0700	170	S-0290	178	S-0710	186
N-0800	157	S-0300	178	S-0720	186
N-0800TC	159	S-0310	178	S-0730	187
N-0810	157	S-0380	179	S-0740	187
N-0810TC	159	S-0382	177	S-0750	187
N-0820	157	S-0390	179	S-0760	187
N-0830	157	S-0400	179	S-0770	188
N-0830TC	159	S-0410	179	S-0780	188
N-0870	158	S-0415	182	S-0780SC	189
		S-0416	182	S-0790	188

VR

VR-0010-23T	192
VR-0030-23T	192
VR-0050-23T	192
VR-0090-23T	192
VR-0110-23T	192
VR-0205-23T	192
VR-0220 -23T	193
VR-0240-23T	193
VR-0250-23T	193
VR-0270	193
VR-0280	193
VR-0300-23TDT	195
VR-0300-23-TIP	195
VR-0310-23TDT	195
VR-0310-23-TIP	195
VR-0400-23TDT	195
VR-0400-23-TIP	195
VR-0410-23TDT	195
VR-0410-23-TIP	195
VR-1600	194
VR-1610	194
VR-1620	194
VR-5000T	195

ALPHABETICAL INDEX

A

Adjustable Eye Speculum 19
Anaesthesia Cannula 225
Anis Needle Holder 166
Anterior Chamber Irrigating Cannula 220, 221
Anterior Chamber Irrigator 219
Anterior Chamber Maintaining Cannula 221
Aspirating Cannula 212, 213, 226, 233, 234
Aspiration Handpiece 230, 232
Aspiration Handpiece only 232
Axis Marker 76
Ayer Chalazion Forceps 141

B

Baby Jones 152
Backhaus towel clamp 153
Baird Chalazion 143
Baird Chalazion Forceps 143
Bangerter Muscle Clamp 147
Bard-Parker Handle 47, 48
Barraquer Cilia Forceps 137, 138
Barraquer Fixation Forceps 137
Barraquer Iris Scissors 182
Barraquer Iris Spatula 55, 59, 64
Barraquer Needle Holder 157, 158, 159, 164, 165
Barraquer Wire Speculum 12, 13, 14
Barrett Eye Speculum 11
Bechert Nucleus Rotator 53, 56, 61, 68
Berke Ptosis Forceps 144
Beyer Rongeur 201
Bimanual I/A 5, 197
Bipolar Cord 153, 154
Bipolar Forceps 5, 153, 154

Bishop-Harmon Anterior Chamber Irrigator 219
Bishop Harmon Fixation Forceps 131
Bishop-Harmon Tissue Forceps 132
Blade Breaker & Holder 170
Blaydes Implantation Forceps 149
Bonaccolto Utility Forceps 133
Bonn Forceps 93, 107, 116, 126, 127
Bores Optic Zone 83, 84
Bores Optic Zone Marker 83, 84
Braunstein Measuring Instrument 241
Brauweiler Aspirating Handpiece 230
Breaker and Holder 169
Brown Triple 21
Bunge 74
Burrs 5, 35

C

Caliper 240
Cannula 111, 212, 213, 214, 215, 216, 217, 218, 219, 220, 221, 223, 224, 225, 226, 227, 228, 233, 234
Cannulas 5, 197
Capsule Polisher 222, 223
Capsulorhexis Forceps 93, 106, 107, 108, 115, 118, 119, 120, 121, 125, 126, 127, 128, 129, 150
Capsulotomy Scissors 180, 181
Carter Sphere Introducer 244
Case 248, 250, 251
Castroviejo Breaker and Holder 169
Castroviejo Caliper 240
Castroviejo Corneoscleral Scissors 175
Castroviejo Eye Speculum 22
Castroviejo Keratoplasty Scissors 176
Castroviejo Lacrimal Dilator 202
Castroviejo Needle Holder 161
Castroviejo Spatula 69

Castroviejo Suturing Forceps 109, 110, 112, 122

Castroviejo Tying Forceps 108

Castroviejo Universal Corneal Scissors 173, 176

Chalazion Forceps 141, 142, 143, 144

Chalazions Curette 75

Chang Phaco Chopper 56, 69

Chisel 205, 206, 207

Chisels 5, 9

Chuck Handle 48

Cilia Forceps 137, 138

Citelli Bone Nibbling Rongeur 201

Colibri Forceps 97, 98, 99, 110, 113, 128

Combo Pre-Chopper 130

Conjunctiva Fixation Forceps 108, 121, 133

Conjunctiva Forceps 116, 134, 140

Conjunctival Scissors 177

Cook Eye Speculum 17

Corneal Dissector 37, 67

Corneal Fixation Forceps 98, 102, 103

Corneal Knife 36

Corneal Rust Ring Remover 41

Corneal Scarifier 37

Corneal Scarifiers 35

Corneal Transplant Scissors 179

Corneal Utility Forceps 97, 98

Corneoscleral Scissors 175

Crawford Retrieving Hook 210

Cyclodialysis Cannula 228

D

DALK Corneal Dissector 67

Depressor and Fixation 34

Derf Needle Holder 170

Desmarres Chalazion Forceps 142

Desmarres cross-action Forceps 144

Desmarres Lid Retractor 25

DeWecker Iridectomy Scissors 181

Dieffenbach 148

Direct Action Forceps 93, 115

Dissector 35, 37, 67

Double Ended Strabismus Hook 33

Double hook 56

Double-level tray 256

Downs Lid Forceps 145

Dressing Forceps 131, 135, 136

Drysdale Nucleus Manipulator 53, 61

E

Elevator 203, 204

Elschnig Skleral Fixation Forceps 133

Endgripping Forceps 192

Endothelial / Descemetorhexis Forceps 117

Entropium Forceps 145

Enucleation Scissors 187

Enucleation Spoon 72

Ergo Grip 92

Ergo Grip Barraquer Needle Holder 157, 158, 159

Eye Scissors 183, 184

Eye Shield 245

Eye Speculum 14, 15, 17, 18, 19, 20, 21, 22

F

Fascia Needle 9, 26

Femto Flaplifter 68

Ferris-Smith-Kerrison Punch 198

Fine Stitch Scissors 174

Fishkind Eye Speculum 22

Fixation Forceps 98, 102, 103, 108, 110, 112, 113, 121, 131, 133, 137

Fixation Hook 29

Fixation Ring 23, 77, 78

Fixation Rings 5, 9

Flieringa Scleral Fixation Ring 23

Foreign Body Forceps 141, 193

Foreign Body Needle 39, 40

Foreign Body Spud 40

Fragment Forceps 114, 125, 127

Francis Chalazion Forceps 143

Freer Lacrimal Chisel 207

Freer Periosteal Elevator 203

G

Gaskin Fragment Forceps 114, 125, 127
Gass Retinal Detachment Hook 32
Gass Scleral Marker 24
GIANNETTI MIC Capsulorhexis Forceps 118, 120
Gills-Vannas Capsulotomy Scissors 180, 181
Gills-Welsh Capsulotomy Scissors 181
Gills Welsh Olive Tip Cannula 223
Ginsberg Eye Speculum 20
Globe Retractor 72
Gradle Cilia Forceps 138
Graefe model 137
Graefe Muscle Hook 31
Graether Collar Button Micro Iris Hook 64
Grasping Forceps 194
Green Muscle Hook 32

H

Haefliger Phaco cleaver 52, 57
Halsted Mosquito Forceps 152
Harms Trabeculotomy Probe 207
Hartmann Mosquito Forceps 152
Helveston Muscle Hook 31
Hokey Knife 36
Hooks 5, 9, 51
Hunt Chalazion Forceps 142
Hydrodissection Cannula 216, 217, 218

I

ICL Manipulator 59
ILM Forceps 195
Implantation Forceps 115, 149
Inamura Capsulorhexis Forceps 150
Infant Lid Speculum 16
Infusions Plug 234
Injection Cannula 216
Instruments 5, 58, 91, 156, 197, 239, 257
IOL Cutter 5, 50

IOL Cutters 35
IOL Implantation Forceps 115
IOL Inserting Forceps 123
IOL-Manipulator 53, 54, 62
Iridectomy Scissors 181
Iris Forceps 100, 101, 122, 135, 136
Iris Hook 52, 58, 59, 61, 64, 65, 66
Iris Retractor 34
Iris Scissors 178, 182
Iris Spatula 55, 59, 64, 70
Irrigating-Aspirating Cannula 234
Irrigating Vectis Cannula 224
Irrigation Handpiece 230, 232
Irrigation Handpiece only 232

J

Jaeger Lid Plate 26
Jaeger Strabismus Hook 33
Jaffe Capsulorhexis Forceps 115, 118, 125
Jaffe Forceps 117
Jaffe Stitch Scissors 174
Jaffe Wire Lid Retractor 25
Jakobi Needle Holder 167
Jameson Caliper 240
Jameson Muscle Forceps 146
Jameson Muscle Hook 32
Jayle Forceps 134
Jensen Capsule Polisher 222
Jeweler's Forceps 139
Joseph Periosteal Elevator 204

K

Kalt Needle Holder 162
Kelman Cannula 228
Kelman-McPherson Forceps 92, 113, 116
Keratometer 5, 88
Keratoplasty Scissors 176
Kerrison Bone Nibbling Rongeur 200
Knapp Strabismus Scissors 187

Knife Tray 249
Knives 5, 35
Knolle-Pearce Irrigating Vectis Cannula 224
Koch Stop & Chop Manipulator 59
Kraff Eye Seculum temporal 20
Kratz-Barraquer Wire Speculum 14
Kratz Cystotome Cannula 227
Krumeich Eye Speculum 11
Kuglen Iris Hook 52, 61, 66

L

Lacrimal Cannula 213, 214, 215
Lacrimal Dilator 202, 203
Lacrimal Instruments 5, 197
Lacrimal Intubations Set 209
Lacrimal Probe 208, 215
Lacrimal Probe Cannula 215
Lacrimal Sac Chisel 205
Lacrimal Sac Retractor 22, 28, 30
Lambert Chalazion Forceps 141
Lancaster Eye Speculum 19
LASEK Alcohol Well 85
LASIK Cannula 225
LASIK Corneal Marker 80
LASIK Flap Forceps 116
LASIK Flap Protector 85
LASIK Spatula 65
Lehner Capsulorhexis Forceps 93, 106, 115, 120, 128
Lens cutting Scissors 194
Lens Folding Forceps 115
Lens Forceps 114, 147, 149
Lens Holding Forceps 114, 115, 125, 151
Lens Loop 73
Lens Manipulating Hook 54, 62, 63, 66, 67
Lens Manipulator 52, 55, 59, 61, 63, 66
Lens Spatula 67
Lester Conjunctiva Fixation Forceps 133
Lester IOL-Manipulator 53, 54, 62
Lid Clamp 146
Lid Forceps 145
Lid Plate 26

Lid Retractor 25
Lid Speculum 16
Lieberman Aspirating Speculum 15
Lieberman Eye Speculum 15, 19, 20, 21
Ligabue Nucleus Manipulator & Splitter 55
Ligature Scissors 188
Lister Conjunctiva Forceps 134
Loop 73, 224
Lordan Chalazion Forceps 144
Luer-Lock Adaptor 237

M

Mallet 204
Manipulators 5, 51
Marking Instruments 239
Marking Pen 5, 88
Masket Phaco Spatula 56, 69
McIntyre-Binkhorst Cannula 218
McNeil-Goldman Scleral Fixation Ring 23
McPherson Needle Holder 160
McPherson-Westcott Stitch Scissors 173
Meibomian Gland Compressor 117
Mellinger Eye Speculum 18
Membrane Hook 86
Membrane Scratcher 86
Membran-Peeler 86
Mendez Degree Gauge 76
Metal Bole 245
Metzenbaum Scissors 188, 189
Meyerhoefer Chalazions Curette 75
MIC Corydon Capsulorhexis Forceps 118, 119
MIC Jaffe Capsulorhexis Forceps 118
Micro Colibri Forceps 98
Micro Finger 65
Micro Gripping Forceps 195
MIC Utrata Capsulorhexis Forceps 118
Miscellaneous 5, 243
Monomanual Phaco Safe-Lock I/A Handpiece 235
Moorfields Conjunctiva Forceps 140
Mosquito Forceps 152
Murdoch Eye Speculum 18

Muscle Clamp 147
Muscle Forceps 146
Muscle Hook 31, 32

N

Nagahara Phaco-Chopper 52, 53, 57
Nasal Speculum 202
Needle & Gouge 27
Needle Holder 155, 157, 158, 159, 160, 161, 162, 163, 164, 165,
166, 167, 170
Needle Holder with TC-Dust 159
Needles 5, 35
Neuhann Nucleus Rotator 58
Nevyas Capsulorhexis 93
Nevyas Capsulorhexis Forceps 93
Noyes Iris Scissors 178
Nucleus Cracker 111
Nucleus Cracking Forceps 124
Nucleus Divider 58
Nucleus Loop 73
Nucleus Manipulator 53, 55, 61, 65
Nucleus Manipulator & Splitter 55
Nucleus Rotator 53, 56, 58, 61, 68
Nucleus Rotator & Dividing Hook 68
Nucleus Rotator (push-pull) & Dividing Hook 56, 68
Nugent Utility Forceps 136

O

Ocular Landers Lens Forceps 147
Oil Pen 244
Optic Zone Marker 83, 84
Orbital Spatula 71
O-Ring Replacement Set 235

P

Paper Filters 252
Paper Label/Indicator 253
Partsch Lacrimal Sac Chisel 205

Partsch Mallet 204
Paufigue Suture Forceps 104, 105
Phaco-Chopper 52, 53, 57, 63
Phaco Knife 44
Phaco Pre-Chopper 123
Phaco Splitter 52, 61, 69
Phaco Trays 248, 249
Pickup & Folding Forceps 124
Pigtail Probe 211
Pisacano Nucleus Rotator 53, 61
Plastic Safety Seals 253
Pressure Plate 87
Pterygium Drill 41
PTFE Teflon Filters 252
Ptosis Forceps 144
Punches 5, 35
Push-pull 52, 56, 61, 66, 68
Putterman Forceps 146

R

Radial Marker 79
Remky Superior Rectus Forceps 134
Removal Cannula Forceps 111
Replacement burr 41
Replacement tip 52, 53, 54, 55, 56
Retractors 5, 9
Retrieving Hook 210
Retrograde Probe 209
Ribbon type 184, 186
Rongeur 200, 201
Rosen Phaco Splitter 52, 61, 69
Rowen Phaco Spatula & Bechert Nucleus Rotator 56, 68
Ruler 241

S

Sac Retractor 22, 28, 30
Safe - Lock I/A 235, 236
Safe - Lock I/A tips 235
Sauer 16, 105
Schepens Scleral Depressor 24

Schocket Scleral Depressor 24

Scissors 5, 171, 172, 173, 174, 175, 176, 177, 178, 179, 180, 181, 182, 183, 184, 185, 186, 187, 188, 189, 190, 193, 194, 195

Scleral Depressors 5, 9

Scleral Knife 37

Scleral Plug Forceps 147

Serrifine 148

Shaaf Cilia Foreign Body Forceps 141

Shepard Lens Forceps 149

SIGNUM GIANETTI MICS Forceps 95

SIGNUM INAMURA MICS Forceps 95

SIGNUM LEHNER MICS Forceps 95

Silicone Finger Mats 255

Silicone Irrigation Sleeve 236

Simcoe Cortex Extractor 211

Simcoe Irrigating - Aspirating Cannula 233

Sinsky II Lens Manipulating Hook 54, 62, 67

Skleral Fixation Forceps 133

Sleeve Spreader 140

Sleeve tip 236

SMILE GRIP IOL Loading Forceps 96

Snapper Hook 65

Snellen Entropium Forceps 145

Spatula 51, 55, 56, 58, 59, 60, 64, 65, 67, 68, 69, 70, 71, 73

Speculum 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 202

Spencer Ligature Scissors 188

Spoon 5, 72, 73

Spud & Gouge 35

Sterilisation Case 250, 251

Sterilisation Case aluminium 251

Sterilisation cassette 256

Steri Oil Pen 244

Stevenson Lacrimal Sac Retractor 22

Stevens Tentotomy Scissors 185, 189

Stitch Scissors 174, 184, 185

Stop & Chop Manipulator 59

Strabismus Hook 33

Strabismus Scissors 186, 187

Super-Cut 189

Superior Rectus Forceps 134

Suturing Forceps 92, 101, 102, 104, 105, 109, 110, 112, 114, 122, 123

Synechiaie Spatula 67

T

„Teaser“ Hook 32

Temporal Approach 13

Tennant Tying Forceps 92

Tenotomy Scissors 173, 177, 179, 186, 190

Textile Filters 252

Thornton Fixation Ring 77

Thorpe Conjunctiva Fixation Forceps 108

Tip only 195

Tissue Forceps 132

Tissue Retractor 27

Titanium 52, 88, 91, 97, 101, 102, 112, 113, 114, 115, 125, 126, 150, 156, 165

Tonometer 242

Towel clamp 152, 153

Trabeculectomy Marker 85

Trabeculectomy Punch 49

Trabeculotomy Probe 207

Trephine Blade 47

Trephine Handle 47

Trephines 5, 35

Troutman-Chris Blade Breaker & Holder 170

Troutman-Katzin Corneal Transplant Scissors 179

Troutman Tying Forceps 123

Tying Forceps 92, 96, 101, 106, 107, 108, 111, 113, 121, 123, 125, 127, 128, 129

U

Universal Corneal Scissors 173, 176

Universal Trephine Handle 47

Uribe-Stern Capsulotomy Scissors 181

Urretts-Zavalía Depressor and Fixation 34

Utility Forceps 97, 98, 99, 133, 136

Utility Scissors 188

Utrata Capsulorhexis 93, 106, 108, 115, 118, 121, 125, 126, 127, 128, 129

Utrata Capsulorhexis Forceps 93, 106, 108, 115, 118, 121, 125, 126, 127, 128, 129

V

Vannas Capsulotomy Scissors 180, 181

Vannas Scissors 175

Viscoelastic Aspirating Cannula 226

Vitrectomy 5, 87, 191, 192, 193

Vitrectomy Forceps 192, 193

Vitrectomy Foreign Body Forceps 193

Vitrectomy Scissors 193

W

Wash Tray 254

Watzke Sleeve Spreader 140

Westcott Tenotomy Scissors 173, 177, 179, 190

Wies Chalazion Forceps 143

Wilder Lacrimal Dilator 203

Williams Eye Speculum 19

Wire Speculum 12, 13, 14, 17

Y

Y-Hook 55, 59

Z

Zone Marker 81, 82, 83, 84



MAR-0001 / v.4.0